

Product datasheet for **SC336426**

DDX52 (NM_001291476) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDX52 (NM_001291476) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDX52
Synonyms:	HUSSY19; ROK1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336426 representing NM_001291476.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCATCTGTAGAAGCAAAGATTGAAGACAAAAAGTTCAGAGAGAAAAGTAACTAACTCCGGAAG
TTGGAGAATCTCAGAAAAGAAAAGATAAACTTCTTGGGAATAAACACAAAATTCACGTCCAAGGAACC
GATCTTCTGACCCAATTGCTACATTTTCAGCAACTTGACCAGGAATATAAAATCAATTCTGACTACTT
CAGAACATTCTAGATGCAGGTTTCCAAATGCCTACGCCAATCCAATGCAAGCCATCCCAGTTATGCTG
CATGGTCGGGAATTCTGGCTTCTGCTCCAACGGATCTGAAAAACATTAGCTTTTAGCATTCTTATT
TTAATGCAGCTGAAACAACCCGCAATAAAGGCTTCAGAGCCCTGATTATACCAACACGAGAATT
GCCAGCCAGATTCACAGAGAGTTAATAAAATTTCTGAGGGAACAGGATTCAGAATACACATGATCCAC
AAAGCAGCAGTGGCAGCCAAGAAATTTGGACCTAAATCATCTAAAAGTTTGATATTCTTGACTACT
CCAAATCGACTAATCTATTTATTAAGCAAGATCCCCCGGAATCGACCTAGCAAGTGTTGAGTGGCTT
GTAGTAGCGAATCAGATAAACTGTTTGAAGATGGCAAACTGGGTTTCAGAGACCAGCTGGCTTCCATT
TTCTGGCTGCACATCCACAAGGTCCGAAGAGCTATGTTCAAGTCAACTTTTGCATATGATGTTGAA
CAGTGGTGCAAACTCAACTGGACAATGTCATCAGTGTGTCCATTGGAGCAAGGAATTCGCAGTAGAA
ACTGTAGAACAAGAGCTTCTCTTTGTTGGATCTGAGACCGGAAAACCTTCTGGCCGTGAGAGAATTGTT
AAAAAGGTTTCAATCCACCTGTTCTTGTTTTGTTCAGTCCATTGAAAGGGCTAAAGAATTTTTCAT
GAGCTCATATGAAGGTATTAATGTGGATGTTATTCATGCAGAGAGAACAACAACAGAGAGATAAC
ACAGTCCACAGTTTCAGAGCAGGAAAACTGGGTTCTGATTTGTACAGCCTTGTAGCAAGAGGGATT
GATTTTAAAGGTGTGAAGTGGTATCAACTATGACTTCCAAGTCAAGTGGATATATCCACAGG
ATAGGTCGAAGTGAAGAGCAGGGAATAAGGGAAAAGCAATTACATTTTCACTGAGGATGATAAGCCA
TTATTAAGAAGCGTTGCTAATGTTATACAGCAGGCTGGGTGTCCTGTACCAGAATACATAAAAGTTTTT
CAGAAACTACTAAGCAAAACAAAAGAAAAGATGATTAAGAAACCATTGAAAGGGAGAGCATTAGTACA
ACTCCAAATGTTTCTTAGAAAAAGCTAAGGATAAACAGAAAAGGCTCACTGGTCAGAACAGCAAGAAG
AAAGTAGCTCTTGAAGACAAAAGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001291476

Insert Size:	1476 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291476.1</u>
RefSeq Size:	6556 bp
RefSeq ORF:	1476 bp
Locus ID:	11056
UniProt ID:	<u>Q9Y2R4</u>
Cytogenetics:	17q12
MW:	55.2 kDa