

Product datasheet for **SC336433**

LETM2 (NM_001286819) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LETM2 (NM_001286819) Human Untagged Clone
Tag:	Tag Free
Symbol:	LETM2
Synonyms:	SLC55A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336433 representing NM_001286819.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTTCTACAGTTATAATTCAGTTCTGGCTATTGCTCGAACAAGATTCCTAGCCATTTTGTCCAT
CCTACCTGCTCTTCTTATTCCCCATCATGTGCATTTCTTCACTTGCCAGATTCCCATTTAAATAAGACA
TGTATGAAGAAGCTATGAGAGCAAGAAGTACTCGGATCCTAGTCAGCCAGGCAATACAGTACTTCACCCA
GGAAGTACTAATAACAAAAGCTACACACATCCACTTGCTGGCTGCAAGAAGTTCCTGGCAAACCTCAG
CTGGAGCAAGCCACAAAACATCCACAGGTGACAAGCCCTCAGGCCACAAAAGAACTGGCATGGAGATT
AAAGAAGGCAAACAATCTTATAGACAAAAATCATGGTGAAGTAAAATATTATTACAATGGATTCTAC
TACTTTGGATTGACGCCAAAGTTGCTGCCAGAATGGTTTGGAGGCTGTTGCATGGACAGGTCTGACC
AGACGAGAGAGACGAAGGCTGTTGAGAAGTGTGTTGATTTCTCCGCCTGGTCCATTTATGGTGTTC
TTAATTGTACCCTTCATGGAATTCATTACCAGTGTTCGAACTCTCCAGAGATGTTGCCATCA
ACTTTTAAAAGTGAATCCAAAAGGAAGAAAACAGAAAAGAAAATGGCTGTAAAGTTGGAAGTACGCA
AAATTTCTTCAAGAAACCATGACAGAAATGGCAAGGAGGAACAGAGCCAAGATGGCGCATGCCTCTACA
CAGCTCTCATCTACGTGAAGCAGGTCCAGACAGGCCACAAGCCAGCACAAAGGAGATAGTTTCGCTTC
TCCAAACTATTTGAGGACCAGCTGGCCCTGGAACACTTAGATCGCCCTCAGCTGGTTGCCCTTGGCAA
CTGCTGGAATTGCAGACATTTGGAACCAACAACCTGCTCCGCTTTCAGCTCCTGATGAACTGAAGTCT
ATAAAAGCAGATGATGAAATAATTGCCAAGGAAGGGGTGACAGCATTGAGTGTATCAGAAGTACAGGCT
GCCTGTAGGGCCCAGGGATGAGTACTGGGTCTCACGGAGGAACAACCTGCGACAACAGCTCACGGAG
TGGCAGGACCTCCACCTGAAGGAGAACGTCCTCCTTCCCTTTTGTCTCCTGTCCCGCACCTCTACCTG
ATAGATGTGAAGCCCAAGCCGATTGAGATACCAGTCACTGAGGAGGCTCCAAAGACTGATATTCTTGTG
GAATTACCTACTTTCAGTGAATCTAAAGAGAACATGGTGGATCTTGCACCTCAACTGAAGGGAAGTAAAG
GATGAAGACTTTATACAGCCGCCACCAGTTACATCATCACCATAACACCATCAACACCTATTTCTTATA
CCTAAAGGACCCATCACTTCTTCTGAAGAAGCTACACTCCAGGCCAAATCACAATGACGGCCAGAAC
AGCAAGGCTAGTTCAAAAGGAGCATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001286819

Insert Size:	1476 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286819.1</u>
RefSeq Size:	2891 bp
RefSeq ORF:	1476 bp
Locus ID:	137994
UniProt ID:	<u>Q2VYF4</u>
Cytogenetics:	8p11.23
MW:	55.9 kDa