

## Product datasheet for **SC336458**

### Estrogen Receptor beta (ESR2) (NM\_001291723) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Estrogen Receptor beta (ESR2) (NM_001291723) Human Untagged Clone
Tag:	Tag Free
Symbol:	ESR2
Synonyms:	ER-BETA; Erb; ESR-BETA; ESRB; ESTRB; NR3A2; ODG8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336458 representing NM\_001291723.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGATATAAAAACTCACCATCTAGCCTTAATTCCTCCTCCTACAACCTGCAGTCAATCCATCTTA
CCCTGGAGCAGGCTCCATATACATACCTTCTCCTATGTAGACAGCCACCATGAATATCCAGCCATG
ACATTCTATAGCCTGCTGTGATGAATTACAGCATTCCAGCAATGTCACTAACTTGGAAGTGGGCT
GGTCGGCAGACCAAGCCAAATGTGTTGTGGCCAACACCTGGGCACCTTTCTCCTTTAGTGGTCCAT
CGCCAGTTATCACATCTGTATGCGGAACCTCAAAAGAGTCCCTGGTGTGAAGCAAGATCGCTAGAACAC
ACCTTACCTGTAACAGAGAGACTGAAAAGGAGGTTAGTGGGAACCGTTGCGCCAGCCCTGTTACT
GGTCCAGGTTCAAAGAGGGATGCTCACTTCTGCGTGTCTGCAGCGATTACGCATCGGGATCACTAT
GGAGTCTGGTGTGAAGGATGAAGGCCTTTTTAAAAGAAGCATTCAAGGACATAATGATTATATT
TGTCAGCTACAATCAGTGTACAATCGATAAAAACCGGCAGAGCTGCCAGGCTGCCGACTTCGG
AAGTGTTACGAAGTGGGAATGGTGAAGTGTGGCTCCCGGAGAGAGAGATGTGGTACCGCCTTGTCGG
AGACAGAGAAGTGCCGACGAGCAGCTGCACTGTGCCGGCAAGGCCAAGAGAAGTGGCGGCCACGCGCC
CGAGTGGCGGAGCTGCTGCTGGACGCCCTGAGCCCCGAGCAGCTAGTGCTCACCTCCTGGAGGCTGAG
CCGCCCATGTGCTGATCAGCCGCCCACTGCGCCCTCACCGAGGCCTCCATGATGATGTCCTGACC
AAGTTGGCCGACAAGGAGTTGGTACACATGATCAGCTGGGCCAAGAAGATTCCCGCTTTGTGGAGCTC
AGCCTGTTGACCAAGTGCAGCTCTTGGAGAGCTGTTGGATGGAGGTGTTAATGATGGGGCTGATGTGG
CGCTCAATTGACCACCCCGCAAGCTCATCTTTGCTCCAGATCTGTTCTGGACAGGGATGAGGGGAAA
TGCGTAGAAGGAATCTGAAATCTTTGACATGCTCCTGGCACTACTTCAAGGTTTCGAGAGTTAAAA
CTCCAACACAAGAATATCTCTGTGCAAGGCCATGATCCTGCTCAATTCCAGTATGACCTCTGGTC
ACAGCGACCCAGGATGCTGACAGCAGCCGGAAGCTGGCTCACTTGAACGCCGTGACCGATGCTTTG
GTTTGGTGATTGCCAAGAGCGGCATCTCTCCAGCAGCAATCCATGCGCCTGGCTAACCTCCTGATG
CTCCTGTCCACGTCAGGCATGCGAGGGCAGAAAAGGCCTCTCAAACACTCACCTCATTTGGAATGAAG
ATGGAGACTCTTTGCTGAAGCAACGATGGAGCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_001291723

<b>Insert Size:</b>	1488 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001291723.1</a>
<b>RefSeq Size:</b>	2745 bp
<b>RefSeq ORF:</b>	1488 bp
<b>Locus ID:</b>	2100
<b>UniProt ID:</b>	<a href="#">Q92731</a>
<b>Cytogenetics:</b>	14q23.2-q23.3
<b>Protein Families:</b>	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
<b>MW:</b>	55.5 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the family of estrogen receptors and superfamily of nuclear receptor transcription factors. The gene product contains an N-terminal DNA binding domain and C-terminal ligand binding domain and is localized to the nucleus, cytoplasm, and mitochondria. Upon binding to 17beta-estradiol or related ligands, the encoded protein forms homo- or hetero-dimers that interact with specific DNA sequences to activate transcription. Some isoforms dominantly inhibit the activity of other estrogen receptor family members. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been fully characterized. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (l) represents an alternative promoter usage and differs in the 5' UTR, 3' coding region and 3' UTR, compared to variant a. The resulting isoform (2) has a distinct and shorter C-terminus, compared to isoform 1. Variants b, k and m encode the same isoform 2.</p>