

## Product datasheet for SC336459

### C3orf23 (TCAIM) (NM\_001282913) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C3orf23 (TCAIM) (NM_001282913) Human Untagged Clone
Tag:	Tag Free
Symbol:	TCAIM
Synonyms:	C3orf23; TOAG-1; TOAG1
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC336459 representing NM_001282913. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGTTTTGCCACCTGAGACCTATGAGGAGGTTATGTCTAGAGAAGATATTTCCACACTGGTTCCCTTT
TCAAGAGCTTTATCGGGAGCTGAAGCAGTCAATGCCTTGAGGCCTTTCTATTTTGCAGTACATCCAGAT
TTCTTTGGACAGCACCCCGTAGAAAGGGAAATCAATGAAAATTCTCTTAAAAGGTTAAGTGTCTACCTA
GAAAACCTCCAGAAACCAGGCTTCAAGTCTTTGAAACCAACTCAGCTTACATTTTATGTAAGAGAAAACA
GACCAGAGTTCCTCCGATGGCCAGGAACCTTTTAGTACTTCCGGATTTTCGAGCAGTCAAATTTACTTTG
CACACCAGAGATCTGCTAAGCACAGTGTATATATTCTCAACTCTGCAGTTTATCTGTTGAACATATC
CAAAGCTTGAATACTAATATGCATACCCAGCCTCTCAAAGAAGCTAAAAGGATGCCTGACAGGCCATC
AAATGGGACAAGTCTTATTACTCTTTACTGGATTCAAGGACCCTGATGAAGACCTTGAACAAGTCTCG
AGAGTGGAAACAACCTCTCACATCCTGGTTAGATAACAATGGGAAAAGTGCTGTTAAAAAGCTAAAGAAC
AGTTTGGCACTTAGAAAAGAAGTATGATCGTTTAAAAGATGAACTGTCTCATCAATTGCAACTCTCAGAT
ATCAGGTGGCAGAGGAGCTGGGGCATCGCCACCCTGTAGCCAGCTGCATAGTTTAAAGCCGCTTAGCA
CAGCAGAAATTTGGAAACACTTAAAAAGCAAAGGGTGTACAATCATATTTACAGACCGTTCTGGCATG
AGTGCAGTGGCCATGTGATGCTAGGAACAATGGATGTCCATCACCCTGGACAAAACCTTTTGAAGA
TTGCCAAGTTATTTTACCTTCAGAGGAGGCTGATGATTTTGAAGACCAATAAGCTATCTTTTAGGT
GGCATAACAAGTTGTTTATATTGAAGAATTACAGCCAGTATTGACACTTGAAGAATATTACTCTCTTCTT
GATGTGTTTTATAATAGACTGTTGAAAAGTAGAATACTATTTACCCTCGAAGTTTTCGCTGGTTTACAA
ATGATCCTTAACAGTGACAGATATGCTCCAAGCTTGCATGAACTCGGGCATTTTAATATTCCAACACTC
TGTGATCCAGCAAATCTCCAGTGGTTTATTCTCACAAAGCTCAGCAGGCAAGAGAGAACATGAAAAGA
AAGGAAGATTAAGGTTATTGAAAATGAATTGATACAGGCATCAACAAAGAAATTTTCTTTGGAGAAG
TTATATAAAGAGCCAGCATTCTAGTACAAATGGTGGATTGTTGTAAGAGACTTCTAGAACAAATCA
CTGCCTTACCTACATGGGATGCACCTCTGCATTTACATTTTTACTCTGTTATGCAAGATGGAGACCTT
TGTATTCCTTGAATTGGAAGAATGGAGAAGCCATTAAGTAA

```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001282913

**Insert Size:** 1491 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001282913.1](#)

**RefSeq Size:** 3520 bp

**RefSeq ORF:** 1491 bp

**Locus ID:** 285343

**UniProt ID:** [Q8N3R3](#)

**Cytogenetics:** 3p21.31

**MW:** 57.9 kDa

**Gene Summary:** May regulate T-cell apoptosis.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Both variants 1 and 4 encode the same isoform (1).