

Product datasheet for **SC336465**

SLC2A14 (NM_001286235) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | SLC2A14 (NM_001286235) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | SLC2A14 |
| Synonyms: | GLUT14; SLC2A3P3 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >SC336465 representing NM_001286235.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGACAACAGACAGAATGTCACCCAGCTCTGATCTTTGCCATCACAGTTGCTACAATCGGCTCTTTTC
CAGTTTGGCTACAACACTGGGGTCATCAATGCTCCTGAGACGATCATAAAGGAATTTATCAATAAACT
TTGACGGACAAGGCAAATGCCCTCCCTCTGAGGTGCTGCTCACGAATCTCTGGTCTTGTCTGTGGCC
ATATTTCCGTCGGGGTATGATCGGCTCTTTCCGTCGGACTCTTTGTTAACCCTTTGGCAGGCGC
AATTCAATGCTGATTGTCAACCTGTTGGCTGCCACTGGTGGTGCCTTATGGGACTGTGTAATAAGCT
GAGTCAGTTGAAATGCTGATCCTGGGCCGCTTGGTTATTGGCCTCTTCTGCGGACTCTGCACAGTTTTT
GTGCCATGTACATTGGAGAGATCTCGCCTACTGCCCTGAGGGGTGCCTTTGGCACTCTCAACCAGCTG
GGCATAGTTATTGGAATCTGGTGGCCAGATCTTTGGTCTGGAATCATCTTGGGTCTGAAGAGCTA
TGGCCGGTCTATTAGGCTTTACCATCCTTCCAGCTATCCTGCAAAGTGCAGCCCTTCCATGTTGCCCT
GAAAGTCCAGATTTTTGCTCATTAAACAGAAAAAAGAGGAGAATGCTACGCGGATCCTCCAGCGGTTG
TGGGGCACCAGGATGTATCCCAAGACATCCAGGAGATGAAAGATGAGAGTGCAAGGATGTCACAAGAA
AAGCAAGTCACCGTGTGGAGCTCTTAGAGTGTCCAGCTACCGACAGCCATCATCATTCCATTGTG
CTCCAGCTCTCTCAGCAGCTCTCTGGGATCAATGCTGTGTTCTATTACTCAACAGGAATCTTCAAGGAT
GCAGGTGTTCAACAGCCCATCTATGCCACCATCAGCGCGGGTGTGGTTAAACTACTTCTTCACTTTACTT
TCTCTATTTCTGGTGGAAAGGCGAGGAAGAAGGACTCTGCATATGATAGGCCTTGGAGGGATGGCTTTT
GGGGCTATCTGGTCTTTGTGGCCTGTTTTGAAATTGGACCAGGCCCATCCCTGGTTTATTGTGGCC
GAACCTTTCAGCCAGGGCCCCCGCCAGCTGCGATGGCAGTGGCCGGCTGCTCCAACGGACCTCCAAC
TTCTAGTCGGATTGCTCTTCCCCTCTGCTGCTTACTATTTAGGAGCCTACGTTTTTTATTATCTTACC
GGTCTCTCATTACCTTCTTGGCCTTACCTTCTCAAAGTCCCTGAGACCCGTGGCAGGACTTTTGAG
GATATCACACGGGCCCTTGAAGGGCAGGCACACGGTGCAGATAGATCTGGGAAGGACGGCCTCATGGGG
ATGAACAGCATCGAGCCTGCTAAGGAGACCACCACCAATGTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCG
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001286235

| | |
|-------------------------------|--|
| Insert Size: | 1494 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001286235.1</u> |
| RefSeq Size: | 3912 bp |
| RefSeq ORF: | 1494 bp |
| Locus ID: | 144195 |
| UniProt ID: | <u>Q8TDB8</u> |
| Cytogenetics: | 12p13.31 |
| Protein Families: | Transmembrane |
| MW: | 54 kDa |
| Gene Summary: | <p>Members of the glucose transporter (GLUT) family, including SLC2A14, are highly conserved integral membrane proteins that transport hexoses such as glucose and fructose into all mammalian cells. GLUTs show tissue and cell-type specific expression (Wu and Freeze, 2002 [PubMed 12504846]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate upstream start codon, compared to variant 1. The encoded isoform (b) has a distinct N-terminus and is shorter than isoform a. Variants 3 and 4 encode the same isoform (b).</p> |