

Product datasheet for **SC336470**

C14orf93 (NM_001282970) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C14orf93 (NM_001282970) Human Untagged Clone
Tag:	Tag Free
Symbol:	C14orf93
Synonyms:	RTFC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336470 representing NM_001282970.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCTTCAGTGCCACCATTCTCTTCTCCCCTCCCAGTGGCAGCGAGGCCAGATGCTGCTGCTGCGCC
TGTAAGAGTGAGACTAATGGAGGCAACACAGGCTCCCAGGGTGGGAATCCTCCTCCAGCACCCCATC
ACAGTGACTGGACATGGCTTGGCTGTTCCAGAGCTCAGAGCAGCTCCTGCATGTTATCTACCAGCGGGTC
GATAAGGCAGTGGGTTTGGCTGAAGCTGCTCTGGGCTTGGCCAGGGCCAACAATGAGTTGTTAAAACGT
CTCCAGGAGGAAGTGGGTGACCTGAGGCAAGGAAAGTGTCCATCCCTGATGAAGATGGGAAAGCCGG
GCACATAGTTCCCACCTGAGGAGCCTGGGCTCTCAAGGAAAGTCCCGGGAAGCCTTTAAGGCTCTG
TCTGCCGTGGAAGAGGAGTGTGACAGCGTGGGCAGCGCGTGCAGGTGGTATTGAGGAGCTGCGGCAG
CTGGGAGCAGCCTCAGTGGGCTGGGCTTTGGGCTTCCAGCAACTCAGAGGGACATGCGGCTCCCA
GGGTGCACGCTGGCTGCCAGCGAGGCGGCCCTTCTCAATCCTCTGGTGGATGATTACGTGGCCTCT
GAGGGTGCAGTACAGCGAGTTCTGGTCCCTGCTTATGCCAAGCAACTCTACCAGCCACACAACCTGGCA
ATCCAGCGGGCAACCCAGAGACAGGACCAGAAAATGGAACCAAGCTGCCACCACCCCGCCTGAGGAC
ATGCTCAATGCCCTGCTGCGCTGGACAGTGCCTTGAAGAGTCAAGCCCTGGGAGCACTGGGGAGCTG
AGACACTCTCTAGGGCTGACCGTTTCCCATGCAGGACCAGGAAAGTGGGCAGAAGAACTCCAGGGCG
AAGCGGATCTTGTACTCTCTAACTGGTCCACAATGTGCATAACCACATACCAATGACAAGAGATTC
AATGGGTCTGAAAGCATCAAGTCTCTTGAATATTTAGTAGTGAAGTTTCTTCTGAAAAGCTCAAG
CAAGAGCTGGTGACCAGTCCCACAATTACACTGATAAGGAGCTAAAAGGAGCCTGTGTGGCCTACTTC
CTTACTAAGAGGCGTGAGTACCGCAACTCCCTGAACCCCTTTAAGGCCTGAAGGAAAAAGAGGAGAAG
AAACTTCGAAGTCGCCGATATCGCTTTTGGCAACCGATCCAGTATCATGAGGCATTTTGGACCTGAG
GACCAACGCTCTGTGGAATGATGTGACAGAGAACTGATGTCAGATGAAGAGGACAGTCTTAAACGAGCCA
GACAGACTGCCTTCTGCTGAAGCCAGCTCCTTCCACCAGAACTTTACAATCCTAATTTCCAAGAAGAG
GAAGATGAGGGAGGGATGAGAATGCACCTGGCTCCCATCTTTTACCAACCCACAAAACCTGCTGT
CCTGACTTGAACCTATTCATTGAAATCAAGGTGAAAAGGATGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282970

Insert Size:	1497 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001282970.1</u>
RefSeq Size:	2310 bp
RefSeq ORF:	1497 bp
Locus ID:	60686
UniProt ID:	<u>Q9H972</u>
Cytogenetics:	14q11.2
Protein Families:	Secreted Protein
MW:	54.2 kDa