

## Product datasheet for **SC336485**

### ANKRD42 (NM\_001300972) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD42 (NM_001300972) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD42
Synonyms:	PPP1R79; SARP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC336485 representing NM\_001300972.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCCGGGGTGGCCAATTCAGGCCCTCCACTTCCTCTAGGGAGACTGCAAACCCCTGTTCCAGGAAG
AAGGTGCATTTTGGCAGCATACATGATGCAGTACGAGCTGGAGATGTAAGCAGCTTTCAGAAATAGTG
GTACGTGGAGCCAGCATTAAATGAACTTGATGTTCTCCATAAGTTTACCCCTTTACATTGGGCAGCACAT
TCTGGAAGTTTGGAGTGTCTTATTGGCTGCTCTGGCATGGAGCTGATATCACACACGTAACAACGAGA
GGTTGGACAGCATCTCACATAGCTGCAATCAGGGGTGAGGATGCTTGTGTACAGGCTCTTATAATGAAT
GGAGCAAATCTGACAGCCAGGATGACCGGGATGCACTCCTTTACATCTTGCTGCAACTCATGGACAT
TCTTTCACTTTACAAATAATGCTCCGAAGTGGAGTGGATCCCAGTGTGACTGATAAGAGAGAATGGAGA
CCTGTGCATTATGCAGCTTTTCATGGCGGCTTGGCTGCTTGAACCTTCTGTAAATGGGGTTGTAGC
ATAGAAGATGTGGACTACAATGAAACCTTCCAGTTCCTTAGCAGCCATGGAAGGCCACCTTCACTGT
TTCAAATTCCTAGTCAGTAGAATGAGCAGTGCAGCAAGTTTTAAAAGCTTTCAATGATAATGGAGAA
AATGTACTGGATTTGGCCAGAGGTTCTCAAGCAGAACATTTTACAGTTTATCCAGGGGCTGAGTAT
GAAGGAAAAGACCTAGAGGATCAGGAACTTTAGCATTTCAGGTCATGTGGCTGCCTTTAAGGGTGAT
TTGGGGATGCTTAAGAAATAGTGGAAGATGGAGTAATCAATATTAATGAGCGTGCTGATAATGGATCA
ACTCCTATGCATAAAGCTGCTGGACAAGGCCACATAGAGTGTTTGCAGTGGTTAATTAATGGGAGCA
GACAGTAATATTACCAACAAAGCAGGGGAGAGACCCAGTATGTGGCAAAGAGGTTTGGCCATTTGGCA
GCAGTGAAGCTGTTAGAGGAGCTACAGAAATATGATATAGATGACGAAAATGAAATTGATGAAAATGAT
GTGAAATATTTTATAAGACATGGTGTGGAGGAAGCACTGATGCCAAGGATGATTTATGTCTGAGTGAC
TTGGATAAAACAGATGCCAGAATGAGAGCTTACAAGAAAATTTAGAATTGAGACACCTCCTGAAAT
GCCGAGAGCAACTATAAACACTTGGGAGGCATAACAGAAGAAGATTTAAAGCAGAAGAAAAGAACAGCTT
GAGTCTGAAAAGACCATCAAAGAACTGCAGGGCCAGCTGGAGTATGAACGACTACGTAGAGAAAAATTA
GAATGTCAGCTTGATGAATATCGAGCAGAAGTTGATCAACTCAGGGAAACACTGGAAAAATTCAGTC
CCAAACTTTGTGGCTATGACAGCGCTTCTTGTGAGTCAAACAAAGAGAAGAGGGCGAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001300972

<b>Insert Size:</b>	1509 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001300972.1</a></u>
<b>RefSeq Size:</b>	2662 bp
<b>RefSeq ORF:</b>	1509 bp
<b>Locus ID:</b>	338699
<b>UniProt ID:</b>	<u><a href="#">Q8N9B4</a></u>
<b>Cytogenetics:</b>	11q14.1
<b>MW:</b>	56.2 kDa