

## Product datasheet for **SC336497**

### ZNF596 (NM\_001287256) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF596 (NM_001287256) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF596
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC336497 representing NM_001287256. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCCATCACCGGATTCCATGACCTTCGAGGATATCATTGTAGACTTCACTCAAGAAGAGTGGGCCCTG
CTGGACACATCCCAGAGAAAAGCTGTTTCAAGATGTGATGTTGGAGAACATCAGTCATCTGGTCTCTATT
GGCAAACAGCTCTGCAAATCAGTTGTGCTTCCCAATTGGAGCAAGTAGAGAACTTTCAACACAAAAGA
ATAAGCTTACTGCAAGGTAGAGAAGTTGGCATTAAACATCAAGAGATACCATTATTCAACATATCTAT
CAGAAGGGCAGTCCACCATCAGCACAATGAGATCTCATACTCAAGAGGATCCTTTTCTATGCAATGAC
TTAGGAGAAGATTTCACTCAACATATAGCATTGACTCAAAATGTGATTACCTACATGAGAACGAAACAC
TTTGTAAGCAAAAAGTTTGGGAAAATCTTCAGTGACTGGTTATCCTTTAATCAACACAAGGAAATTCAC
ACCAAATGTAAATCATATGGAAGTCATCTATTTGATTATGCCTTTATCCAAAACCTGCCCCTTAGACCA
CACAGTGTGACTCACACTAGAGAGATAACATTGGAATGTCGTGTGTGGGAAAACCTTTAGCAAAAAT
TCTAATCTTAGGCGACATGAGATGATTCACACTGGAGAGAAACCACCGGATGTCATCTATGTGGGAAA
GCCTTTACTCATTGCTCTGATCTTCGAAAACATGAGAGAACTCACACTGGAGAGAAGCCATATGGATGT
CATCTATGTGGGAAAGCCTTCAGTAAAAGTTCTAACCTTAGACGACATGAGATGATTCACACTAGAGAA
AAAGCACAGATATGCCATCTATGTGGGAAAGCCTTCACCTATTGCTCTGACCTTAGAAAACATGAGAGA
ACTCACTTAGGAGATAAACCATATGGATGTCTCCTATGTGGGAAGGCTTTCAGTAAATGTTCTTACCTT
AGACAACATGAAAGAACCAATGGAGAGAAACCATATGAATGTCATCTATGTGGAAAAGCCTTCTCT
CATTGTTCTCACCTTAGACAACATGAGCGAAGTCAATGGAGAGAAACCACATGGATGTCATCTATGT
GGGAAAGCATTCACTGAATCTTCTGTGCTTAAACGACATGAGAGAATTCACACTGGAGAGAAACCATAT
GAGTGCCATGTATGTGGGAAAGCCTTCACTGAATCTTCTGACCTCAGACGACATGAGAGAACTCACACT
GGAGAAAACCATATGAATGCCATCTATGCGGAAAAGCCTTCAATCACTCTTCTGTCTTAGACGACAT
GAGAGAACTCACACTGGAGAGAAACCATATGAATGCAATATATGTGGTAAAGCCTTCAATAGAAGTTAC
AACTTTAGACTTCATAGAAGAGTTCACTGGAGAGAAACCATATGTATGTCCTCTATGTGGGAAAGCC
TTTAGTAAATTTTTAACCTTAGACAACATGAGAGAACTCACACTAAAAAAGCAATGAATATGTAA
```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001287256

**Insert Size:** 1515 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001287256.1](#)

**RefSeq Size:** 2692 bp

**RefSeq ORF:** 1515 bp

**Locus ID:** 169270

**UniProt ID:** [Q8TC21](#)

**Cytogenetics:** 8p23.3

**Protein Families:** Druggable Genome, Transcription Factors

**MW:** 58.5 kDa

**Gene Summary:** May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (6) differs in the 5' UTR compared to variant 1. Variants 1 through 6 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.