

Product datasheet for SC336500

Annexin A11 (ANXA11) (NM_001278408) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Annexin A11 (ANXA11) (NM_001278408) Human Untagged Clone
Tag: Tag Free
Symbol: ANXA11
Synonyms: ALS23; ANX11; CAP-50; CAP50
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC336500 representing NM_001278408.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGAGCTACCCTGGCTATCCCCGCCCCAGGTGGCTACCCACCAGCTGCACCAGGTGGTGGTCCCTGG
GGAGGTGCTGCCTACCCTCCTCCGCCAGCATGCCCCCATCGGGCTGGATAACGTGGCCACCTATGCG
GGGCAGTTCACCCAGGACTATCTCTCGGGAATGGCGCCAACATGTCTGGGACATTTGGAGGAGCCAAC
ATGCCAACCTGTACCCTGGGGCCCTGGGGCTGGCTACCCACCAGTGCCCCCTGGCGGCTTTGGGCAG
CCCCCTCTGCCAGCAGCCTGTTCTCCCTATGGGATGTATCCACCCAGGAGAAACCCACCCTCC
AGGATGCCCTCATATCCGCCATACCCAGGGCCCTGTGCGGGCCAGCCATGCCACCCCGGACAG
CAGCCCCAGGGGCTACCCTGGGACGCCACCAGTGACCTACCCTGGTCAGCCTCCAGTGCCACTCCCT
GGGCAGCAGCAGCAGTCCGAGCTACCCAGGATACCCGGGGTCTGGGACTGTCACCCCGCTGTGCC
CCAACCCAGTTTGAAGCCGAGGCACCATCACTGATGCTCCCGGCTTTGACCCCTGCGAGATGCCGAG
GTCCTGCGGAAGGCCATGAAAGGCTTCGGGACGGATGAGCAGGCCATCATTGACTGCCTGGGGAGTCGC
TCCAACAAGCAGCGGCAGCAGATCCTACTTTCTTCAAGACGGCTTACGGCAAGGATTTGATCAAAGAT
CTGAAATCTGAACTGTCAGGAACTTTGAGAAGACAATCTTGGCTCTGATGAAGACCCAGTCTCTTT
GACATTTATGAGATAAAGGAAGCCATCAAGGGGTTGGCACTGATGAAGCCTGCCTGATTGAGATCCTC
GTTCCCGCAGCAATGAGCACATCCGAGAATTAACAGAGCCTACAAAGCAGAATTCAAAAGACCCTG
GAAGAGGCCATTGAAAGCAGACATCAGGGCACTTCCAGCGGCTCCTCATCTCTCTCAGGGAAC
CGTGATGAAAGCACAACGTGGACATGTCACTCGCCAGAGAGATGCCAGGAGCTGTATGCGGCCGGG
GAGAACCCTGGGAACAGACGAGTCCAAGTTCAATGCGGTTCTGTGCTCCCGGAGCCGGGCCACCTG
GTAGCAGTTTTCAATGAGTACCAGAGAATGACAGGCCGGGACATTGAGAAGAGCATCTGCCGGGAGATG
TCCGGGGACCTGGAGGAGGCATGCTGGCCGTGGTAAATGTCTCAAGAATACCCAGCCTTCTTTGCG
GAGAGGCTCAACAAGCCATGAGGGGGCAGGAACAAGGACCGACCCTGATTCGCATCATGGTGTCT
CGCAGCGAGACCGACCTCCTGGACATCAGATCAGAGTATAAGCGGATGTACGGCAAGTCGCTGTACCAC
GACATCTCGGGAGATACTTCAGGGGATTACCGGAAGATTCTGCTGAAGATCTGTGGTGGCAATGACTGA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001278408

Insert Size: 1518 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278408.1](#)

RefSeq Size: 3073 bp

RefSeq ORF: 1518 bp

Locus ID: 311

UniProt ID: [P50995](#)

Cytogenetics: 10q22.3

MW: 54.4 kDa

Gene Summary: This gene encodes a member of the annexin family, a group of calcium-dependent phospholipid-binding proteins. Annexins have unique N-terminal domains and conserved C-terminal domains, which contain calcium-dependent phospholipid-binding sites. The encoded protein is a 56-kD antigen recognized by sera from patients with various autoimmune diseases. Several transcript variants encoding two different isoforms have been identified. [provided by RefSeq, Dec 2015]
Transcript Variant: This variant (e) differs in the 5' UTR compared to variant a. Variants a, b, c, d, and e all encode the same isoform (1).