

Product datasheet for **SC336521**

FSD1L (NM_001287191) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSD1L (NM_001287191) Human Untagged Clone
Tag:	Tag Free
Symbol:	FSD1L
Synonyms:	CCDC10; CSDUFD1; FSD1CL; FSD1NL; MIR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336521 representing NM_001287191.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACTCCAGAAAGAAGCTCTACAGAGGATCATTTCAACTCTGGCAAATAAAAATGATGAAATTCAG
AACTTTATTGATACACTACATCATACTAAAAGGAGTTCAGAAAATTCGTCCAACATACTCTCAGAG
TTAGATGAAGAATTTGATAGTTTATACTCTATACTGGATGAAGTAAAAGAAAGTATGATTAAGTATC
AAGCAGGAACAAGCTCGTAAATCCCAAGAGTTACAGAGTCAGATTAGTCAATGTAATAATGCCCTGGAG
AACTCTGAAGAACTATTAGAATTTGCAACAAGGTCATTAGATATAAAGGAACCTGAAGAATTTCAAAG
GCTGCCAGACAGATCAAGGATAGAGTCACAATGGCTTCAGCCTTTCGCCTTTCCTTGAACCAAAGGTC
AGTGACAACATGACTCATTTAATGGTGGATTTCTCACAGGAAAGACAGATGCTGCAAACCTTGAAGTTT
TTGCCAGTCCCCAAAGCTCCAGAGATAGATCCAGTAGAGTGTGGTGGCAGATAACTCTGTAACAGTG
GCTTGGAGAATGCCAGAAGAAGATAATAAGATTGACATTTTATACTGGAACATAGGAAGACTAATTTT
GATGGACTTCCAGTGTAAGGATGAGCGATGCTGGGAGATAATTGATAATTAAGGGTACTGAATAT
ACACTATCAGGCTTAAAATTTGATTCAAAGTATATGAATTTAGAGTGGCAGCTTGTAAACAAGGCTGTG
GCTGGAGAGTATTCTGATCCAGTGACTCTAGAGACCAAAGCACTTAACTCAATTTGGATAACTCCTCA
TCCCATTTGAACCTGAAAGTTGAAGATACATGTGTAGAGTGGGATCCTACTGGAGAAAAGGTCAGAA
AGTAAAATTAAGGAAAAGAGAACAAGGCAGTGTCCATGTTACTTCACTGAAGAAAACATACAAGTGGT
ACACCATCCCCAAAACGAACATCTGTAGGCTCCAGGCCACCAGCAGTAAGAGGCAGTAGAGATCGTTTT
ACTGGAGAATCATACACAGTGTGGGAGACACTGCTATTGAAAGTGGACAACATTATTGGGAGGTCAAG
GCCCAGAAGGATTGTAATCCTACAGTGTGGGAGTAGCATACAAAACGTTGGGGAATTTGACCAATTG
GGAAAGACAAACACTAGTTGGTGTATCCATGTCACAACTGGCTACAAAACACATTTGCAGCAAAGCAT
AATAATAAAGTCAAAGCTTTGGATGTTACTGTTCTGAAAAAATAGGTGATTTTGTGATTTGATGGG
GGTCAACTTTCTATGATGCAAATCTAAACAGTTGCTATATTCTTTAAGACAAAATTTACTCAG
CCAGTACTACCTGGTTTCATGGTATGGTGTGGTGGACTTCTTTGAGTACTGGGATGCAGGTTCCAAGT
GCTGTGAGAACACTTCAGAAAAGTAAAATGGAATGACTGGTTCAGTACGAGCCTGAACAATGTTGTT
ACTCAATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001287191

Insert Size:	1527 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001287191.1</u>
RefSeq Size:	7684 bp
RefSeq ORF:	1527 bp
Locus ID:	83856
UniProt ID:	<u>Q9BXM9</u>
Cytogenetics:	9q31.2
MW:	57 kDa