

Product datasheet for **SC336523**

ZNF343 (NM_001282498) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF343 (NM_001282498) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF343
Synonyms:	dj734P14.5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336523 representing NM_001282498.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCTGGAGAATTACAGGAATCTTCTCATTGGCAGAACCAAGCCAGAAATCTACACTTGTTCCTCC
TGCCTTCTGGCCTTCTCCTGTCCAGGATTCCTCAGTCAACATGTACTTCAGATCTTCTGGGCTTATGT
GCAGAAAATCACTTCCATCCAGGAATTCTAGCCAGGGCATTGGAACAGCAGGGGCAGCAGTATTCC
CATGTAAGCTGTTGGTTTAAAAATGCAGAAGGTCAAGGAGAGAGGAGGTGGCTCCAAAACCTGGTCTGCA
AGGACAGAGGAGAGAGAAACCTCAAGGGCATTCCCAGCCACTCCAAAGACAGTCAGCAAGTCCTAGA
AAAGGCAACATGGTGGTAGAAACAGAGCCAGCTCAGCCCAAAGACCAAACCTGTGCAGCTAGACAAA
GGCTTGAAGGAATTAGAAACCTTGAGATTTGGAGCAATCAACTGTAGAGAGTATGAACCGGACCATAAC
CTGGAATCAAACCTTATTACAAACCCGAGGACCCTTTAGGGAAGAAGCCCTACATTTGCAGTGATTGT
GGGCGAAGCTTTAAAGATAGATCAACCCTCATCAGACACCATCGTATACACTCGATGGAGAAGCCTTAT
GTGTGCAGTGAGTGGGCGAGGTTTTAGCCAGAAGTCCAACCTCAGCAGACACCAGAGAACACATTCA
GAAGAGAAGCCTTATTTGTGCAGGGAGTGTGGCAAAAGCTTTAGAAGTAAGTCCATCCTCAATAGACAT
CAGTGGACTCACTCAGAGGAGAAGCCCTATGTTTGCAGCGAGTGTGGGCGAGGCTTTAGCGAGAAGTCA
TCCTTCATCAGACACCAGAGGACACACTCCGGTGAGAAACCTATGTGTGCCTGGAGTGTGGACGAAGC
TTTTGTGATAAGTCAACCCTCAGAAAACACCAGAGGATACACTCAGGGGAGAAGCCTTATGTTTGCAGG
GAGTGTGGGCGAGGCTTTAGCCAGAACTCAGATCTCATCAAACACCAGAGGACACACTTGGATGAGAAG
CCTTATGTTTGCAGGGAGTGTGGGCGAGGCTTTGTGACAAGTCAACCCTCATCATAACAGAGCGGACG
CACTCTGGAGAGAAGCCTTATGTGTGGTGTGAGTGTGGCCGAGGCTTTAGTCGGAATCACTCCTCCTT
GTCCACCAGAGGACACACTCAGGGGAGAAGCATTATGTCTGCAGGAGTGTAGGCGAGGTTTTAGCCAG
AAGTCAAATCTCATCAGACACCAGAGGACGCACTCAAATGAGAAGCCTTATATTTGCAGGAATGTGGG
CGAGGCTTTTGTGACAAGTCAACCCTATTGTACATGAGAGGACACACTCAGGAGAGAAGCCTTACGTG
TGCAGTGAGTGTGGCCGAGGCTTTAGCCGAAATCACTCCTCCTTGTCCACCAGAGGACACACTCAGGG
GAGAAGCATTATGTTTGTAGGGAGTGTGGGCGAGGCTTTAGTCATAAGTCAAATCTCATCAGACACCAG
AGGACACACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001282498

Insert Size:	1530 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001282498.1</u>
RefSeq Size:	3318 bp
RefSeq ORF:	1530 bp
Locus ID:	79175
UniProt ID:	<u>Q6P1L6</u>
Cytogenetics:	20p13
Protein Families:	Transcription Factors
MW:	58.8 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (5) differs in its 5' UTR and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (3) is shorter, compared to isoform 1.