

## Product datasheet for **SC336530**

### ZNF92 (NM\_001287533) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF92 (NM_001287533) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF92
Synonyms:	HEL-203; HPF12; HTF12; TF12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC336530 representing NM\_001287533.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTGTTCTCATTTTGCCTAAGATGTTTGGCCAGAGCACAGCATAAAAGATTCTTTCCAAAAAGTGATA
CTGAGAACATATGGAAAATATGGACATGAGAATTTACAGCTAAGAAAAGACCATAAAAAGTGGATGCA
TGTAAAGGTGTACAAAGGAGTTATAATGGACTTAACAGTGTGGTGAACAATACTGACAGCAAGATATTT
CAGTGTGATAAATATGTGAAAGTCTTTTCATAAATTTCCAATGTAATAGAAAATAAGATAAGACATACT
GGAAAAGAACCTTTCAAATGTAAAAACCGTGGCAAATCATTTTGCATGCTTTACAATTAACCAACAT
AAGAAAATTCATACTAGAGAGTATTCTTACAAATGTGAAGAATGTGGTAAAGCCTTTAAGTGGTCTCA
ACCCTTACTAACATAAGATAATTCATACTGGAGAAAAACCTACAAATGTGAAGAATGTGGCAAAGCT
TTTAAACCGGTCTCAAATCTTACTAACATAAAAATAATTCATACTGGAGAGAAACCTACAAATGTGA
GAATGTGGCAAAGCTTTTAAACCGGTCTCAACCTTACTAAACATAAAAAGAATTCATACAGAAGAGAAA
CCCTACAAATGTGAAGAATGTGGCAAGGCCTTTAACCAGTCTCGATTCTAATAACATAAGAGAATT
CATATGGAAGATAAACCCCTACAAATGTGAAGAATGTGGCAAGGCCTTTAGAGTATTCTCAATTCTAAA
AAACATAAGATAATCCATACTGGGGAAAAACCATACAAATGTGAAGAATGTGGCAAGGCCTTTAACCAG
TTCTCAAACCTTACTAACATAAGATAATTCATACTGGAGAGAAACCTACAAATGTGATGAATGTGGC
AAAGCCTTTAACCAGTCTCAACCTTACTAAACATAAAAAGAATTCATACGGGAGAAAAACCTACAAA
TGTGAAGAATGTGGCAAAGCTTTTAAACAGTCTCAACCTTACTGAACATAAGATAATTCATACTGGA
GAGAAACCTACAAATGTGAAAATGTGGCAAGGCCTTTAGCTGGTCTCAGCTTTTACTAAACATAAG
AGAAATCATATGGAAGATAAACCCCTACAAATGTGAAGAATGTGGCAAAGCCTTTAGTGTATTCTCAAC
CTTACTAACATAAAAATAATTCATACTAGAGAAAAACCTACAAATGTGAAGAATGTGGCAAAGCCTTT
AACAGTCTCAATTTTACTAACATAAGATAATTCACACTGAAGGGAATCCTACAAATGTGAAAAA
TGTGGCAATGCTTTTAAACAGTCTCAACCTTACTGCACGTAAGATAATTTATACTGGAGAGAAACCC
TACAAATATGAAGAATGTGACAAAGCCTTTAACAAGTTCTCAACCTTATTACACATCAGATAATTTAT
ACTGGAGAGAAACCTGCAAACATGAATGTGGCAGAGCCTTTAACAAATCCTCAAATTATACTAAAGAG
AAACTACAAACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001287533

<b>Insert Size:</b>	1533 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001287533.1</a></u>
<b>RefSeq Size:</b>	3342 bp
<b>RefSeq ORF:</b>	1533 bp
<b>Locus ID:</b>	168374
<b>UniProt ID:</b>	<u><a href="#">Q03936</a></u>
<b>Cytogenetics:</b>	7q11.21
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	59.4 kDa
<b>Gene Summary:</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR, lacks a portion of the coding region and initiates translation at a downstream start codon, compared to variant 2. Variants 4 and 5 encode the same isoform (4), which is shorter at the N-terminus, compared to isoform 2.