

Product datasheet for **SC336539**

ADA2 (NM_001282225) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADA2 (NM_001282225) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADA2
Synonyms:	ADGF; CECR1; IDGFL; PAN; SNEDS; VAIHS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336539 representing NM_001282225.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTGGTGGATGGCCATCTGAGCGGCCAGCCCTGTGCTTCTTGCTGTTGGCTGTGGCAATGTCTTTC
TTCGGCTCAGCTCTATCCATAGATGAAACACGGGCGCATCTGTTGTTGAAAGAAAAGATGATGCGGCTG
GGGGGGCGGCTGGTCTGAACACCAAGGAGGAGCTGGCCAATGAGAGGCTCATGACGCTCAAATCGCT
GAGATGAAGGAGGCCATGAGGACCTGATATCCACCCAGCATGCACTTTTTCCAGGCCAAGCATCTC
ATTGAGAGAAGTCAAGTGTAAATATTCTAAGGATGATGCCAAAAGGGGCTGCCTTGACCTCCATGAC
ATTGGCATCGTGACTATGGACTGGCTGGTGAAGATGTCACCTACAGGCCTCACTGCCACATCTGTTTC
ACCCAAAGGGGGATCATGCAGTTCAGATTTGCTCACCCAACTCCCGTCCATCAGAAAAATGTTCCAAG
TGGATTCTGCTGGAGGATTATCGGAAGCGGGTGCAGAACGTCCTGAGTTTGTGACAGCTTGTGAGG
AATTTCACTCTGGTGACCCAGCACCCGGAGGTGATTTACAAAACCAAAATGTTGTCTGGTGGAAATTT
GAAACCATCTTCTCACCATCTCTGGTCTCATCCATTACGCACCAGTGTTCAGAGACTATGTCTCCGG
AGCATGCAGGAGTTCTACGAGGACAACGTGCTCTACATGGAGATCAGAGCCAGGCTGCTGCCGGTGTAT
GAGCTCAGTGGAGAGCACCATGACGAAGAGTGGTCACTGAAGACTTACCAGGAAGTAGCTCAGAAGTTT
GTGAAAACACCCCTGAGTTTATTGGAATCAAAATCATTATTCCGGATCAGAGTCCAAAGATGTGGCT
GTCATCGCAGAATCCATCCGAATGGCCATGGGGCTCCGAATCAAGTTCCACCGTGGTGGCAGGGTTT
GACCTGGTGGGGCATGAGGACACTGGCCACTCCTTGATGACTACAAGGAAGCTCTGATGATCCCGCC
AAGGATGGCGTTAAGCTGCCTTACTTCTCCACGCCGGAGAAACAGACTGGCAGGGTACTTCCATAGAC
AGGAACATTCTGGATGCTCTGATGCTGAACACTACCAGAATCGGCCATGGATTTGCTTTGAGCAAACAC
CCCGCAGTCAGGACTTACTCTGAAAAAAGGACATCCCATAGAAGTCTGTCCATCTCAACAGGTG
CTGAAACTGGTGTCTGACTTGGGAACCCCTGTAGCCACTCTGATGGCCACTGGGCACCCCATGGTG
ATCAGCTCTGATGACCCAGCTATGTTTGGTGCCAAAGGCTTGTCTATGATTTCTATGAGGTCTTCATG
GGCATTGGGGGATGAAGGCTGACCTGAGGACCCTCAAACAGCTGGCCATGAACTCTATCAAGTACAGT
ACCCTGTTGGAGAGTGAGAAAAATACTTTCATGGAAATCTGGAAGAAGAGATGGGATAAGTTCATAGCA
GATGTGGCTACAAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282225

Insert Size:	1536 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282225.1
RefSeq Size:	4510 bp
RefSeq ORF:	1536 bp
Locus ID:	51816
UniProt ID:	Q9NZK5
Cytogenetics:	22q11.1
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	58.9 kDa
Gene Summary:	<p>This gene encodes a member of a subfamily of the adenosine deaminase protein family. The encoded protein is one of two adenosine deaminases found in humans, which regulate levels of the signaling molecule, adenosine. The encoded protein is secreted from monocytes undergoing differentiation and may regulate cell proliferation and differentiation. This gene may be responsible for some of the phenotypic features associated with cat eye syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]</p> <p>Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (a). Both variants 3 and 4 encode the same isoform.</p>