

Product datasheet for **SC336541**

ZNF286 (ZNF286A) (NM_001288643) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF286 (ZNF286A) (NM_001288643) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF286A
Synonyms:	ZNF286
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336541 representing NM_001288643.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAGCTCTGTCTTCCAGGATTCTCCCATTTCCAAGAGAAGAGCACAGAAGAGGGAGAAGTGGCT
GCTCTGCCCTCACGGCAGATCCAGGAAACAGTGACATTCAAGGATGTGGCCATGGACTTTACACCA
GAGGAGTGGGGGAAGCTGGATCCTGCACAAAGGATGTGATGCTGGAGAACTATAGGAACCTAGTCTCA
CTTTGGCTTCCAGTTTCAAACCTGAGAGCTACAACCTGGAGAATGAAAAAGAACCATTGAAGCTTGAG
AGAAAAGCCCCAAAAGCAGCTATTAGACATGGAGACTAGACCACAGAGCAAGGATTCAACTTCAGTG
CAAGATTTTTCAAAGCAGAATCATGCAAAGTTGCAATAATAGACAGACTGACACGGAATAGTGTCTAT
GACTCTAACTTGAAGCAGCCCTTGAATGTGAAAATTGGTTAGAGAATCAGCAAGGAAATCAGGAGAGA
CATTGAGAGAAATGTTCACTCACATGAATCACTCTCTGAGGAAACAGACCATAAGCATGATGTATAC
TGGAAAAGCTTCAATCAGAAATCTGTCCTTACTCTGAAGACAGAGTTCCCAAAGGATCTTATGCCTTC
CATACACTTGAAAAAGCTTGAAAACAAAATCAAACCTTAATGAAAAGCAGAGAACTTATAAAGAGAAA
AAACCTCATAAATGTAATGATTGTGGTGAACCTTTCACCTACCATTCAAGTCTTATTCGACACCAGAGA
GTCCATACTGGAGAGAAAACCTATACCTGCAATGAATGTGGGAAATCTTTTAGCCACAGAGCTAATTTA
ACTAAACACCAGAGAACTCATACTAGAATTCTCTTTGAATGCAGTGAATGCAAGAAAACCTTCACAGAA
AGCTCATCCCTTGCAACACATCAGAGAATTCACGTTGGAGAGAGACCTTATGAATGCAATGAATGTGGG
AAAGGTTTTAATCGAAGTACACATCTTGTGCAGCATCAGTTGATTCACTGGAGTGAAGCCTTATGAA
TGTAATGAATGTGATAAAGCTTTTATTCATTATCAGCACTCATTAAACATCAAAGAACTCATACTGGA
GAGAAACCTTATAAATGTCAAGAATGTGGGAAAGCCTTTAGCCATTGCTCATCCCTAACTAAGCATCAG
AGAGTTCATACTGGAGAAAAGCCATATGAATGCAGTGAATGTGGAAAAACTTTTAGTCAGAGCACACAT
CTTGTTCAACATCAGAGAATTCACACTGGAGAGAAAACCTATGAGTGAATGAATGTGGGAAAACCTTC
AGCCGGAGCTCCAATTTTGTAAACATCAAAGAATTCATATTGGAAGAAAACCGTACAAATGTAGCGAG
TGTGGAAAAGCCTTATTATTATCAGCTCTCATTCAACATCAGAGAACTCATACCGGAGAGAAAACCC
TTTAGATGTAATGAGTGTGGGAAAAGCTTTAAGTGCAGTTCATCTCTCATCAGACATCAAAGAGTTCAC
ACTGAAGAGCAACCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001288643

Insert Size:	1536 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001288643.1</u>
RefSeq Size:	5240 bp
RefSeq ORF:	1536 bp
Locus ID:	57335
UniProt ID:	<u>Q9HBT8</u>
Cytogenetics:	17p12
Protein Families:	Transcription Factors
MW:	59 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) lacks an alternate exon in the 5' coding region, compared to variant 3. It encodes isoform 3, which is shorter than isoform 1.