

Product datasheet for **SC336552**

APPBP2 (NM_001282476) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	APPBP2 (NM_001282476) Human Untagged Clone
Tag:	Tag Free
Symbol:	APPBP2
Synonyms:	APP-BP2; HS.84084; PAT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336552 representing NM_001282476.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTACTACAAGACATTTGCTTCATCATTGTTTTTCAGGCTTTGATGGATCATGGTGTAAAGTTGCT
TCAGCTTTGGCCACTCATTTCAGTAGGGCGTCTTATATAGCAGAATCAGATGCTGCAGTAAAGGAA
AAAGCCATTCAGTTGGCTTTGTTTTAGTGGCTTTCTTTCAGATGCAGGCTGGTACAGTATGCTGAG
AAAGTTTTTCTGTCTGCCTTCAGTTGTGTACTCTACACGATGAGATGCTTCATTGGTTTCGTGCAGTA
GAATGTTGTGTGAGGTTGCTTCATGTGCGAAATGGAAACTGCAAATATCATTGGGTGAAGAAACATTT
AAATTAGCTCAGACATATATGGATAAACTATCAAAACATGGCCAGCAAGCAAATAAAGCTGCACTCTAT
GGAGAACTGTGTGCACTCCTATTTGCAAAAAGTCACTATGATGAGGCATACAAATGGTGCATCGAGGCA
ATGAAAGAAATTACAGCAGGCTTACCAGTGAAAGTTGGTGGATGTCTTAAGACAAGCTTCTAAGGCT
TGTGTAGTAAACGTGAATTTAAGAAGGCAGAACAGTTAATTAACATGCAGTGTATTTGGCACGGGAT
CATTTTGGATCCAAACACCCAAAATATTCTGATACACTGCTAGATTATGGTTCTACTTACTCAATGTA
GATAATATCTGTCACTGTTGCAATTTATCAGGCAGCCCTTGACATTAGACAGTCACTGTTTGGTGGC
AAAAATATCCACGTAGCAACAGCTCATGAAGATTTGGCCTACTCTTCTTATGTCCACCAGTATAGCTCT
GGGAAATTTGACAATGCACTATTTTCATGCAGAAAGAGCTATTGGTATCATTACCCACATCCTACCTGAA
GATCATCTTCTTTGGCTTCTTCAAAGAGGGTGAAGCACTTATTTTAGAGGAGATTGCAATTGATTGT
CATAATAAGGAACTGAACAGAGGCTGCTTCAAGAAGCTCATGATTTGCACCTGTCTTCACTCCAATA
GCTAAAAAGCTTTTGGGAATTTAATGTACAGACTGCAAAACTATGGAAACCTTGAAGACTTTAT
CAGTCAATGAGAAAATTAAGGAAGCTGAAGAAATGCACATCAAAGCAATTCAGATTAAGAACAACCTT
CTTGGTCAAGAAGATTATGAAGTAGCCCTTTCAGTGGGACATCTGGCTTCTTTATATAATTATGACATG
AATCAGTATGAAAATGCTGAGAACTTTATTTGCGATCTATAGCAATTGGGAAGAAAATTTTGGTGGAG
GGCTACAGTGGACTAGAATATGATTATCGAGGCTCATTAACTTTACAACCTTGGAAATTACGAG
AAAGTGTTTGAATATCACAATGTTCTGTCTAACTGGAACCGGTTGCGAGATCGGCAATATTCAGTGACA
GATGCTCTTGAAGATGTCAGCACCAGCCCCAGTCCACTGAAGAAGTGGTGCAGTCTTCTGATTTCT
CAGAAATGTCGAGGGACCGAGCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282476

Insert Size:	1545 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282476.1
RefSeq Size:	6410 bp
RefSeq ORF:	1545 bp
Locus ID:	10513
UniProt ID:	Q92624
Cytogenetics:	17q23.2
Protein Families:	Druggable Genome
MW:	58.5 kDa
Gene Summary:	<p>The protein encoded by this gene interacts with microtubules and is functionally associated with beta-amyloid precursor protein transport and/or processing. The beta-amyloid precursor protein is a cell surface protein with signal-transducing properties, and it is thought to play a role in the pathogenesis of Alzheimer's disease. The encoded protein may be involved in regulating cell death. This gene has been found to be highly expressed in breast cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (2) lacks an in-frame exon in the 5' coding region, compared to variant (1). The encoded isoform (2) is shorter, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>