

Product datasheet for **SC336561**

PXK (NM_001289098) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PXK (NM_001289098) Human Untagged Clone
Tag:	Tag Free
Symbol:	PXK
Synonyms:	MONAKA; SLOB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336561 representing NM_001289098.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTTCATGGAGAAGCCGCCAGCCGCAAGGTGCTGCTGGACGACACGGTGCCGCTGACAGCAGCC
ATCGAGGCGAGCCAGACCTGCAGTCCCACACGGAATATATTATTCGAGTGCAAAGAGGAATTTCTGTG
GAAAACAGCTGGCAGATTGTTAGAAGATACAGTGACTTTGATTTGCTTAACAACAGCTTACAGATTGCA
GGCCTAAGTCTACCTCTTCTCCCAAAAATTGATTGGTAACATGGATCGTGAATTCATAGCTGAAAGG
CAGAAAGGTCTTCAGAACTATCTCAACGTGATCACAACAATCATATCTTGTCTAATTGTGAGCTGGTT
AAGAAGTTTTAGATCCAAACAATATTCGCAAACTATACTGAGATTGCCTTGCAACAGGTTTCCATG
TTCTTCCGATCAGAACCAAAGTGGGAGGTGGTGAACCTTTGAAAGACATAGGTTGGAGAATAAGGAAG
AAATATTTCTTGATGAAGATTAATAATCAGCCAAAGAACGGCTAGTGTTAAGCTGGGCTGACCTTGGC
CCAGACAAGTATTTGTCAGATAAAGATTTTCAGTGTCTAATCAAACCTTCTGCCTTCTGTTTGCACCT
TACATCTATCGGGTTACCTTTGCCACAGCTAATGAATCCTCAGCGTTGCTAATTAGGATGTTAACGAA
AAGGGAACATTGAAGGATCTGATCTACAAGGCAAACCAAAGACCCATTCTAAAGAAGTACTGCAAC
CCTAAGAAGATTCAGGGCTGGAACCTCCAGCAAATAAAAACATATGGACGGCAAATATTAGAGGACTG
AAGTTTCTTATGACAAGGGATCCCTTATGGGCATCTTCACGCCTCAAATGTGATGCTCGATGGGGAC
ACTTGGCCGGTGTGGACCTTGAGAATTCCTTATTGGGCTGCCTTCTTCTACCGATCTTATTTTTCA
CAATTCAGGAAAATCAATACATTGGAAGTGTGGATGTCCACTGCTTTGGCCACTTACTGTATGAAATG
ACTTATGGACGACCGCCAGACTCGGTGCCTGTGGACTCCTTCCCTCCTGCCCGTCCATGGCTGTGGTG
GCCGTGTTGGAGTCTACGCTGTCTTGTGAAGCCTGTAAAAATGGCATGCCTACCATCTCCCGCTCTTA
CAGATGCCATTATTCAGCGATGTTTTACTAACCCTTCTGAAAAACCACAGTTTAAGATCCCTACAAG
TTAAAAGAGGCATTGAGAATTGCCAAAGAATGTATAGAGAAGAGACTAATTGAGGAACAGAAACAGATT
CACCAGCATCGAAGACTGACAAGAGCTCAGTCCCACCATGGATCTGAGGAGGAAAGAAAAAAGAAAG
ATTTTAGCTCGAAAGAAGTCAAACGATCTGCTCTTGAATAAGTGAAGAGCATTAGCGAAGTACAGC
AACTCCAATAATTCAGCAGGATCTGGGGCCAGCTCACCTCTCACGTCGCCGTCATCGCCAACTCCACCC
TCTACATCAGTAGAGCATGCACCATTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001289098

Insert Size:	1548 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289098.1</u>
RefSeq Size:	3053 bp
RefSeq ORF:	1548 bp
Locus ID:	54899
UniProt ID:	<u>Q7Z7A4</u>
Cytogenetics:	3p14.3
Protein Families:	Druggable Genome, Protein Kinase
MW:	58.7 kDa
Gene Summary:	<p>This gene encodes a phox (PX) domain-containing protein which may be involved in synaptic transmission and the ligand-induced internalization and degradation of epidermal growth factors. Variations in this gene may be associated with susceptibility to systemic lupus erythematosus (SLE). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (4) contains an alternate exon in the 3' region, resulting in an early stop codon, compared to variant 1. It encodes isoform d (also known as isoform 4), which is shorter and has a distinct C-terminus, compared to isoform a. Variants 4 and 8 both encode the same isoform (d).</p>