

## Product datasheet for **SC336583**

### CLUL1 (NM\_001289036) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CLUL1 (NM_001289036) Human Untagged Clone
Tag:	Tag Free
Symbol:	CLUL1
Synonyms:	RA337M
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336583 representing NM\_001289036.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAGGACCTGGGACTACAGATGGAGTCTTGTCTCGTTGCCAGACTGGAGTGCCTGCTGCAATCTC
AGCTCACTGCAACCTCTACCTCCAGGTTCAAGCGATTCTCTGCCTCAGCCTCTCGAGTGGCTGGGAC
TATAGTAACAGCGGGAACATGAAGCCGCCACTCTTGGTGTATTATTGTGTCTGCTGTGGTTGAAAGAC
AGTCACTGCGCACCCACTTGGAAGGACAAAACCTGCTATCAGTGAAAACCTGAAGATTTTTCTGAGGTG
GGGAGATAGATGCAGATGAAGAGGTGAAGAAGGCTTTGACTGGTATTAAGCAAATGAAAATCATGATG
GAAAGAAAAGAGAAGGAACACCAATCTAATGAGCACCTGAAGAAATGCAGAGAAGAAAAGCAGGAG
GCCCTGAAACTTCTGAATGAAGTTCAAGAACATCTGGAGGAAGAAGAAAGGCTATGCCGGGAGTCTTTG
GCAGATTCTGGGTGAATGCAGGTCTTGCCTGGAAAATAACTGCATGAGAATTTATAACAACCTGCCAA
CCTAGCTGGTCTCTGTGAAAAATAAGATTGAACGTTTTTCAGGAAGATATATCAATTTCTATTTCTT
TTCCATGAAGATAATGAAAAGATCTCCCATCAGTGAAAAGCTCATTGAGGAAGATGCACAATTGACC
CAATGGAGGATGTGTTACGCCAGTTGACTGTGGATGTGAATTTCTCTTTAACAGGAGTTTAAAGTC
TTCAGACAGATGCAGCAAGAGTTTGACCAGACTTTTCAATCACATTTTATATCAGATACAGACCTAACT
GAGCCTTACTTTTTCCAGCTTTCTCTAAAGAGCCGATGACAAAAGCAGATCTTGAGCAATGTTGGGAC
ATTCCAACTTCTCCAGCTGTTTTGTAATTTCAAGTGTCTCTATTTATGAAAGTGTCAAGTAAACAATT
ACTAAGATGCTGAAGGCAATAGAAGATTACCAAAAACAAGACAAAGCTCCTGACCACGGAGGCTGATT
TCAAAGATGTTACCTGGGAGGACAGAGGACTGTGTGGGAACTTGACCAGAATTTGTCAAGATGTTTC
AAATTTTCAAAAAATGCCAAAAATGTCAGGCTCACCTATCTGAAGACTGTCTGATGTACCTGCTCTG
CACACAGAATTAGACGAGGCGATCAGGTTGGTCAATGTATCCAATCAGCAGTATGGCCAGATTCTCCAG
ATGACCCGGAAGCACTTGAGGACACCGCCTATCTGGTGGAGAAGATGAGAGGGCAATTTGGCTGGGTG
TCTGAACTGGCAAACCAGGCCCCAGAAACAGAGATCATCTTTAATTCAATACAGGTAGTTCCAAGGATT
CATGAAGGAAATATTTCCAAACAAGATGAAACAATGATGACAGACTTAAGCATTCTGCCTTCTCTAAT
TTCACACTCAAGATCCCTCTTGAAGAAAGTGTGAGAGTTCTAACTTCATTGGCTACGTAGTGGCAAAA
GCTCTACAGCATTTTAAGGAACATTTTAAACCTGGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



**ACCN:**

NM\_001289036

<b>Insert Size:</b>	1557 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001289036.1</a></u>
<b>RefSeq Size:</b>	2196 bp
<b>RefSeq ORF:</b>	1557 bp
<b>Locus ID:</b>	27098
<b>UniProt ID:</b>	<u><a href="#">Q15846</a></u>
<b>Cytogenetics:</b>	18p11.32
<b>Protein Families:</b>	Secreted Protein
<b>MW:</b>	60.1 kDa