

Product datasheet for **SC336589**

TUBA6 (TUBA1C) (NM_001303114) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TUBA6 (TUBA1C) (NM_001303114) Human Untagged Clone
Tag:	Tag Free
Symbol:	TUBA1C
Synonyms:	bcm948; TUBA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336589 representing NM_001303114.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACAGCCTGGTTCAATGGGGTGGAGGAGCTAGGGAGGGATGAGTGCTTTGTGTGCTTGAATTAGAT
CTTTCAAATGGATCCTTTCTGAATGCAAACTGTACATCTCTAACTGGATTCTTATTACTTCACCAGG
ACTCTTCAGTCCCTGCGCCTTTTAACACATGCACATCCAGCAAAGCAGAGGAGAACCTGGCTGTGAT
TCAAAGCGTGAGTGCATCTCCATCCACGTTGGCCAGGCTGGTGTCCAGATTGGCAATGCCTGCTGGGAG
CTCTACTGCCTGGAACACGGCATCCAGCCCGATGGCCAGATGCCAAGTGACAAGACCATTGGGGGAGGA
GATGATTCTTCAACACCTTCTTCAGTGAACGGGTGCTGGCAAGCATGTGCCCCGGGCAGTGTTTGTGA
GACTTGAACCCACAGTCATTGATGAAGTTCGCACTGGCACTTACCGCCAGCTCTCCACCCTGAGCAA
CTCATCACAGGCAAGGAAGATGCTGCAATAACTATGCCCGAGGGCACTACACCATTGGCAAGGAGATC
ATTGACCTCGTGTGGACCGAATTCGAAGCTGGCTGACCAGTGCACCGGTCTTCAGGGCTTCTGGTT
TCCACAGCTTTGGTGGGGAACTGGTTCTGGGTTACCTCGCTGCTCATGGAACGTCTCTCAGTTGAT
TATGGCAAGAAGTCCAAGCTGGAGTTCTCCATTTACCGGCGCCCCAGGTTTCCACAGCTGTAGTTGAG
CCCTACAACCTCCATCCTACCACCACACCACCCTGGAGCACTCTGATTGTGCCTTCATGGTAGACAAT
GAGGCCATCTATGACATCTGTCTAGAAAACCTCGATATCGAGCGCCCAACCTACACTAACCTAACCCGC
CTTATTAGCCAGATTGTGCTCCTCCATCACTGCTTCCCTGAGATTTGATGGAGCCCTGAATGTTGACCTG
ACAGAATTCAGACCAACCTGGTGCCCTACCCCGCATCCACTTCCCTCTGGCCACATATGCCCTGTG
ATCTCTGTGAGAAAAGCCTACCAGAACAGCTTACTGTAGCAGAGATCACCAATGCTTGTGAGCCA
GCCAACCAGATGGTAAAATGTGACCCTCGCCATGGTAAATACATGGCTTGCTGCCTGTTATACCGTGGT
GACGTGGTCCCAAAGATGTCAATGCTGCCATTGCCACCATCAAACCAAGCGTACCATCCAGTTTGTG
GATTGGTGCCCGCCACTGGCTTCAAGGTTGGCATTAAATACCAGCCTCCCACTGTGGTGCCTGGCGGAGAC
CTGGCCAAGGTACAGAGAGCTGTGTGCATGCTGAGCAATACCACAGCTGTTGCCGAGGCCTGGGCTCGC
CTGGACCACAAGTTTACCTGATGTATGCCAAGCGTGCCTTTGTTCACTGGTACGTGGGTGAGGGGATG
GAGGAAGGCGAGTTTTCAGAGGCCCGTGGAGCATGGCTGCCCTTGAAGGATTATGAGGAGTTGGA
GCAGATAGTGCTGACGGAGAGGATGAGGGTGAAGAGTATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001303114

Insert Size:	1560 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001303114.1</u>
RefSeq Size:	1707 bp
RefSeq ORF:	1560 bp
Locus ID:	84790
Cytogenetics:	12q13.12
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection
MW:	57.7 kDa
Gene Summary:	<p>Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (a).</p>