

## Product datasheet for **SC336602**

### **ATL3 (NM\_001290048) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ATL3 (NM_001290048) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATL3
Synonyms:	HSN1F
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336602 representing NM\_001290048.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGAGCAGCAAGCCTGGTCCAGTGCAGTGTGTTTGGTTTCAGAAAGATCAACATTCCTTTGAGCTA
GATGAGAAGCCTTGGCCAGCATCCTCTGCAGGACCACATCCGAGATCTTGATGTGGTGGTGGTTTCA
GTGGCTGGTGCCTTCCGAAAGGGCAAGTCTTCATTCTGGATTTTATGCTACGATACTTATATTCTCAG
AAGGAAAGTGCCATTCAAATTGGTTGGGTGACCCAGAAGAACCCTTAACAGATTTTCTGGAGAGGG
GGATCTGATCCAGAAACCACTGGGATCAAATCTGGAGTGAAGTTTTCACTGTGGAGAAGCCAGGTGGG
AAGAAGTTGCAGTTGTTCTGATGGATACCCAGGGGCCATTTGACAGCCAGTCAACTGTAAAGACTGT
GCTACCATCTTTGCTCTAAGCACTATGACTAGTCTGTTTCAGATTTATAATTTATCTCAGAACATTCAA
GAAGATGATCTTCAACAGCTGCAGCTCTCACAGAATACGGTCGTCTGGCAATGGATGAAATTTTCCAA
AAGCCTTCCAGACTGATGTTTTGGTTAGAGATTGGAGTTCCCTTATGAATATAGCTATGGACTC
CAAGGAGGAATGGCATTGTTTGGATAAGCGTTTACAGGTGAAGGAACATCAACATGAAGAAATTCAGAA
GTTTCAAATCACATTCATGTTTCTCGATGTCACCTGCTTTCTTACCACATCCAGGACTCCAG
GTGGCCACAAGCCCTGACTTTGATGGGAAATTAAGATATTGCTGGTGAATTCAGAGCAGTTACAG
GCACTGATACCGTATGTATTAACCCATCTAAGTTAATGGAAAAGGAGATCAATGGCTCAAAGGTCACC
TGTCGGGGACTACTGGAGTATTTAAGGCATATATTAATTTATCAAGGAGAAGATCTGCCTCACCCC
AAGTCCATGCTTCAGGCCACTGCTGAAGCCAACAACCTAGCAGCTGCAGCCTCTGCCAAGGACATTTAT
TATAACAACATGGAAGAGTTTGTGGGGGAGAGAAACCTATTTGTCTCCAGACATTCTAGAGGAGAAG
CACTGTGAATTCAAACAACCTGCTCTGGACATTTAAGAAGACCAAGAAGATGGGTGGGAAGGATTTT
AGCTTTCGTTACCAGCAGGAGCTGGAGGAGGAAATCAAGGAATTATATGAGAACTCTGCAAGCACAA
GGTAGCAAGAACGTCTTACGACCTCCGAACCCCTGCAGTGTGTTACGGGCATTGTAGCTTTGTAC
ATAGCCTCAGGCCTCACTGGCTTCATAGGTCTTGGAGTTGTAGCCAGTTGTTCAACTGTATGGTTGGA
CTACTGTTAATAGCACTCCTCACCTGGGGCTACATCAGGTATTCTGGTCAATATCGTGAGCTGGCGGA
GCTATTGATTTTGGTGCCGCATATGTGTTGGAGCAGGCTTCTTCTCATATCGGTAATTCCTCAGGCC
ACTGTGAGGGATGCAGTTGTTGGAAGACCATCCATGGATAAAAAAGCTCAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



**ACCN:**

NM\_001290048

<b>Insert Size:</b>	1572 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001290048.1</a></u>
<b>RefSeq Size:</b>	7156 bp
<b>RefSeq ORF:</b>	1572 bp
<b>Locus ID:</b>	25923
<b>UniProt ID:</b>	<u><a href="#">Q6DD88</a></u>
<b>Cytogenetics:</b>	11q13.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	58.8 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of a family of dynamin-like, integral membrane GTPases. The encoded protein is required for the proper formation of the network of interconnected tubules of the endoplasmic reticulum. Mutations in this gene may be associated with hereditary sensory neuropathy type IF. Alternatively spliced transcript variants that encode distinct isoforms have been described. [provided by RefSeq, Feb 2014]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at an alternate start codon compared to variant 1. The encoded protein (isoform 2) has a distinct N-terminus and is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>