

Product datasheet for **SC336619**

SWAP70 (NM_001297714) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SWAP70 (NM_001297714) Human Untagged Clone
Tag:	Tag Free
Symbol:	SWAP70
Synonyms:	HSPC321; SWAP-70
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336619 representing NM_001297714.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGGAGCTTGAAGGAGGAGCTGCTCAAAGCCATCTGGCAGCCTTACCCGCACTCGACCAGGACCAC
AGCGGCAAGGCTCCAAGTCCCAGCTCAAGTCTTTCCATAACCTGTGCACGGTCTGAAGTTCTCT
CATGACCCAGTTGCCCTTGAAGAGCACTTCAGGGATGATGATGAGGGTCCAGTGTCCAACAGGGCTAC
ATGCCTTATTTAAACAGGTTTATTTTGGAAAAGATTGAATACCTGCTTAAGAAGCTTACAGAAGCTATG
GGAGGAGTTGGCAGCAAGAACAATTTGAACATTATAAAATCAACTTTGATGACAGTAAAAATGGCCTT
TCTGCATGGGAATTATTGAGCTTATTGAAAATGGACAGTTTAGCAAAGGCATGGACCGGCAGACTGTG
TCTATGGCAATTAATGAAGTCTTTAATGAACCTTATTAGATGTGTTAAAGCAGGGTTACATGATGAAA
AAGGGCCACAGACGAAAAACTGGACTGAAAGATGGTTTGTACTAAAACCAACATAATTTCTTACTAT
GTGAGTGAGGATCTGAAGGATAAGAAAGGAGACATTCTCTGGATGAAAATTGCTGTGTAGAGTCTTTG
CCTGACAAAGATGAAAAGAAATGCCTTTTCTCGTAAAATGTTTTGATAAGACTTTTGAATCAGTGCT
TCAGATAAGAAGAAGAAACAGGAGTGGATTCAAGCCATTCTACTATTTCATCTGTTGAAGCTGGGC
AGCCCTCCACCACAAAAGAAGCCGCCAGCGTCGGAAAAGAACTCCGGAAGAAGCAGCTGGCTGAACAA
GAGGAAGTGGAGCGACAATGAAGGAAGTCCAGGCCGCCAACGAAAGCAAGCAGCAGGAGCTGGAGGCC
GTGCGGAAGAAACTGGAGGAAGCAGCATCTCGTGCAGCAGAAGAGGAAAAGAAACGCCTTACAGCTCAA
GTGGAAGTTCAGGCCAGGTTTCCAGCAGAGCTGGAAAAGAGAGAAGCTTATCAGACAGCAGATGGAAGAA
CAGGTTGCTCAAAGTCTCTGAACTGGAACAGTATTTACAGCGAGTACGGGAGCTGGAAGACATGTAC
CTAAAGCTGCAGGAGGCTCTTGAAGATGAGAGACAGGCCCGCAAGATGAAGAGACAGTGGCGAAGCTT
CAGGCCAGGTTGTTGGAGGAAGAGTCTTCCAAGAGGGCTGAACTAGAAAAGTGGCACTTGGAGCAGCAG
CAGGCCATTTCAGACAACCGAGGCGGAGAAGCAGGAGTTGGAGAATCAGCGTGTCCCTGAAGGAACAGGCC
CTGCAGGAGGCCATGGAGCAGCTGGAGCAGCTTGAAGTTAGAACGGAAGCAAGCACTTGGAGCAGTACGAG
GAAGTTAAAAAGAAGCTGGAGATGGCAACTAATAAGACCAAGAGCTGGAAGGACAAAAGTGGCCCATCAT
GAAGGATTAATTCGACTGATAGAACCAGGTTCAAAGAACCCTCACCTGATCACTAACTGGGGACCTGCA
GCTTTCAGTGGGAGGAGCAACTTGAAGAGAGAGAGAAGAACTGGAAAAGAGAAAAAGACCACGGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001297714

Insert Size:	1584 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001297714.1</u>
RefSeq Size:	4708 bp
RefSeq ORF:	1584 bp
Locus ID:	23075
UniProt ID:	<u>Q9UH65</u>
Cytogenetics:	11p15.4
MW:	62 kDa
Gene Summary:	<p>Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which, independently of RAS, transduces signals from tyrosine kinase receptors to RAC. It also mediates signaling of membrane ruffling. Regulates the actin cytoskeleton as an effector or adapter protein in response to agonist stimulated phosphatidylinositol (3,4)-bisphosphate production and cell protrusion (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>