

Product datasheet for **SC336620**

RCBTB2 (NM_001286831) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RCBTB2 (NM_001286831) Human Untagged Clone
Tag:	Tag Free
Symbol:	RCBTB2
Synonyms:	CHC1L; RLG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336620 representing NM_001286831.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTAGATGTGGAAAGTGGCCAATTTTTCCCTTTGTTCTGAAGAAGAACTACAGTTAATTCGTGAG
GCTTGTGCTTTTGGCAGTGTGGCAATGAAGTTTTATACACTACAGTAAATGATGAGATTTTTGTGCTT
GGCAGAACTGCTGTGGCTGTTTGGGGTTAGGTGACGTCCAGAGCACCATTGAACCTCGGAGACTGGAT
TCTTTAAATGGCAAAAAAATAGCCTGCCTCAGCTATGGGAGTGGTCCACATATTGCTTCTGCAACAACA
GAAGGAGAAGTCTTTACCTGGGGTCATAATGCTTATAGCCAGCTGGGCAATGGGACAATAATCATGGT
TTAGTGCCTGTGATCTCTACTAATCTGTCAAACAACAAGTCATTGAAGTTGCCTGTGGGTCTTAC
CATTCTTTGGTGTAACTCTGATGGAGAGGATTTGCCTGGGGTTATAAATACTCTGGGCAGGTAGGA
TCTGGATCAACAGTTAATCAGCCAATCCCTCGAAGAGTCACTGGTGCCTACAAAATAAAGTAGTTGTG
ACCATAGCATGTGGCAGATGTGCTGCATGGCAGTAGTAGACACGGGGAGGTCTATGTCTGGGGTTAC
AACGAAAACGGGCAGCTTGGACTCGGCAACAGTGGCAACCAGCCAACCCCTTGCAGAGTGGCAGCTTTG
CAAGGCATCCGTGTCAGAGGGTCGCTGTGGCTACGCACACACATTAGTATTAACAGATGAAGGCCAA
GTGATGCTTGGGGCGCCAATTCTTATGGGCAGTTGGGCACTGGCAATAAAAGCAACCAGTCCTATCCT
ACTCCTGTCACTGTGAAAAGGACAGGATTATCGAGATTGCAGCCTGTCACTCCACACACAGTCTGCG
GCCAAGACGCAGGGTGGGCACGTGTACATGTGGGGCCAGTGCCGGGGTCAGTCCGTGATCCTCCCGCAC
CTCACCCACTTCTCTGCACTGACGACGTGTTGCCTGCTTGGCCAGCCCGCCGTACGTGGCGCCTC
CTCTCCGTGGAACCTGATGACCCTCACAGTGGCTGAGTCACTGAAGAGGGAATTTGACAACCCGGAC
ACTGCAGACCTGAAGTTCTAGTTGATGGAAAGTACATTTATGCACATAAAGTCCCTTCTCAAGATTCGA
TGTGAGCATTTCGTTTCGTCATTGGAAGATAACGAGGATGATATTGTAGAAATGAGTGAATTTTCATAT
CCTGTTTACCGGGCCTTCTGGAATACCTATACACAGACAGCATCAGCCTTCTCCTGAGGAGGCAGTA
GGACTGCTAGACTTGGCTACATTTTATAGAGAAAATCGTTTAAAAAGCTCTGCCAACAATAATCAAG
CAAGGCATCTGCGAGGAGAATGCCATCGCTCTGCTCTCGGCTGCGGTGAAGTATGATGCACAGGATTTA
GAAGAATTCTGCTTCAGGTTTTGCATAAACCATCTGACTGTAGTAACACAACATCAGGTTTTGCAGAA
ATGGACCATGATCTCTGAAGAATTTATCAGCAAAGCAAGCAGAGTTGGAGCCTTTAAAAATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001286831

Insert Size:	1584 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286831.1</u>
RefSeq Size:	3020 bp
RefSeq ORF:	1584 bp
Locus ID:	1102
UniProt ID:	<u>O95199</u>
Cytogenetics:	13q14.2
MW:	57.8 kDa
Gene Summary:	<p>This gene encodes a protein containing two C-terminal BTB/POZ domains that is related to regulator of chromosome condensation (RCC). The encoded protein may act as a guanine nucleotide exchange factor. This gene is observed to be lost or underexpressed in prostate cancers. There is a pseudogene of this gene on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2013]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (3) has a shorter N-terminus than isoform 1. Variants 3 and 7 both encode the same isoform (3).</p>