

Product datasheet for **SC336659**

YTHDF3 (NM_001277816) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	YTHDF3 (NM_001277816) Human Untagged Clone
Tag:	Tag Free
Symbol:	YTHDF3
Synonyms:	DF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336659 representing NM_001277816.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCAGATCCATACATGCCTAGTACTATGCTCCATCCATTGGATTCCATATTCTCTTGGGAAGCA
GCGTGGTCCACAGCTGGAGACCAGCCTATGCCATATCTGACAACCTATGGACAAATGAGTAATGGAGAA
CATCACTATATACCAGATGGTGTATTTAGTCAACCTGGGGCATTAGGAAATACCCCTCCATTTCTGGT
CAACATGGATTTAACTTTTTTCTGTTAATGCTGATTTCTCTACATGGGGACAAGTGGATCTCAGGGA
CAATCAACACAAAGTTCTGCTTATAGTAGCAGTTATGGCTATCCACCTAGTTCTCTTGGGAGAGCTATT
ACTGATGGACAGGCTGGATTTGGCAATGATACTTTGAGTAAGGTGCCTGGCATTAGCAGTATTGAGCAA
GGCATGACTGGACTGAAAATTGGTGGTACCTGACAGCTGCAGTGACAAAACTGTAGGTACAGCTTTG
AGCAGCAGTGGTATGACTAGCATTGCAACCAATAGTGTGCCCCAGTTAGCAGTGCAGCACCTAAACCA
ACCTCTGGGCTGCCATTGCCAGAAAGCCTGCCAAACCTCAACCGAAACTTAAACCAAGGGCAATGTG
GGAATTGGGGTTCTGCTGTACCACCACCTCTATAAAACACAACATGAATATTGGAACCTGGGATGAA
AAAGGGTCAGTGGTAAAGGCTCCACCAACCAACCAAGTTCTGCCTCTCAAACTATAATCCAGCAGCCT
CAGCCATTAATCAACCACCACCATTGGTGCAAAGCCAACCTGCCTCAACAGCAGCCTCAACCACCACAA
CCACAGCAGCAACAAGGACCTCAGCCACAGGCCAGCCTCACCAAGTGCAGCCTCAACAGCAGCAGCTG
CAGAATCGCTGGGTAGCTCCTCGTAACAGGGGAGCAGGCTTCAACCAGAACAAATGGAGCGGGCAGTGAA
AACTTTGGTTTAGGTGTTGTACCTGTGCTGCTTACCTTCTAGTGTAGAAGTGCATCCCCTGCTGGAA
AAGCTAAAGGCCATAAACTATAATCCCAAAGACTTTGATTGGAATCTGAAGAATGGACGTGTGTTT
ATAATTAAGCTACTCTGAGGATGACATACATCGTCCATTAATACTCTATCTGGTGTAGTACTGAG
CATGGTAATAAGCGTTTGGATGCAGCTTACCGTTCCTGATGGGAAAGGCCACTCTATTTACTCTTC
AGTGTGAATGGCAGTGGACATTTTTGTGGAGTGGCTGAAATGAAGTCTGTTGTGGACTATAATGCGTAT
GCTGGTGTCTGGTCTCAGGATAAGTGGAAAGGGCAAATTTGAAGTTAAATGGATCTTTGTCAAAGATGTT
CCCAATAACCAATTACGGCATATTGCTTAGAAAAAATGACAACAAACCGGTTACCAATTAAGGGGAC
ACTCAAGAGGTACCCCTAGAAAAAGCTAAGCAAGTGTAAAATAATTGCTACTTTCAAGCATACCACC
TCAATCTTTGATGACTTTGCACATTATGAAAAGCGTCAAGAAGAGGAGGAAGCCATGCGTAGGGAGAGA
AATAGAAACAAACAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001277816
Insert Size:	1605 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001277816.1</u>
RefSeq Size:	5085 bp
RefSeq ORF:	1605 bp
Locus ID:	253943
UniProt ID:	<u>Q7Z739</u>
Cytogenetics:	8q12.3
MW:	58.3 kDa
Gene Summary:	<p>This gene encodes a member of the YTH (YT521-B homology) domain protein family. The YTH domain is common in eukaryotes, is often found in the middle of the protein sequence, and may function in binding to RNA. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2013]</p> <p>Transcript Variant: This variant (5) differs in its 5' UTR and uses a downstream in-frame start codon, compared to variant 1. The encoded isoform (c) has a shorter N-terminus, compared to isoform a. Variants 4, 5, 6, and 7 encode the same isoform.</p>