

## Product datasheet for **SC336668**

### SLC2A14 (NM\_001286237) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SLC2A14 (NM_001286237) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC2A14
Synonyms:	GLUT14; SLC2A3P3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336668 representing NM\_001286237.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAAAGACTCCAACCTGTTGAGAGTGAAGTACTGCTGGGTGTGAAACAAGGAGATGAAATGCGGCAT
TTCTTTTTTTCTTCTCAGACCTCCACTCTGGAGAAGAGTCAGAATGGGGCGTCGGGGAGGAGTCAACC
CCAGCTCTGATCTTTGCCATCACAGTTGCTACAATCGGCTCTTCCAGTTTGCTACAACACTGGGGTC
ATCAATGCTCCTGAGACGATCATAAAGGAATTTATCAATAAACTTTGACGGACAAGGCAATGCCCT
CCCTCTGAGGTGCTGCTCACGAATCTCTGGTCTTGTCTGTGGCCATATTTCCGTGGGGGTATGATC
GGCTCCTTTCCGTGGACTCTTTGTTAACCGCTTTGGCAGGCGCAATCAATGCTGATTGTCAACCTG
TTGGCTGCCACTGGTGGCTGCCTTATGGGACTGTGAAAATAGCTGAGTCAGTTGAAATGCTGATCTG
GGCCGCTTGGTTATTGGCTCTTCTGCGGACTCTGCACAGTTTTGTGCCATGTACATTGGAGAGATC
TCGCCTACTGCCCTGAGGGTGCCTTTGGCACTCTCAACCAGCTGGGCATAGTTATTGGAATTCTGGTG
GCCAGATCTTTGGTCTGGAACATCCTTGGGTCTGAAGAGCTATGGCCGGTGTATTAGGCTTTACC
ATCCTTCCAGCTATCCTGCAAAGTGCAGCCCTCCATGTTGCCCTGAAAGTCCCAGATTTTTGCTCATT
AACAGAAAAAAGAGGAGAATGCTACGCGGATCCTCCAGCGGTTGTGGGGCACCCAGGATGTATCCAA
GACATCCAGGAGATGAAAGATGAGAGTGAAGGATGTCACAAGAAAAGCAAGTCAACCGTGTGGAGCTC
TTAGAGTGTCCAGCTACCGACAGCCATCATCTTTCCATTGTGCTCCAGCTCTCTCAGCAGCTCTCT
GGGATCAATGCTGTGTTCTATTACTCAACAGGAATCTTCAAGGATGCAGGTGTTCAACAGCCCATAT
GCCACCATCAGCGCGGTGTGGTTAATACTATCTTCACTTTACTTTCTCTATTTCTGGTGGAAAGGGCA
GGAAGAAGGACTCTGCATATGATAGCCCTTGGAGGGATGGCTTTTTGTTCCAGCTCATGACTGTTTCT
TTGTTATTAAGAATCACTATAATGGGATGAGCTTTGTCTGTATTGGGGCTATCTTGGTCTTTGGCC
TGTTTTGAAATTGGACCAGGCCCATTCCTGGTTTATTGTGGCCGAACCTTTCAGCCAGGGCCCCCGC
CCAGCTGCGATGGCAGTGGCCGGCTGCTCCAACCTGGACCTCCAACCTTCTAGTCGGATTGCTCTTCCCC
TCTGCTGTTACTATTTAGGAGCTACGTTTTTATTATCTTACCAGGCTTCTCATTACCTTCTTGGCC
TTTACCTTCTTCAAAGTCCCTGAGACCCGTGGCAGGACTTTTGGAGATATCACACGGCCTTTGAAGGG
CAGGCACAGGTGCAGATAGATCTGGGAAGGACGGCTCATGGGGATGAACAGCATCGAGCTGCTAAG
GAGACCACCAATGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



<b>ACCN:</b>	NM_001286237
<b>Insert Size:</b>	1608 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001286237.1</a></u>
<b>RefSeq Size:</b>	3507 bp
<b>RefSeq ORF:</b>	1608 bp
<b>Locus ID:</b>	144195
<b>UniProt ID:</b>	<u><a href="#">Q8TDB8</a></u>
<b>Cytogenetics:</b>	12p13.31
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	58.2 kDa
<b>Gene Summary:</b>	<p>Members of the glucose transporter (GLUT) family, including SLC2A14, are highly conserved integral membrane proteins that transport hexoses such as glucose and fructose into all mammalian cells. GLUTs show tissue and cell-type specific expression (Wu and Freeze, 2002 [PubMed 12504846]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation from an alternate start codon, compared to variant 1. The encoded isoform (d) has a distinct N-terminus and is longer than isoform a.</p>