

## Product datasheet for **SC336678**

### FMO1 (NM\_001282692) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FMO1 (NM_001282692) Human Untagged Clone
Tag:	Tag Free
Symbol:	FMO1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336678 representing NM\_001282692.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGGAGAACATGGCCAAGCGAGTTGCCATTTGTGGAGCTGGGGTCAAGCGCCTGGCCTCCATCAAG
TGCTGTCTGGAAGAAGGACTGGAGCCACCTGCTTTGAGAGGAGCGATGACCTTGGGGGGCTGTGGAGA
TTCACCGAACATGTTGAAGAAGGCAGAGCCAGTCTCTACAAGTCTGTGGTTTCCAACAGCTGCAAGGAG
ATGCTTGTACTCAGACTTTCCATTCCAGAAGATTATCCAAACTATGTGCCAAATTCTCAATTCCTG
GAATATCTCAAAATGTATGCAAACCACTTTGACCTTCTGAAACACATTCAATTCAAGACCAAAGTCTGC
AGTGTAACAAAATGCTCAGATTCTGCTGTCTCTGGCCAATGGGAGGTGGTCACTATGCATGAAGAGAAG
CAAGAGTCAGCCATCTTTGATGCTGTCTGCTGCTGACTGGCTTTCTTACTAATCCTTATTTGCCACTG
GATTCCTTTCCAGGTATTAATGCCTTTAAAGGCCAGTACTTTCATAGCCGGCAATATAAGCATCCAGAT
ATATTTAAGGACAAGAGAGTCTTGTGATTGGAATGGGAAATTCGGCACAGACATTGCTGTGGAGGCC
AGCCACCTGGCGAAAAGGTGTTCTCAGCACCACCGAGGGGGATGGGTGATCAGCCGAATCTTTGAC
TCGGGCTACCCATGGGACATGGTGTTCATGACAGCTTTTCAGAACATGTTGAGAAATTCCTCCCAACC
CCAATTGTGACTTGGTTGATGGAGCGAAAGATAAACAACCTGGCTCAATCATGCAAATTACGGCTTAATA
CCAGAAGACAGGACTCAGCTGAAAGAGTTTGTGCTAAATGATGAGCTCCCAGGACGCATCATCACTGGG
AAAGTGTTCATCAGGCCAAGCATAAAAGAGGTAAGGAAAACCTCTGTCATATTTAACAATACTTCAAAG
GAAGAGCCTATTGACATCATTGTCTTTGCCACTGGATACACATTTGCTTTCCCCTTCTTGATGAGTCT
GTAGTAAAGTTGAAGATGGCCAGGCCTCACTGTACAAGTATATCTTCCCTGCACATCTGCAAAAAGCCA
ACCCTGGCCATTATTGGCCTCATCAAACCCTTGGGCTCCATGATACCTACAGGAGAAACACAAGCTCGG
TGGGCTGTTTCGAGTCTGAAAGGTGTAATAAGTTACCACCACCAAGTGCATGATAGAGGAAATTAAT
GCAAGGAAAGAAAACAAGCCAGTTGGTTTGGCTTGTGCTACTGCAAGGCTTTACAATCAGATTATATC
ACATACATAGATGAACTCCTGACCTATATCAATGCAAAAACCAACCTGTTCTCTATGCTCCTAACGGAT
CCACATCTGGCTCTGACCGTCTTCTTTGGCCATGCTCACCATACCAGTCCGCTTGACTGGCCAGGA
AAATGGGAAGGAGCCAGAAATGCCATCATGACCCAGTGGGACCGAACATTCAAGGTATCAAAAGCTCGA
GTTGTACAAGAGTCTCCATCTCCCTTTGAAAGTTTCTTAAAGTCTTTAGCTTTCTGGCTTTGCTTGTG
GCTATTTTCTGATTTTCTATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_001282692
<b>Insert Size:</b>	1611 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001282692.1</a>
<b>RefSeq Size:</b>	2198 bp
<b>RefSeq ORF:</b>	1611 bp
<b>Locus ID:</b>	2326
<b>UniProt ID:</b>	<a href="#">Q01740</a>
<b>Cytogenetics:</b>	1q24.3
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Drug metabolism - cytochrome P450
<b>MW:</b>	60.8 kDa
<b>Gene Summary:</b>	<p>Metabolic N-oxidation of the diet-derived amino-trimethylamine (TMA) is mediated by flavin-containing monooxygenase and is subject to an inherited FMO3 polymorphism in man resulting in a small subpopulation with reduced TMA N-oxidation capacity resulting in fish odor syndrome Trimethylaminuria. Three forms of the enzyme, FMO1 found in fetal liver, FMO2 found in adult liver, and FMO3 are encoded by genes clustered in the 1q23-q25 region. Flavin-containing monooxygenases are NADPH-dependent flavoenzymes that catalyzes the oxidation of soft nucleophilic heteroatom centers in drugs, pesticides, and xenobiotics. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).</p>