

Product datasheet for **SC336686**

ZNF75A (NM_001302109) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF75A (NM_001302109) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF75A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336686 representing NM_001302109.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATGATGGTAGATCTGAAAGTGGCTGCGTACTTGGACCCTCAGATCAGGGCTTTGTGGGAGACCAAG
GGCCTGCAAGAGAGAGCTCCGGTCAGAGTAAAAATCTCCTCAATGGACTGTCTCGATCCTAAGAGC
TCTTGCTGGCACTTCCGGAATTTACCTATGATGAAGCAGGTGGACCCCGTGAGGCTGTAGCAAACTT
CAAGAATTATGTCATCTATGGCTGAAGCCAGAGATCCACTCAAAAGAGCAGATACTGGAAGTCTGGTG
CTGGAGCAGTTCTGACTATTCTGCCAGGGAGACACAGACCCAGATGCAGAAGCACCATCCACAGAGC
ATTGAGGAGGCTGTGGCTCTGGTAGAACACTTGCAGAGGGAATCTGGTCAAACATGGAATGGGGTTGCA
GTCCATGAGCTGGGAAAGGAGGCAGTGTCTTGGGAGAAAACAGCAGAGGCCTCAAGTTTCGGGCTGAAG
CCAACAGAGTCCCAACAGTGGGCGTATCCAAGATGAAGAATTTTGAATACATACGAGGGTCTGCAA
GAACAGCTCAGCAGGAATACTCATAAAGAGACTGAGCCTGTGTATGAGAGGGCTGTGCCTACTCAACAG
ATTCTAGCTTTTCTGAGCAAACAAACACCAAGACTGGACAGTGACACCTGAGCAGCTTGCCTGAG
TCCAGAGCTTGTGACATTTGAAGAAGTGCCATGTATTTTCCAGGAAGAATGGGAGTTATTGGAT
CCCACTCAGAAGGCCCTCTACAATGATGTAATGCAGGAAAATATGAGACTGTACTCTCTAGCATTG
TTTGTGCTCCCCAAACCTAAAGTATCTCTGTCTAGAGCAAGGGGAAGAGCCATGGGTTCAAGTATCC
CCGGAGTTAAGGATAGTGCCGAAAAATCTCTACAGGGTTAAAGCTCAAAAACGACACTGAAAATCAT
CAGCCTGTGTCTTTCTGACTTAGAAAACAAGCATCAGCAGGCGTCATATCAAAAAGGCCAAAGTA
AAAGTTCCCCAGAAAACAGCAGGCAAAGAAAATCATTTTGATATGCACAGAGTGGGAAAATGGCACCAA
GATTTTCCAGTGAAGAAAAGAAAACCTTTCAACCTGGAACAAGAGCTGTCAAACCTTATGGATCGT
CACAAGAAAGATTGTGCAAGAGAGAAGCCTTTAAATGTCAGGAATGTGGGAAAACCTTCAGAGTTAGC
TCTGACCTTATTAAGCACCAAAGAATTCACACTGAAGAGAAAACCTATAAATGTCAACAGTGTGATAAG
AGGTTTAGATGGAGTTCAGATCTTAATAAGCACTTAACAACACACCAAGGAATAAAACCATATAAATGT
TCATGGTGTGGGAAAAGCTTCAGTCAAAATACAAATTTACATACACACCAAGAAGTATACAGGAGAA
AAGCCCTTCACATGTCATGAATGTGAAAAAAATTCAGTCAGAACTCCACCTTATTAACACCGGAGA
ACCCACACAGGTGAGCAGCCATACTTGTAGCATATGCAGGAGAAAACCTCAGCAGGCGGTCAAGCCTT
CTTAGACACCAGAACTCCACCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001302109

Insert Size: 1614 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001302109.1](#)

RefSeq Size: 2692 bp

RefSeq ORF: 1614 bp

Locus ID: 7627

UniProt ID: [Q96N20](#)

Cytogenetics: 16p13.3

Protein Families: Transcription Factors

MW: 62.1 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]