

Product datasheet for **SC336693**

SPARCL1 (NM_001291977) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPARCL1 (NM_001291977) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPARCL1
Synonyms:	MAST 9; MAST9; PIG33; SC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336693 representing NM_001291977.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTGAGCCTCAGGAGAAAAAACTCTCAGAGAACACTGATTTTTGGCTCCTGGTGTAGTTCCTTC
ACAGATTCTAACCAACAAGAAAGTATCACAAGAGAGAGGAAAACCAAGAACAACCTAGAAATTATTCA
CATCATCAGTTGAACAGGAGCAGTAAACATAGCCAAGGCCTAAGGGATCAAGGAAACCAAGAGCAGGAT
CCAAATATTTCCAATGGAGAAGAGGAAGAAGAAAAAGAGCCAGGTGAAGTTGGTACCCACATGATAAC
CAAGAAAGAAAGACAGAATTGCCAGGGAGCATGCTAACAGCAAGCAGGAGGAAGACAATACCCAATCT
GATGATATTTTGAAGAGTCTGATCAACCAACTCAAGTAAGCAAGATGCAGGAGGATGAATTTGATCAG
GGTAACCAAGAACAAGAAGATAACTCCAATGCAGAAATGGAAGAGGAAAAGCATCGAACGTCAATAAG
CACATTCAAGAACTGAATGGCAGAGTCAAGAGGGTAAAACTGGCCTAGAAGCTATCAGCAACCACAAA
GAGACAGAAGAAAAGACTGTTTCTGAGGCTCTGCTCATGGAACCTACTGATGATGGAATACCACGCCC
AGAAATCATGGAGTTGATGATGATGGCGATGATGATGGCGATGATGGCGGCACTGATGGCCCCAGGCAC
AGTGCAAGTGATGACTACTTCATCCAAGCCAGGCCTTTCTGGAGCCGAGAGAGCTCAATCCATTGCC
TATCACCTCAAAATTTGAGGAGCAAGAGAAAAAGTACATGAAAATGAAAATATAGGTACCACTGAGCCT
GGAGAGCACAAGAGGCCAAGAAAGCAGAGAACTCATCAATGAGGAGGAAACGTCAAGTGAAGGCAAC
ATGAGGGTGATGCTGTGGATTCTTGATGAGCTTCCAGTGTAAAAAGAGGCCACATCTGTAAGGCAGAC
CAACAGGGAAAACCTCACTGTGTCTGCCAGGATCCAGTGACTTGTCTCCAACAAAACCCCTTGATCAA
GTTTGTGGCACTGACAATCAGACCTATGCTAGTTCTCTGCATCTATTCGCTACTAAATGCAGACTGGAG
GGGACCAAAAAGGGGCATCAACTCCAGCTGGATTATTTTGGAGCCTGCAAAATCTATTCTACTTGTACG
GACTTTGAAGTGATTCAGTTTCCTCTACGGATGAGAGACTGGCTCAAGAATATCCTCATGCAGCTTTAT
GAAGCCAACTCTGAACACGCTGGTTATCTAAATGAGAAGCAGAGAAATAAAGTCAAGAAAATTTACCTG
GATGAAAAGAGGCTTTTGGCTGGGGACCATCCATTGATCTTCTCTTAAGGGACTTAAAGAAAACTAC
CACATGTATGTATCCTGTGCACTGGCAGTTTAGTGAACTTGACCAACACCCTATGGATAGAGTCTTG
ACACATTCTGAACTTGCTCCTCTGCGAGCATCTCTGGTCCCATGGAACACTGCATAACCCGTTTCTTT
GAGGAGTGTGACCCCAACAAGGATAAGCACATCACCTGAAGGAGTGGGGCCACTGCTTTGGAATAAA
GAAGAGGACATAGATGAAAATCTCTTGTTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001291977
Insert Size:	1620 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291977.1</u>
RefSeq Size:	2670 bp
RefSeq ORF:	1620 bp
Locus ID:	8404
UniProt ID:	<u>Q14515</u>
Cytogenetics:	4q22.1
Protein Families:	Druggable Genome, Secreted Protein
MW:	61.8 kDa