

Product datasheet for **SC336700**

MPP2 (NM_001278371) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MPP2 (NM_001278371) Human Untagged Clone
Tag:	Tag Free
Symbol:	MPP2
Synonyms:	DLG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336700 representing NM_001278371.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGCAAGTCCTGGACAACCTGGGATCCCTCCCAGTGCCACGGGGGCTGCAGAGCTGGACCTGATC
TTCCCTTCGAGGCATTATGGAAAGTCCCATAGTAAGATCCCTGGCCAAGGCCCATGAGAGGCTGGAGGAG
ACGAAGCTGGAGGCCGTGAGAGACAACAACCTGGAGCTGGTGCAGGAGATCCTGCCGGACCTGGCGCAG
CTGGCTGAGCAGAGCAGCACAGCCGCCGAGCTGGCCACATCCTCCAGGAGCCCCACTTCCAGTCCCTC
CTGGAGACGCACGACTCTGTGGCCTCAAAGACCTATGAGACACCACCCCCAGCCCTGGCCTGGACCCT
ACATTCAGCAACCAGCCTGTACCTCCCGATGCTGTGCGCATGGTGGGCATCCGCAAGACAGCCGGAGAA
CATCTGGTGTAACGTTCCGCGTGGAGGGCGCGAGCTGGTGATCGCGGCATTCTGCATGGGGCATG
GTGGCTCAACAAGCCTGCTGCATGTGGTGACATCATCAAGGAGGTGAACGGGCAGCCAGTGGGCAGT
GACCCCGCGCACTGCAGGAGCTCCTGCGCAATGCCAGTGGCAGTGTATCCTCAAGATCCTGCCAGC
TACCAGGAGCCCCATCTGCCCGCCAGGATTTTGTGAAATGTCACCTTTGACTATGACCCGGCCCGAGAC
AGCCTCATCCCCTGCAAGGAAGCAGGCCTGCGCTTCAACGCCGGGACTTGTCTCCAGATCGTAAACCAG
GATGATGCCAACTGGTGGCAGGCATGCCATGTCGAAGGGGGCAGTGTGGGCTCATTCCCAGCCAGCTG
CTGGAGGAGAAGCGAAAGCATTGTCAAGAGGGACCTGGAGCTGACACCAAACCTCAGGGACCCTATGC
GGCAGCCTTTTCAGGAAAGAAAAGAAGCGAATGATGATTTGACCACCAAGAATGCAGAGTTTGACCGT
CATGAGCTGCTCATTTATGAGGAGGTGGCCCGCATGCCCCGTTCCGCCGAAAACCTGGTACTGATT
GGGGCTCAGGGCGTGGGACGGCGCAGCCTGAAGAAACAGCTCATCATGTGGGATCCAGATCGCTATGGC
ACCACGGTGCCCTACACCTCCCGCGGCCGAAAAGACTCAGAGCGGAAGGTGAGGGTTACAGCTTTGTG
TCCCGTGGGGAGATGGAGGCTGACGTCCTGTGCTGGGCGCTACTGGAGCATGGCGAATACGAGGGCAAC
CTGTATGGCACACGTAATTGACTCCATCCGGGGCGTGGTGCCTGCTGGGAAGGTGTGCGTGTGATGTC
AACCCCGAGCGGTGAAGGTGCTACGAACGGCCGAGTTTGTCCCTTACGTGGTGTTCATCGAGGCCCA
GACTTCGAGACCTGCGGGCCATGAACAGGGCTGCGCTGGAGAGTGAATATCCACCAAGCAGCTCAGG
GAGGGGACCTGAGACGGACAGTGGAGGAGAGCAGCCGCATCCAGCGGGGCTACGGGCACTACTTTGAC
CTCTGCCTGGTCAATAGCAACCTGGAGAGGACCTTCCGCGAGCTCCAGACAGCCATGGAGAAGCTACGG
ACAGAGCCCCAGTGGGTGCTGTCAGCTGGGTGTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001278371
Insert Size:	1626 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001278371.1
RefSeq Size:	4282 bp
RefSeq ORF:	1626 bp
Locus ID:	4355
UniProt ID:	Q14168
Cytogenetics:	17q21.31
Protein Families:	Druggable Genome
MW:	60.5 kDa
Gene Summary:	<p>Palmitoylated membrane protein 2 is a member of a family of membrane-associated proteins termed MAGUKs (membrane-associated guanylate kinase homologs). MAGUKs interact with the cytoskeleton and regulate cell proliferation, signaling pathways, and intracellular junctions. Palmitoylated membrane protein 2 contains a conserved sequence, called the SH3 (src homology 3) motif, found in several other proteins that associate with the cytoskeleton and are suspected to play important roles in signal transduction. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 5' coding region and uses a downstream start codon compared to variant 1. It encodes isoform 2 which has a shorter N-terminus compared to isoform 1. Variants 2, 5, and 7 encode the same protein.</p>