

Product datasheet for **SC336715**

ZNF502 (NM_001282880) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: ZNF502 (NM_001282880) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF502
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC336715 representing NM_001282880.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTTGAATATGCAAGGAGCTGAAGAGAGAGACATTAGAAGAGAGACTTGTCCAGGCTGGGTAAACAAG
AACAAAGCTGCTCTGGAGCAGGATGTCTGTAAAATTGACTCATCAGGGATAGTAGTAAAGAGGTTCCAA
GAGGATGAATACCAAGATTCTACATTTGAAGAAAAATATGCATGTGAGGGCATGAAGGAAAACTCTCCT
AGGGAGATTGCTGAATCATGCCTTTTCCAGGAAGGAGGTTTTGGGAGAATAACTTTCATCCACAAAGAA
GCACCCCTGAAATTATTAGTCAAGGATATAATTTTGAGAAAAGCTTGTCTTTCACCTCAAGCCTTGTT
ACACGCTCAGGGTTTTACAGAAGAGAGTCTGCATCAGTGGGAAACAAGTAATATACAAACCAATGAT
ATTTTCAGACCAAAGTAAATGTCCAACCTCTCTGCACACAGAAAAAATCTTGAAATGTAATGAATGTGGA
AAAACCTTTACTCAGAGCTCATCCCTTACCCAACATCAGAGAACTCATACTGGAGAGAGACCCTACACA
TGTGAGGAATGTGGGAAAGCCTTTAGTCGTAGTTCATTCTTGTCAACATCAAAGAATTCACACTGGA
GTGAAACCATATGGATGTGAGCAGTGTGGGAAACATTTTCGATGTCGATCATTCTTACTCAGCATCAA
AGAATTCACACTGGAGAGAAACCTTATAAATGCAATGAATGTGGGAATTCCTTCCGCAATCACTCACAT
CTCACTGAACACCAGAGAATTCACACTGGAGAGAAACCTTATAAATGCAATAGGTGTGGGAAAGGCATTC
AATCAGAATACACACCTTATTCATCATCAGAGAATTCACACTGGTGAAGCCTTACATATGCAGTGAA
TGTGGCTCTTCTTTTCGAAAACACTCAAATCTTACGCAACATCAGAGAATTCACACTGGGAAAAACCC
CATAAATGTGACGAATGTGGGAAAACCTTTCAAAACAAAGGCAAACCTCTCTCAGCATCAGAGAATTCAT
AGTGGAGAGAAAACCTTATAAATGTAAGAATGTGGCAAAGCCTTTTGTGAGGCCATCTCTTATTTAAA
CACCAGCGAATTCATACTGGAGAAAAACCATATAAGTGTAAAGAATGTGGCAAAGCGTTTACTCAGAGC
ACCCCACTCACTAAACATCAGAGAATACATACAGGGGAGAGACCCTACAAATGCAGTGAATGTGGTAAA
GCCTTCATTTCAGAGCATTGCTTATTCGGCACCAGAGAAGTACACTGGAGAAAAACCTTATAAATGC
AATGAATGTGGAAAGGCTTTAATCAGAACACCTGCCTCACTCAGCATATGAGAATTCATACTGGAGAG
AAGCCCTATAAATGTAAGAATGTGGGAAAGCCTTTGCTCATAGCTCATCTTACTGAACATCATAGA
ACTCACACTGGTGAAGCTCTATAAATGTAGTGAGTGTGAGAAAACCTTCGCAAGTATGCACACCTT
AGTGAACATTACAGAATTCACACTGGTGAAGCCTTATGAGTGTATTGAGTGTGAAAGTCTTTCAGA
CATAGTTCAGTCTTTTCAGACATCAGAAAACCTTCACAGTGGTGACTAA
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001282880

Insert Size: 1635 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282880.1](#)

RefSeq Size: 3407 bp

RefSeq ORF: 1635 bp

Locus ID: 91392

UniProt ID: [Q8TBZ5](#)

Cytogenetics: 3p21.31

Protein Families: Transcription Factors

MW: 62.9 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (5) represents the longest transcript. Variants 1-5 encode the same protein.