

Product datasheet for **SC336717**

POTEB (NM_001277304) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POTEB (NM_001277304) Human Untagged Clone
Tag:	Tag Free
Symbol:	POTEB
Synonyms:	A26B1; CT104.5; POTE-15; POTE15; POTEB2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336717 representing NM_001277304.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGTGGCTGAGGTTTGTCAATGCCCGCTGCCTCTGCTGTGAAGAAGCCATTTGATCTCAGGAGCAAG
ATGGGCAAGTGGTGCCACCACCGCTTCCCTGCTGCAGGGGGAGCGGCACGAGCAATGTGGCAGCTTCT
GGAGACCATGACGACTCCTTTATGAAGACACTCAGGAGCAAGATGGGCAAGTGGTGTCTACTGCTTC
CCCTGCTGCAGGGGGAGCGGCAAGAGCAACGTGGGCACTTGGGGAGACTACGACGACAGCGCCTTCATG
GAGCCGAGGTACCACGTCGTCGAGAAGATCTGGACAAGCTCCACAGAGCTGCCTGGTGGGTAAAGTC
CCCAGAAAGGATCTCATCGTCATGCTCAGGGACACTGACATGAACAAGAGGGACAAGCAAAAGAGGACT
GCTCTACATTTGGCCTCTGCCAATGAAATTCAGAAGTAGTACAACCTCTGCTGGACAGACGATGTCAA
CTTAACGTCCTTGACAACAAAAAAGGACAGCTCTGATAAAGGCCGTACAATGCCAGGAAGATGAATGT
GTGTTAATGTTGCTGGAACATGGCGCTGATGGAATATTCAAGATGAGTATGGAATACCGCTCTACAC
TATGCTATCTACAATGAAGATAAATTAATGGCCAAAGCACTGCTCTTATATGGTGTGATATTGAATCA
AAAAACAAGTGTGGCCTCACACCACTTTTGCTTGGCGTACATGAACAAAAACAGGAAGTGGTGAATTT
TTAATCAAGAAAAAGCTAATTTAAATGCACTTGATAGATATGGAAGAACTGCCCTCATACTTGCTGTA
TGTTGTGGATCAGCAAGTATAGTCAATCTTCTACTTGAGCAAAATGTTGATGTATCTTCTCAAGATCTA
TCTGGACAGACAGCCAGAGAGTATGCTGTTTCTAGTCATCATGTAAATTTGTGAATTACTTTCTGAC
TATAAAGAAAAACAGATGCTAAAAATCTTCTGAAAACAGCAATCCAGAACAGACTTAAAGCTGACA
TCAGAGGAAGAGTCACAAAGGCTTAAAGTCAGTGAAAAATAGCCAGCCAGAGAAAAATGTCTCAAGACCA
GAAATAAATAAGGACTGTGATAGAGAGGTTGAAGAAGAAATAAAGAAGCATGGAAGTAATCTGTGGGA
TTACCAGAAAACCTGACTAATGGTGCCAGTCTGGCAATGGTGTGATGATGATTAATTCACAAAGGAAG
AGCAGAAAACCTGAAAATCAGCAATTTCTGACACTGAGAATGAAGAGTATCACAGTGACGAACAAAAT
GATACCCAGAAACAACCTTTCTGAAGAACAGAACTGGAATATCACAAGATGAGATTCTGACTAATAAAA
CAAAAGCAGATAGAAGTGGCTGAAAAGGAAATGAATTCTGAGCTTTCTCTTAGTCATAAGAAAGAAGAA
GATCTCTTGCGTGAAAACAGCATGTTGCGGGAAGAAATTGCCAAGCTAAGACTGGAAGTATGAAACA
AAACATCAGAACAGCTAAGGGAAAATAAAATTTTGGAGGAAATTGAAAGTGTAAGAAAGAACTTCTA
AAGACTATACAACTGAATGAAGAAGCATTAAACGAAAACAGTATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001277304

Insert Size: 1635 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001277304.1

RefSeq Size: 1687 bp

RefSeq ORF: 1635 bp

Locus ID: 100996331

UniProt ID: H3BUK9

Cytogenetics: 15q11.2

MW: 61.7 kDa