

Product datasheet for **SC336718**

ZNF891 (NM_001277291) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF891 (NM_001277291) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF891
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336718 representing NM_001277291.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCAGTTATGGACCTATCCTCCCATGGGCTTTAACTAAACAGGACTCTGCCTGTTTCCATCTGAGA
AATGCTGAAGAGGAAAGGATGATTGCTGATTTCTGACAACCTGGTTACAGGAACCAATGACTTTCAAA
GATGTAGCTGTGGAGTTCACTCAGGAGGAGTGGATGATGCTGGATTCTGCTCAAAGAAAGTCTGTACAGA
GATGTGATGTTGAAAACTATAGAAATCTCACCTCCGTGGAATATCAGTTGTACAGGCTTACTGTGATC
TCCCATTTGGATCAAGAAGAGATAAGGAATATGAAAAAGAGAATTCCTCAAGCCATCTGTCCAGACCAG
AAGATTCAACCTAAAACCAAGAATCAACTGTGCAGAAAATTTTGTGGGAGGAACCATCCAATGCAGTG
AAAATGATAAACTCACAATGCATAATTGGTCTCCACATTAAGAGAAGACTGGGAATGCCATAAGATC
AGAAAACAGCATAAGATTCCAGGGGACATTGGAGGCAATGATATATGCCCCAAAGAAAACAGTACCT
CAGGAGCTATTCCGTGACTATCATGAATTAGAGGAAAACCTCTAACTTGGATCAAACTTATTTTTTCA
CAGAGCATTTTACCAGTAAACATTGCCAAAAATGTTACTCAGAGATTGGATGTTTGAAACACAATTCA
ATCATAAAACAATTATGTGAAAACTCTATAAGTGAGAAGCTCTATGAAAGTCATGAATGTGATACAACT
CTATGGCATTTCAGAGAAATCAACAGTACAAAAAGAGTACACATATTCTAAACATGGAATGCACTTC
ACACATAATATGTTTCTGTACCTAACAAATTTGCATATGGCACAGAATGCCTGTGAATGTAATAAGAT
GAAACGCTGTGTCATCAATCATCCCTTAAGAAACAAGGACAACTCATACTGAAAAGAAACATGAATGC
AATCAATGTGGAAGGCTTCAAAGGATTTCTAATCTTACTTTTACAAGAAAAGTACATGGGTGAG
AAACAATATGAATGTAAGAATGTGGTAAAGTGTCAATGATTCTCAACCTTAAGGAGACATGTAAGA
ACTCACACTGGAGAGAAACCTATGAATGTAATCAATGTGGAAGGCTTTCAGTCAAAAAACATCCCTT
AAGGCTCATGAGAACTCACACTGGAGAGAAACCTTATGAATGTAATCAGTGTGAAAAATCCTTTGGC
ACAAGCTCTTACCTTATAGTGCACAAGAGAATACACACTGGAGAAAACTCTATGAATGCAGTGAGTGT
GGAAAAGCCTTCAACACGAGCTCTCACCTTAAAGTTCATAAGAAAATTCATACCGAGAGAAATGTTTAT
GAATGCAGTGACTGTGGGAAAGTTTTCAGTGGGTCTCATCCCTAGAATGCATATAAGAACTCACACT
GGAGAGAAACCTATGAATGTAAGGAATGTAGGAAAGCCTTCAGTGTTCCTCTTCCCTTAGGAGGCAT
GTGAGAATTACACTGGAGAGAAACCTATGAATGTATTCAATGTGGAAGGCTTTCAGTCAGAGCTCT
TCACTTATTATACACAAGAGAATTCATACTGAGAGAGAAACCTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001277291

Insert Size: 1635 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001277291.1
RefSeq Size:	2303 bp
RefSeq ORF:	1635 bp
Locus ID:	101060200
UniProt ID:	A8MT65
Cytogenetics:	12q24.33
MW:	63.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]