

## Product datasheet for **SC336722**

### **ABLIM3 (NM\_001301028) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ABLIM3 (NM_001301028) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABLIM3
Synonyms:	HMFN1661
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336722 representing NM\_001301028.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACACTAGCATTCCTTATCAGCAGAATCCTTACAATCCACGGGGCAGCTCCAATGTCATCCAGTGC
TACCGCTGTGGAGACACCTGCAAAGGGGAAGTGGTCCGCGTGCACAACAACCACTCCACATCAGATGC
TTCACCTGTCAAGTATGTGGCTGTGGCCTGGCCAGTCAGGCTTCTTCTCAAGAACCAGGAGTACATC
TGCACCCAGGACTACCAGCAACTCTATGGCACCCGCTGTGACAGCTGCCGGGACTTCATCACAGGGCAA
GTCATCTCGGCCCTGGGCCGCACTTACCACCCCAAGTGCCTTCGTGTGACAGCTTGTGCAGGAAGCCTTTC
CCCATTGGAGACAAGGTGACCTTCAGCGGTAAGAATGTGTGCCAAACGTGCTCCAGTCCATGGCC
AGCAGTAAGCCCATCAAGATTCGTGGACCAAGCCACTGTCCGGGTGCAAGGAGGAGATCAAGCACGGC
CAGTCACTCCTGGCTCTGGACAAGCAGTGGCACGTGAGTCTCAAGTCCAGACTGCAGCGTCATC
CTCACCGGGGAGTATATCAGCAAGGATGGTGTCCATACTGTGAGTCCGACTACCATGCCAGTTTGGC
ATTAATGTGAGACTTGTGACCGATACATCAGTGGCAGAGTCTTGAGGGCAGGAGGAAGCACTACCAC
CCAACCTGTGCCAGTGTGTACGCTGCCACCAGATGTTACCCGAAGGAGAGGAAATGTACCTCACAGGT
TCCGAGGTTTGGCACCCCATCTGCAAACAGGCAGCCCGGGCAGAGAAGAAGTTAAAGCATAGACGGACA
TCTGAAACCTCCATCTCACCCCTGGATCCAGCATTGGGTCAACCAACCGAGTCACTGCGACATCTAC
GAGAACCTGGACCTCCGGCAGAGACGGGCTCCAGCCCGGGTACATAGACTCCCCACCTACAGCCGG
CAGGGCATGTCCCCACCTTCTCCCGCTCACCTCACCCTACTACCGCTCTGGTGATTTGTCTACAGCA
ACCAAGAGCAAAACAAGTGAAGACATCAGCCAGACCTCCAAGTACAGTCCCATCTACTCGCCAGACCCC
TACTATGCTTCGGAGTCTGAGTACTGGACTACCATGGGTCCCCAAAGTGCCCGAGCCAGAAGGTTTC
TCGTCTGGAGGAGAGGAGGATGATTTTGACCCGACATGCACAAGCTCCAAGTGGAAATGGCCGGCTG
ATTCTGAAGGAAGAAATGAAGGCCCGGTGAGCTCCTATGCAGATCCCTGGACCCCTCCCGGAGCTCC
ACCAGCAGCCGGGAAGCCCTGCACACAGCTGGCTATGAGATGTCCTCAATGGCTCCCTCGGTGCGAC
TACCTGGCTGACAGTGATCCTCTCATCTCAAATCTGCCTCCCTGCCTACCAGGAAGAAATGGGCTG
CACAGGACACCCAGCGCAGACCTTCCACTACGACAGCATGAACGCAGTCAACTGGGGCATGCGAGAG
TACAAGATCTACCTTATGAACTGCTGCTGGTACTACAAGAGGAAGAAACCGACTGCCCAAGGATGTA
GACAGGACCCGTTAGAGGGAACTTTTGAAGAGTGGCTGCTTATGA
ACGGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_001301028
<b>Insert Size:</b>	1635 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001301028.1</a>
<b>RefSeq Size:</b>	3131 bp
<b>RefSeq ORF:</b>	1635 bp
<b>Locus ID:</b>	22885
<b>UniProt ID:</b>	<a href="#">O94929</a>
<b>Cytogenetics:</b>	5q32
<b>Protein Pathways:</b>	Axon guidance
<b>MW:</b>	61.6 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the actin-binding LIM (abLIM) family of proteins. These proteins are characterized by an N-terminal LIM domain and a C-terminal dematin-like domain. The encoded protein interacts with actin filaments and may be a component of adherens junctions in several cell types. A variant of this gene may be associated with pain sensitivity in male human patients. [provided by RefSeq, Sep 2016]</p> <p>Transcript Variant: This variant (5) lacks several in-frame exons and has an alternate splice site in the 3' coding region, compared to variant 1. The resulting isoform (4) lacks two internal segments and has a shorter and distinct C-terminus, compared to isoform 1.</p>