

Product datasheet for **SC336727**

BTD (NM_001281723) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BTD (NM_001281723) Human Untagged Clone
Tag:	Tag Free
Symbol:	BTD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336727 representing NM_001281723.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGCCAGAGGGAGCGGGACTAGCAGGAGATTGCTGCCTATGCAAAGCAGATTTGTGGTCTGCATTATG
TCTGGAGCCAGAAGTAAGCTTGCTCTTTTCTCTGCGGCTGTTACGTGGTTGCCCTGGGAGCCACACC
GGGGAGGAGAGCGTGGCTGACCATCACGAGGCTGAATATTATGTGGCTGCCGTGTATGAGCATCCATCC
ATCCTGAGTCTGAACCTCTGGCTCTCATCAGCCGCCAAGAGGCCTTGGAGCTCATGAACCAGAACCTT
GACATCTATGAACAGCAAGTGACTGACGCCAAAAGGATGTACAGATTATAGTGTTCAGAAAGAT
GGCATTTCATGGATTCAACTTTACAAGAACATCCATTTATCCATTTTGGACTTCATGCCGCTCCCCAG
GTGGTCAGGTGGAACCCATGCCTGGAGCCTCACCGCTTCAATGACACAGAGGTGCTCCAGCGCTGAGT
TGTATGGCCATCAGGGGAGATAGTTCTTGGTGGCCAATCTTGGGACAAAGGAGCCTTGTCATAGCAGT
GACCAAGGTGCCAAAAGATGGGAGATACCAGTTC AACACAAAATGTCGTGTT CAGCAATAATGGAACC
CTTGTTGACCGCTACCGTAAACACAACCTCTACTTTGAGGCAGCATTGATGTTCTCTTAAAGTGGAT
CTCATCACCTTTGATACCCCTTTGCTGGCAGGTTTGGCATCTTCACATGCTTTGATATATTGTTCTTT
GACCCTGCCATCAGAGTCTCAGAGACTACAAGTGAAGCATGTTGTGTACCCAACCTGCCTGGATGAAC
CAGCTCCCCTCTTGGCAGCAATTGAGATT CAGAAAAGCTTTTGCTGTTGCCTTTGGCATCAACGTTCTG
GCAGCTAATGTCCACCACCCAGTTCTGGGGATGACAGGAAGTGGCATAACACCCCTCTGGAGTCCTTT
TGGTACCATGACATGGAAAATCCAAAAGTCACTTATAATTGCCAGGTGGCCAAAATCCAGTGGGT
CTCATTGGTGCAGAGAATGCAACAGGTGAAACGGACCCATCCCATAGTAAGTTTTAAAAATTTGTCA
GGCGATCCGTA CTGTGAGAAGGATGCTCAGGAAGTCCACTGTGATGAGGCCACCAAGTGGAAACGTGAAT
GCTCTCTCCACATTTCACTCTGAGATGATGATGACAATTTACCCTGGTCCCTGTCTGGGAAAAGGAA
GGCTATCTCCACGTCTGTTCCAATGGCCCTCTGCTGTTATTTACTTTACGAGAGGCCACCTTATCCAAA
GAGCTGTATGCCCTGGGGTCTTTGATGGGCTTACACAGTACATGGCACTTACTACATCCAAGTGTGT
GCCCTGTCAGGTGTGGGGTCTTGGCTTGCACACCTGTGGACAGGAAATCACAGAGGCCACGGGGATA
TTTGAGTTTACCTGTGGGGCACTTCACTACTTCTATATCTTTCTTTGTTTCTGACCTCAGGGATG
ACCCTAGAAGTCCCTGACCAGCTTGGCTGGGAGAATGACCACTATTTCTGAGGAAAAGTAGGCTGTCC
TCTGGGCTGGTGACGGCGCTCTATGGGCGCTTGTATGAGAGGGACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN:	NM_001281723
Insert Size:	1638 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001281723.2
RefSeq Size:	4169 bp
RefSeq ORF:	1638 bp
Locus ID:	686
UniProt ID:	P43251
Cytogenetics:	3p25.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Biotin metabolism, Metabolic pathways
MW:	61.3 kDa
Gene Summary:	<p>The protein encoded by this gene functions to recycle protein-bound biotin by cleaving biocytin (biotin-epsilon-lysine), a normal product of carboxylase degradation, resulting in regeneration of free biotin. The encoded protein has also been shown to have biotinyl transferase activity. Mutations in this gene are associated with biotinidase deficiency. Multiple transcript variants encoding different isoforms have been described. [provided by RefSeq, Aug 2013]</p> <p>Transcript Variant: This variant (1) encodes isoform 1. Variants 1-4 and 6 encode the same isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>