

Product datasheet for **SC336735**

TAP1 (NM_001292022) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAP1 (NM_001292022) Human Untagged Clone
Tag:	Tag Free
Symbol:	TAP1
Synonyms:	ABC17; ABCB2; APT1; D6S114E; PSF-1; PSF1; RING4; TAP1*0102N; TAP1N
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336735 representing NM_001292022.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCATTCCATTCTTTACGGGCCGCCTCACTGACTGGATTCTACAAGATGGCTCAGCCGATACCTTC
ACTCGAAACTTAACTCTCATGTCCATTCTCACCATAGCCAGTGCAGTGCCTGGAGTTCGTGGGTGACGGG
ATCTATAACAACACCATGGGCCACGTGCACAGCCACTTGCAGGGAGAGGTGTTTGGGGCTGCTCTGCGC
CAGGAGACGGAGTTTTTCCAACAGAACCAGACAGGTAACATCATGTCTCGGGTAACAGAGGACACGTCC
ACCCTGAGTGATTCTCTGAGTGAGAATCTGAGCTTATTCTGTGGTACCTGGTGCAGGCCTATGTCTC
TTGGGGATCATGCTCTGGGGATCAGTGTCCCTCACCATGGTCACCCTGATCACCTGCCTCTGCTTTTC
CTTCTGCCAAGAAGGTGGGAAATGTTACCAGTTGCTGGAAGTGCAGGTGCGGGAATCTCTGGCAAAG
TCCAGCCAGGTGGCCATTGAGGCTCTGTCGGCCATGCCTACAGTTCGAAGCTTTGCCAACGAGGAGGGC
GAAGCCAGAAGTTTAGGAAAAGCTGCAAGAAATAAGACACTCAACCAGAAGGAGGCTGTGCCCTAT
GCAGTCAACTCCTGGACCACTAGTATTTCAAGTATGCTGCTGAAAGTGGGAATCCTCTACATTGGTGGG
CAGCTGGTGACCAAGTGGGGCTGTAAAGCAGTGGGAACCTTGTACATTTGTCTCTACCAGATGCAGTTC
ACCCAGGCTGTGGAGGACTGCTCTCCATCTACCCAGAGTACAGAAGGCTGTGGGCTCCTCAGAGAAA
ATATTTGAGTACCTGGACCGCACCCCTCGCTGCCACCCAGTGGTCTGTTGACTCCCTTACACTTGAG
GGCCTTGTCCAGTCCAAGATGTCTCCTTTGCCTACCCAAACCGCCAGATGTCTTAGTGCTACAGGGG
CTGACATTCACCCACGCCCTGGCGAGGTGACGGCGCTGGTGGGACCAATGGGTCTGGGAAGAGCACA
GTGGCTGCCCTGCTGCAGAATCTGTACCAGCCACCGGGGACAGCTGCTGTTGGATGGGAAGCCCTT
CCCCAATATGAGCACCCTACCTGCACAGGCAGGTGGCTGCAGTGGGACAAGAGCCACAGGTATTTGGA
AGAAGCTTCAAGAAAATATTGCCTATGGCCTGACCCAGAAGCCAATATGGAGGAAATCACAGCTGCT
GCAGTAAAGTCTGGGGCCATAGTTTTCATCTCTGGACTCCCTCAGGGCTATGACACAGAGGTAGACGAG
GCTGGGAGCCAGCTGTCAGGGGGTCAAGGACAGGCAGTGGCGTTGGCCCGAGCATTGATCCGGAAACCG
TGTGTAATTCCTGGATGATGCCACCAGTGCCTGGATGCAAACAGCCAGTTACAGGTGGAGCAGCTC
CTGTACGAAAGCCCTGAGCGGTAATCCCCTCAGTGTCTCTCATACCCAGCACCTCAGCCTGGTGGAG
CAGGCTGACCACATCCTCTTTCTGGAAGGAGGCGCTATCCGGGAGGGGGAAACCACCAGCAGCTCATG
GAGAAAAGGGGTGCTACTGGGCCATGGTGCAGGCTCCTGCAGATGCTCCAGAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001292022
Insert Size:	1644 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001292022.1
RefSeq Size:	2234 bp
RefSeq ORF:	1644 bp
Locus ID:	6890
UniProt ID:	Q03518
Cytogenetics:	6p21.32
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	ABC transporters, Antigen processing and presentation, Primary immunodeficiency
MW:	60 kDa
Gene Summary:	<p>The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. The protein encoded by this gene is involved in the pumping of degraded cytosolic peptides across the endoplasmic reticulum into the membrane-bound compartment where class I molecules assemble. Mutations in this gene may be associated with ankylosing spondylitis, insulin-dependent diabetes mellitus, and celiac disease. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2014]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (2) is shorter at the N-terminus compared to isoform 1.</p>