

Product datasheet for **SC336736**

C17orf80 (NM_001288771) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C17orf80 (NM_001288771) Human Untagged Clone
Tag:	Tag Free
Symbol:	C17orf80
Synonyms:	HLC-8; MIG3; SPEP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336736 representing NM_001288771.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTGATAATCCACCAGAATGGAAGTGTGCCTTACTGTAGAAGCCATTTAAACGATTAATAATCC
CACTTGCCATACTGTAGATGATAGGACCAACCATACTACTGATCAAAAAGTTTATCAGTCCAAGCCA
GCTACACTCCCACGTGCTAAAAAGATGAAAGGACCAATCAAAGATTTAATTAAAGCTAAAGGGAAAGAG
TTAGAGACAGAGAATGAAGAAAAGAAATTCTAAGTTGGTGGTGGACAAACCAGAACAGACAGTGAAGACC
TTTCCACTGCCAGCTGTTGGTTTGGAAAGAGCAGCTACTACAAAGGCAGATAAAGACATCAAGAATCCA
ATCCAACCATCCTTCAAATGTTAAAAATACTAAACCAATGACTACTTTCCAAGAAGAAACCAAGGCT
CAGTTTTACGCATCAGAGAAAACCTCTCCTAAAAGAGAACTTGCCAAAGATTTGCCTAAATCAGGAGAA
AGTCGATGTAATCCTTCAGAAGCTGGAGCGTCTTACTGGTTGGCTCAATAGAACCTTCTTGTCAAT
CAAGATAGAAAATATTCCTCAACTCTACCTAATGATGTACAACTACCTCTGGTGTCTCAAATTTGGAC
AAAATTGATCCCCAAAGACAGGAACCTTCTAGTAAAATTACTAGATGTGCCTACTGGTGTGTCATATT
TCTCCAAAGAATGTCAGTGATGGGGTAAAAGGGTAAGAACATTATTAAGCAATGAGAGAGATTCCTAAA
GGCAGGGATCACCTCTCAGGAGTCCCTACTGATGTTACAGTTACTGAGACTCCAGAAAAGAACACAGAA
TCCCTCATTTTAAAGCCTTAAAATGAGCTCATTAGGTAAAATCCAAGTCATGGAGAAAACAAGAGAAAGGA
CTTACCCTGGGAGTAGAGACGTGTGGGAGCAAAGGAAATGCAGAGAAAAGTATGTCTGCAACAGAAAAG
CAGGAACGGACTGTCATGAGCCATGGCTGTGAGAACTTCAACACCAGGGATTTCAGTCACAGGAAAGGAG
TCTCAAGGGGAAAGACCACATTTAAGTTTGTTCATTCCGAGGGAGACGACTTACCAGTTTCATTCTGTA
TCGCAGTCAAGTAGTCAAAGTCTTGCCTCTAGCTACAACATTTCTTCAAGAAAAGAAAGCAGAAGCC
CAGAATCATAATTGTGTCCCTGATGTAAGGCATTAATGGAGAGTCCCAGGGACAGTTATCTCTGGAG
CCCAAATCTGATAGTCAGTTCCAAGCATCACACTGGGTGCCAGAGCCCTTTATGTTTCAGCCCAGCGT
CACACTCCTCAGAGCCCTTACCAATCATGCTGACGCTGCTGGCAGGAAGACTTTCGCAGCTGCATG
GGGCTGGAGTGGTTTCCAGAGCTCTATCCTGGTTACCTTGGACTAGGGGTGTTGCCAGGGAAGCCTCAG
TGTTGGAATGCAATGACCCAGAAGCCACAACCTTATCAGTCCCCAGGGGAAAGACTCTCACAGGCTGG
ATCAGGTGCAACACCACCATAAGGAAGAGTGGATTCGGTGGCATCAGGATGCTTTCACAGGATACTTC
GTCCTGTGTTGTAGCTGGAGTTTCAGACGCTCTGAAGCTGCAGCGGTGGCGTAAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001288771
Insert Size:	1644 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001288771.1</u>
RefSeq Size:	2168 bp
RefSeq ORF:	1644 bp
Locus ID:	55028
Cytogenetics:	17q25.1
Protein Families:	Transmembrane
MW:	60.8 kDa