

Product datasheet for **SC336743**

RNF214 (NM_001278249) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RNF214 (NM_001278249) Human Untagged Clone
Tag:	Tag Free
Symbol:	RNF214
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336743 representing NM_001278249.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGCGTCTGAGGTGTGTTGGTGTGGCCAATGCCCCAGTCTCCGGAATCTTCTAGTTTATGT
GCTTCCAATCAGACGAAGGTCTCCAGATGGTCTAAGCACAAAGACTCTGCACAGAAGCAGAAGAAC
TCGCCTCTGTTGAGTACTGACTTTAAGACAGCTGATTCAGAGGTAACACAGATCAAGATATTGAAAAG
AATTTGGATAAAATGATGACAGAGAGAACCCTGTTGAAAGAGCGTTACCAGGAGTCTGGACAAACAG
AGGCAAGTGGAGAATCAGCTCCAAGTGAATTAAGCAGCTTCAGCAAAGGAGAGAAGAGGAAATGAAG
AATCACCAGGAGATATTAAGGCTATTCAGGATGTGACAAATAAGCGGGAAGAAACAAAGAAGAAGATA
GAGAAAGAGAAGAAGGAGTTTTTGCAGAAGGAGCAGGATCTGAAAGCTGAAATTGAGAAGCTTTGTGAG
AAGGGCAGAAGAGAGGTGTGGGAAATGGAAGTGGATAGACTCAAGAATCAGGATGGCGAAATAAATAGG
AACATTATGGAAGAGACTGAACGGGCTGGAAGGCAGAGATCTTACTAGAGAGCCGGAAGAGTTA
CTGGTACTGAAACTAGAAGAAGCAGAAAAGAGGCAGAATTGCACCTTACTTACCTCAAGTCAACTCCC
CCAACACTGGAGACAGTTCGTTCCAACAGGAGTGGGAGACGAGACTGAAATGGAGTTCGGATAATGAAA
AAGAATGTTCTGTGACCAATTTAATAGTCATATCCAGTTAGTGAGGAACGGAGCCAAGCTGAGCAGCCTT
CCTCAATCCCTACTCCCACTTTACCTCCACCCCATCAGAGACAGACTTCATGCTTCAGGTGTTTCAA
CCCAGTCCCTCTCTGGCTCCTCGGATGCCCTTCTCCATTGGGCAGGTCACAATGCCCATGGTTATGCC
AGTGCAGATCCCCGCTCCTTGTCTTTCCAATCCTGAACCCTGCCCTTCCCAGCCAGCCAGCCTTCC
TCACCCCTTCTGGCTCCCATGGCAGAAATAGCCCTGGCTGGGTTCCCTTGTGACCCCCACGGTCCA
CACATGCCCCCTGCCCTCCATCCCACCTCCCCAGGCTTGGGGGTGTTAAGGCTTCTGCTGAACT
CCCCGGCCCCAACCAAGTAGACAAACTGGAGAAGATCCTGGAGAAGCTGCTGACCCGGTCCCACAGTGC
AATAAGGCCCAGATGACCAACATTCTTCAGCAGATCAAGACAGCACGTACCACCATGGCAGGCCTGACC
ATGGAGGAACCTTATCCAGTTGGTTGCTGCACGACTGGCAGAACATGAGCGGGTGGCAGCAAGTACTCAG
CCACTTGGTCGCATCCGGGCTTGTTCCTGCTCCACTGGCCCAATCAGTACCCCAATGTTCTTGCCT
TCTGCCAAGTTTCATATCCTGGAAGGTCTTACATGCTCCAGCCACCTGTAAGCTATGTCTAATGTGC
CAGAAACTCGTCCAGCCAGTGGCTGCATCCAATGGCGTGTACCCATGTATTGCACAAGGAGTGTATC
AAATTCTGGGCCAGACCAACACAAATGACACTGTCCCTTTGTCCAACCTTAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN:	NM_001278249
Insert Size:	1647 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001278249.1</u>
RefSeq Size:	2245 bp
RefSeq ORF:	1647 bp
Locus ID:	257160
UniProt ID:	<u>Q8ND24</u>
Cytogenetics:	11q23.3
MW:	61.3 kDa