

Product datasheet for **SC336744**

ZNF18 (NM_001303281) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF18 (NM_001303281) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF18
Synonyms: HDSG1; KOX11; Zfp535; ZKSCAN6; ZNF535; ZSCAN38
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC336744 representing NM_001303281.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCCGTTGACTTGGGGCAGGCCCTAGGCCTGCTGCCATCGCTGGCGAAGGCCGAGGACTCCCAGTTCC
TCAGAATCAGATGCTGCCCTTCAAGAGGAAGCTCTCCAGCCCTGAGACCGCACGCCAGCTTTTCAGGCAG
TTCCGTTACCAGGTGATGTCTGGGCCTCATGAGACCTTGAAGCAACTTCGGAAGCTCTGTTTCCAGTGG
CTACAGCCAGAGGTTACACCAAAGAGCAGATCCTAGAGATCCTCATGTTGGAGCAGTTTCTGACCATC
CTGCCTGGGAGATCCAGATGTGGGTGCGGAAACAGTGTCCAGGAAGTGGAGAAGAGGCAGTGACCCTT
GTGGAAAGCTTGAAGGGGACCCCAAGAGACTGTGGCAATGGATCAGTATCCAGTTCTAGGACAGGAC
ATCTTATCAGAGAAGATGGAATCTCAAGCTGCCAAGTGGGGAAAGTGGAGCCCCATCTTGAAGTGGTG
CCTCAGGAGTTGGGACTTGAGAATTCATCCTCAGGGCCTGGGGAGCTTCTGAGCCACATCGTGAAGAG
GAATCTGACACAGAAGCAGAAGTACCCCTGGCTGCCTCCCAGCCTGCCCGACTGGAGGAAAGGCTGATC
AGAGACCAGGACCTCGGAGCCTCACTGCTCCCAGCAGCACCTCAGGAACAGTGGAGACAAGTGGATTCC
ACTCAAAGGAGCAATACTGGGATCTCATGCTGGAGACCTATGGGAAAATGGTCTCAGGAGCAGGCATT
TCCCATCCCAAATCTGACCTGACTAATTAATAGAATTTGGGGAAGAGCTGGCAGGAATATACCTTCAT
GTCAATGAGAAGATCCCAAGACCACCTGCATAGGAGATAGACAAGAGAATGACAAGGAGAACCTAAAT
TTGGAGAATCACAGGGACCAGGAGCTCCTGCATGCTTCTGCAAGCTTCAGGAGAGTTTCTTCTCAG
GCTTCTTGGAGGGCTTCTTCACTGAGGATGAGCCAGGATGCTTTGGAGAAGGAGAGAATCTCCCTGAG
GCTCTGCAAAAACATTCAGGATGAGGGAACAGGGGAACAGCTGTCTCCTCAAGAAAGGATTCTGAGAAA
CAACTAGGTCAGCATTTGCCTAATCCTCATTAGGAGAAAATGTCCACCATGTGGCTTGGAGAGAAGAGA
GAGACCTCCCAGAAGGGGACGCAAGAGCCCCATGGCCAGAAGCTCCCACCTGCAGGGAGTGTGGG
AAGACCTTTTATAGGAATTCAGCTTATTTTTACCAAAGAAGTACACCCGGAGAGACATACTTTCAG
TGCACCATCTGCAAAAAGCCTTCTGCGGAGTTCAGACTTTGTGAAGCATCAGAGAAGTACACCGGA
GAGAAGCCCTGTAATGTGATTACTGTGGAAAGGCTTTAGTGACTTCTCAGGATTGCGCCACCACGAG
AAAATCCACACAGGAGAGAAACCTATAAATGTCTATCTGTGAGAAAAGTTTCATTAGAGATCAAAC
TTAATAGACATCAGAGGTTTCACTGGAGAGAAACCTATAAATGTTGCACTGTGGGAAAAGTTTC
AGCTGGAGCTCGAGCCTTGACAAACATCAAAGATCCCCTTAGGAAAGAAGCCCTTTCAATAG
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001303281



[View online »](#)

Insert Size:	1650 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001303281.1</u>
RefSeq Size:	2344 bp
RefSeq ORF:	1650 bp
Locus ID:	7566
UniProt ID:	<u>P17022</u>
Cytogenetics:	17p12
Protein Families:	Transcription Factors
MW:	62.3 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode isoform 1.