

Product datasheet for **SC336745**

TBC1D3G (NM_001291462) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBC1D3G (NM_001291462) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBC1D3G
Synonyms:	TBC1D3C; TBC1D3D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336745 representing NM_001291462.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGACGTGGTAGAGGTCGCGGGCAGTTGGTGGGCACAAGAGCGAGAGGACATCATTATGAAATACGAA
AAGGGACACCGAGCTGGGCTGCCAGAGGACAAGGGGCCTAAGCCTTTTCGAAGCTACAACAACAACGTC
GATCATTGGGGATTGTACATGAGACGGAGCTGCCTCCTCTGACTGCGCGGGAGGCGAAGCAAATTCGG
CGGGAGATCAGCCGAAAGAGCAAGTGGGTGGATATGCTGGGAGACTGGGAGAAATACAAAAGCAGCAGA
AAGCTCATAGATCGAGCGTACAAGGGAATGCCCATGAACATCCGGGGCCCGATGTGGTCAGTCCTCTG
AACATTGAGGAAATGAAGTTGAAAAACCCCGGAAGATACCAGATCATGAAGGAGAAGGGCAAGAGGTCA
TCTGAGCACATCCAGCGCATCGACCGGGACATAAGCGGGACATTAAGGAAGCATATGTTCTTCAGGGAT
CGATACGGAACCAAGCAGCGGAACTACTCCACATCCTCTGGCATATGAGGAGTATAACCCGGAGGTG
GGCTACTGCAGGGACCTGAGCCACATCGCCGCCTTGTTCTCTCTATCTTCTGAGGAGGATGCATTC
TGGGCACTGGTGCAGCTGCTGGCCAGTGAGAGGCACTCCCTGCAGGGATTTACAGCCCAAATGGCGGG
ACCGTCCAGGGGCTCCAAGACCAACAGGAGCATGTGGTAGCCACGTCACAACCAAGACCATGGGGCAT
CAGGACAAGAAAGATCTATGTGGGCAGTGTCCCGTTAGGCTGCCTCATCCGGATATTGATTGACGGG
ATCTCTCTCGGGCTCACCTGCGCTGTGGGACGTGTATCTGGTAGAAGGCGAACAGGCGTTGATGCCG
ATAACAAGAAATCGCCTTTAAGGTTACAGAGAAGCGCCTCACGAAGACGTCCAGGTGTGGCCCGTGGGCA
CGTTTTTGCAACCGGTTCTGTGATACCTGGGCCAGGGATGAGGACACTGTGCTCAAGCATCTTAGGGCC
TCTATGAAGAACTAACAAGAAAGCAGGGGGACCTGCCACCCCGAGCCAAACCCGAGCAAGGGTCGTG
GCATCCAGGCCTGTGCCGGCTTCACGTGGCGGGAAGACCCTCTGCAAGGGGGACAGGCAGGCCCTCCA
GGCCACACAGCCCGTTCCCGGCCCATTTGGTCAGTTCCCGGCCACGGGCACCTCGTTCTTCCACA
CCCTGTCTGTGGGGCTGTCCGGGAAGACACCTACCCTGTGGGCACTCAGGGGTGTGCCAGCCCGGCC
CTGGCTCAGGGAGGACCTCAGGGTTCCTGGAGATTCTGCACTGGAACCTCATGCCCGCCTCCCAACG
GACCTGGACGTAGAGGGCCCTTGTTCCGCCATTATGATTTAGACAGAGCTGCTGGGTCCGTGCCATA
TCCAGGAGGACAGCTGGCCCTGCTGGCAGGCTGAACACCCTGCGGAGCGGGTGAGATCGGCTTTC
GCTGCACCCAGCACTGATTCCGACAGGGCACCCCTTCAGAGCTAGGGACGAACAGCAGTGTGCTCCC
ACCTCAGGGCCTTGCTCTGCGGCCTCCACTTGAAAGTTCTCAGTTCCCTCCAGGCTTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001291462

Insert Size: 1650 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291462.1</u>
RefSeq Size:	2090 bp
RefSeq ORF:	1650 bp
Locus ID:	101060321
UniProt ID:	<u>Q6IPX1</u>
Cytogenetics:	17q12
MW:	62.2 kDa
Gene Summary:	Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11 (By similarity). [UniProtKB/Swiss-Prot Function]