

Product datasheet for **SC336749**

TBC1D3E (NM_001291466) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBC1D3E (NM_001291466) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBC1D3E
Synonyms:	PRC17; TBC1D3; TBC1D3A; TBC1DE3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336749 representing NM_001291466.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGACGTGGTAGAGGTCGCGGGCAGTTGGTGGGCACAAGAGCGAGAGGACATCATTATGAAATACGAA
AAGGGACACCGAGCTGGGCTGCCAGAGGACAAGGGGCCTAAGCCTTTTCGAAGCTACAACAACAACGTC
GATCATTGGGGATTGTACATGAGACGGAGCTGCCTCCTCTGACTGCGCGGGAGGCGAAGCAAATTCGG
CGGGAGATCAGCCGAAAGAGCAAGTGGGTGGATATGCTGGGAGACTGGGAGAAATACAAAAGCAGCAGA
AAGCTCATAGATCAAGCGTACAAGGGAATGCCCATGAACATCCGGGGCCCGATGTGGTCAGTCCTCTG
AACACTGAGGAAATGAAGTTGAAAAACCCGGAAGATACCAGATCATGAAGGAGAAGGGCAAGAAGTCA
TCTGAGCACATCCAGCGCATCGACCGGGACGTAAGCGGGACATTAAGGAAGCATATATTCTTCAGGGAT
CGATACGGAACCAAGCAGCGGAACTACTCCACATCCTCTGGCATATGAGGAGTACAACCCGGAGGTG
GGCTACTGCAGGGACCTGAGCCACATCGCCGCCTTGTTCTCTCTATCTTCTGAGGAGGATGCATTC
TGGGCACTGGTGCAGCTGCTGGCCAGTGAGAGGCACTCCCTGCAGGGATTCACAGCCCAAATGGCGGG
ACCGTCCAGGGGCTCCAAGACCAACAGGAGCATGTGGTAGCCACGTCACAACCAAGACCATGGGGCAT
CAGGACAAGAAAGATCTATGTGGGCACTGTTCCCGTTAGGCTGCCTCATCCGGATATTGATTGACGGG
ATCTCTCTCGGGCTCACCTGCGCTGTGGGACGTGTATCTGGTAGAAGGCGAACAGGCGCTGATGCCG
ATAACAAGAAATCGCCTTTAAGGTTACAGAGAAGCGCCTCACGAAGACGTCCAGGTGTGGCCCGTGGGCA
CGTTTTTGCAACCGGTTCTGTGATACCTGGGCCAGGGATGAGGACACTGTGCTCAAGCATCTTAGGGCC
TCTATGAAGAACTAACAAGAAAGAAAGGGGACCTGCCACCCCGAGCAAACCCGAGCAAGGGTCGTG
GCATCCAGGCCTGTGCCGGCTTCACGTGGCGGGAAGACCCTCTGCAAGGGGGACAGGCAGGCCCTCCA
GGCCACCAAGCCCGTTCCCGGCCCATTTGGTCAGTTCCCGGCCACGGGCACCTCGTTCTTCCACA
CCCTGTCTGTTGGGGCTGTCCGGGAAGACACCTACCCTGTGGCACTCAGGGGTGTGCCAGCCCGGCC
CTGGCTCAGGGAGGACCTCAGGGTTCCTGGAGATTCTGCAGTGGAACCTCATGCCCGCCTCCCAACG
GACCTGGACGTAGAGGGCCCTTGTTCCGCCATTATGATTTAGACAGAGCTGCTGGGTCCGTGCCATA
TCCAGGAGGACAGCTGGCCCTGCTGGCAGGCTGAACACCCTGCGGAGCGGGTGAGATCGGCTTTC
GCTGCACCCAGCACTGATTCCGACAGGGCACCCCTTCAGAGCTAGGGACGAACAGCAGTGTGCTCCC
ACCTCAGGGCCTTGCTCTGCGGCTCCACTTGGAAGTTCTCAGTTCCCTCCAGGCTTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001291466

Insert Size: 1650 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291466.1</u>
RefSeq Size:	2063 bp
RefSeq ORF:	1650 bp
Locus ID:	102723859
UniProt ID:	<u>A0A087X179</u>
Cytogenetics:	17q12
MW:	62.1 kDa
Gene Summary:	Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11 (By similarity). [UniProtKB/Swiss-Prot Function]