

## Product datasheet for **SC336753**

### TTC39A (NM\_001297664) Human Untagged Clone

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids                        |
| Product Name:             | TTC39A (NM_001297664) Human Untagged Clone |
| Tag:                      | Tag Free                                   |
| Symbol:                   | TTC39A                                     |
| Synonyms:                 | C1orf34; DEME-6                            |
| Mammalian Cell Selection: | Neomycin                                   |
| Vector:                   | pCMV6-Entry (PS100001)                     |
| E. coli Selection:        | Kanamycin (25 ug/mL)                       |



[View online »](#)

**Fully Sequenced ORF:** >SC336753 representing NM\_001297664.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCGCCCTGGACCTTTCCTCACCAACCGTTCTCAGAAGCACTCAGCTACCTCAAGCCGAGAACC
AAGGAAAGCATGTACCACCTCACTGACATATGCCACCATCCTGGAGATGCAGGCCATGATGACCTTTGAC
CCTCAGGACATCCTGCTTGCCGGCAACATGATGAAGGAGGCACAGATGCTGTGTGACAGGGCACCGGAGG
AAGTCTTCTGTAACAGATTCTTTCAGCAGCCTGGTGAACCGCCACGCTGGGCCAATTCCTGAAAGAG
GAAATCCACGCTGAGGTCTGCTATGCAGAGTGCCTGCTGCAGCGAGCAGCCCTGACCTTCTGCAGGAC
GAGAACATGGTGAGCTTCATCAAAGGGCGCATCAAAGTTCGAAACAGCTACCAGACCTACAAGGAGCTG
GACAGCCTTGTTCAGTCTCACAACTGCAAGGGTGAGAACCACCCGCACTTTGAAGGAGGAGTGAAG
CTTGGTGTAGGGCCCTTCAACCTGACACTGTCCATGCTTCTACTAGGATCCTGAGGCTGTTGGAGTTT
GTGGGTTTTTCAGGAAACAAGGACTATGGGCTGCTGCAGCTGGAGAGGGAGCGTCAGGGCACAGCTTC
CGCTCTGTGCTCTGTGTCATGCTCCTGCTGTGCTACCACACCTTCTCACCTTCGTGCTCGGTACTGGG
AACGTCAACATCGAGGAGCCGAGAAGCTTTGAAGCCCTACCTGAACCGGTACCCTAAGGGTGCCATC
TTCTGTCTTTGCAGGGAGGATTGAAGTCATTAAGGCAACATTGATGCAGCCATCCGGCGTTTCGAG
GAGTGTGTGAGGCCAGCAGCACTGGAAGCAGTTCACCCACATGTGCTACTGGGAGCTGATGTGGTGC
TTCACCTACAAGGGCCAGTGAAGATGTCTACTTCTACGCCGACCTGCTCAGCAAGGAGAACTGCTGG
TCCAAGGCCACCTACATTTACATGAAGGCCGCTACCTCAGCATGTTTGGGAAGGAGGACCACAAGCCG
TTCGGGGACGACGAAGTGAATATTTTCGAGCTGTGCCAGGCCTGAAGCTCAAGATTGCTGGGAAATCT
CTACCCACAGAGAAGTTGCCATCCGGAAGTCCCGGCGCTACTTCTCCTCAACCCTATCTCGCTGCCA
GTGCTGCTCTGGAATGATGTACATCTGGAACGGCTACGCCGTGATTGGGAAGCAGCCGAAACTCAGC
GATGGGATACTTGAGATTATCACTAAGGCTGAAGAGATGCTGGAGAAAGGCCAGAGAACGAGTACTCA
GTGGATGACGAGTGTGGTGAATTTGTTGAAAGCCTGTGCTGAAATACCTGGGCCGTGTCCAGGAG
GCCGAGGAGAATTTAGGAGCATCTCTGCAATGAAAAGAAGATTAATATGACCACTACTTGATCCCA
AACGCCCTGTGGAGCTGGCCCTGCTGCTTATGGAGCAAGACAGAAACGAAGAGGCCATCAAACCTTTG
GAATCTGCCAAGCAAACACTACAAGAATTACTCCATGGAGTCAAGGACACACTTTCGAATCCAGGCAGCC
ACACTCCAAGCCAAGTCTTCCCTAGAGAACAGCAGCAGATCCATGGTCTCATCAGTGTCTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



|                        |   |
|------------------------|---|
| ACCN:                  | NM_001297664  |
| Insert Size:           | 1653 bp   |
| OTI Disclaimer:        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| Components:            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| Reconstitution Method: | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| RefSeq:                | <u><a href="#">NM_001297664.1</a></u>   |
| RefSeq Size:           | 2774 bp   |
| RefSeq ORF:            | 1653 bp   |
| Locus ID:              | 22996   |
| UniProt ID:            | <u><a href="#">Q5SRH9</a></u>   |
| Cytogenetics:          | 1p32.3  |
| MW:                    | 63.2 kDa  |