

Product datasheet for **SC336757**

PTGS1 (NM_001271164) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTGS1 (NM_001271164) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTGS1
Synonyms:	COX1; COX3; PCOX1; PES-1; PGG/HS; PGHS-1; PGHS1; PHS1; PTGHS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336757 representing NM_001271164.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGAGCCGGAGTCTTTGCTCTGGTCTTGTGTTCTCTGCTCCTGCTCCCGCCGCTCCCCGTCTGCTC
GCGGACCCAGGGCGGCCACGCCAGTGAATCCCTGTTGTTACTATCCATGCCAGCACCAGGGCATCTGT
GTCCGCTTCGGCCTTGACCGCTACCCAGTGTGACTGCACCCGCACGGGCTATTCCGGCCCCAACTGCACC
ATCCCTGGCCTGTGGACCTGGCTCCGGAATCACTGCGGCCAGCCCTCTTTACCCCACTTCTGCTC
ACTCACGGGCGTGGTTCTGGGAGTTTGTCAATGCCACCTTCATCCGAGAGATGCTCATGCGCTGGTA
CTCACAGGAAGAAGCAGTTGCCAGATGCCAGCTCCTGGCCCGCGCTTCTGCTCAGGAGGAAGTTC
ATACCTGACCCCCAAGGCACCAACCTCATGTTTGCCTTCTTGCACAACACTTACCCACCAGTTCTTC
AAAATCTCGGAAGATGGTCTGGCTTACCAAGGCCTTGGGCCATGGGGTAGACCTCGGCCACATT
TATGGAGACAATCTGGAGCGTCAGTATCAACTGCGGCTCTTAAGGATGGGAACTCAAGTACCAGGTG
CTGGATGGAGAAATGTACCCGCCCTCGGTAGAAGAGGCGCTGTGTTGATGCACTACCCCGAGGCATC
CCGCCCAAGCCAGATGGCTGTGGGCCAGGAGGTGTTGGGCTGCTTCCGCGCTCATGCTGTATGCC
ACGCTCTGGCTACGTGAGCACAACCGTGTGTGACCTGCTGAAGGCTGAGCACCCACCTGGGGCGAT
GAGCAGCTTTTCCAGACGACCCGCTCATCCTCATAGGGGAGACCATCAAGATTGTCATCGAGGAGTAC
GTGACGAGCTGAGTGGCTATTTCTGCAGCTGAAATTTGACCCAGAGCTGCTGTTCCGGTGTCCAGTTC
CAATACCGCAACCGCATTGCCATGGAGTTCAACCATCTCTACCACTGGCACCCCTCATGCCTGACTCC
TTCAAGGTGGGCTCCCAGGAGTACAGCTACGAGCAGTTCTTGTTCACACCTCCATGTTGGTGGACTAT
GGGGTTGAGGCCCTGGTGGATGCCTTCTCTGCCAGATTGCTGGCCGGATCGGTGGGGCAGGAACATG
GACCACCACATCCTGCATGTGGCTGGATGTCATCAGGGAGTCTCGGGAGATGCGGCTGCAGCCCTTC
AATGAGTACCGCAAGAGTTTGGCATGAAACCCTACACCTCCTTCCAGGAGCTCGTAGGAGAGAAGGAG
ATGGCAGCAGAGTTGGAGGAATTGTATGGAGACATTGATGCGTTGGAGTTCTACCCTGGACTGCTTCTT
GAAAAGTGCCATCCAAACTCTATCTTTGGGGAGAGTATGATAGAGATTGGGGCTCCCTTTTCCCTCAAG
GGTCTCTAGGGAATCCCATCTGTTCTCCGGAGTACTGGAAGCCGAGCACATTTGGCGGCGAGGTGGGC
TTTAACATTGTCAAGACGGCCCACTGAAGAAGCTGGTCTGCCTCAACACCAAGACCTGTCCCTACGTT
TCCTTCCGTGTGCCGGATGCCAGTCAGGATGATGGCCTGCTGTGGAGCGACCATCCACAGAGCTCTGA

ACCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001271164

Insert Size: 1656 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271164.1](#)

RefSeq Size: 4901 bp

RefSeq ORF: 1656 bp

Locus ID: 5742

Cytogenetics: 9q33.2

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Arachidonic acid metabolism, Metabolic pathways

MW: 63.3 kDa

Gene Summary:

This is one of two genes encoding similar enzymes that catalyze the conversion of arachinodate to prostaglandin. The encoded protein regulates angiogenesis in endothelial cells, and is inhibited by nonsteroidal anti-inflammatory drugs such as aspirin. Based on its ability to function as both a cyclooxygenase and as a peroxidase, the encoded protein has been identified as a moonlighting protein. The protein may promote cell proliferation during tumor progression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Transcript Variant: This variant (3) lacks an alternate exon in the coding region, compared to isoform 1. The encoded isoform (3) is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.