

Product datasheet for **SC336777**

ETFDH (NM_001281738) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ETFDH (NM_001281738) Human Untagged Clone
Tag:	Tag Free
Symbol:	ETFDH
Synonyms:	ETFQO; MADD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336777 representing NM_001281738.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAAGGTTTGCAGAAGAAGCAGATGTTGTAATAGTTGGTGCAGGCCCTGCAGGGCTCTCTCAGCT
GTTTCGTCTAAAACAGTTGGCTGTGGCAGATGAAAAGGACATCCGTGTGTCTAGTGGAGAAAGCTGCC
CAGATAGGAGCTCATACTCTCTCAGGGCTTGCCTTGATCCAGGTGCTTTTAAAGAACTCTCCAGAC
TGGAAAGAGAAGGGGGCTCCACTTAACACTCCTGTAACAGAAGACAGATTTGGAATTTTAAACAGAGAAA
TACAGAATTCCTGTGCCAATTCTCCAGGGCTTCCAATGAATAATCATGGCAATTACATTGTACGCTTG
GGACATTTAGTGAGCTGGATGGGCGAACAAGCAGAAGCCCTTGGTGTGAAGTATACCCTGGTTATGCA
GCTGCTGAGGTCCTTTTTCATGATGATGGTAGTGTAAAAGGAATTGCCACTAACGATGTAGGATACAA
AAGGATGGTGCACCAAAGGCAACATTTGAGAGAGGACTGGAACACATGCTAAAGTCACAATTTTTCGA
GAAGTTGCCATGGACATCTAGCCAAGCAACTATAAAGAAGTTTGATTTGAGAGCAAATGTGAACCT
CAAACCTACGGGATTGGACTGAAGGAGTTATGGGTTATTGATGAAAAGAAGCTGGAACCTGGGAGAGTA
GATCACACTGTTGGTTGGCCCTTGACAGACATACCTATGGAGGATCTTCTCTATCATTGAATGAA
GGTGAACCCCTAGTAGCTCTGGTCTTGTGGTGGTCTAGACTATCAGAATCCATACCTGAGTCCATTT
AGAGAGTTCAAAAGGTGGAAACACCATCCTAGCATTCCGGCAACCTTGGAAAGGTGAAAAAGGATTGCA
TACGGAGCCAGAGCTCTCAATGAAGGTGGCTTTCAGTCTATACCAAACTCACCTTCTGGTGGTTTA
CTAATTGGTTGTAGTCTGGTTTTATGAATGTTCCCAAGATCAAAGTACTCACACAGCAATGAAAAGT
GGAATTTTAGCAGCAGAATCTATTTTAACTCACTAAGTGTAAAATCTCCAATCAAAGACAATAGGA
CTCCATGTAAGTGAATATGAGGACAATTTGAAGAACTCATGGGTATGGAAAGAGCTATATTCTGTTAGA
AATATAAGACCGTCCGCCAGGAGTACTGGGTGTATATGGAGGGATGATTTACACTGGAATCTTTTAC
TGGATATTGAGAGGAATGGAGCCGTGGACTCTGAAACATAAAGGTTCTGACTTTGAACGGCTCAAGCCA
GCCAAGGATTGCACACCTATTGAGTATCCAAAACCCGATGGACAGATCAGTTTTGACCTCTTGTCTCT
GTGGCTCTGAGTGGTACTAATCATGAACATGACCAGCCGGCACACTTAACCTTAAGGGATGACAGTATA
CCTGTAATAGAAATCTGTCGATATATGATGGGCCGAGCAGCGATTCTGTCTGCAGGAGTTTATGAA
TTTGTACCTGTGGAACAAGGTGATGGATTTCCGTTACAGATAAATGCTCAGAAGTGTGTACATTGTAAA
ACATGTGATATTAAGATCCAAGTCAGAATATTAAGTGGTGGTACCTGAAGGTGGAGGAGACCTGCT
TACAATGGAATGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001281738

Insert Size: 1671 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001281738.1](#)

RefSeq Size: 1952 bp

RefSeq ORF: 1671 bp

Locus ID: 2110

Cytogenetics: 4q32.1

MW: 61.4 kDa

Gene Summary:

This gene encodes a component of the electron-transfer system in mitochondria and is essential for electron transfer from a number of mitochondrial flavin-containing dehydrogenases to the main respiratory chain. Mutations in this gene are associated with glutaric acidemia. Alternatively spliced transcript variants that encode distinct isoforms have been observed. [provided by RefSeq, Aug 2013]

Transcript Variant: This variant (3) represents use of an alternate promoter, differs in the 5' UTR, and uses a downstream start codon compared to variant 1. The resulting protein (isoform 3) is shorter and has a distinct N-terminus compared to isoform 1.