

Product datasheet for **SC336783**

CNOT6 (NM_001303241) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CNOT6 (NM_001303241) Human Untagged Clone
Tag:	Tag Free
Symbol:	CNOT6
Synonyms:	CCR4; Ccr4a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336783 representing NM_001303241.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCAAAGAAAAATACGAGCCCCCTGACCCTCGGAGGATGTATACAATTATGTCTTCTGAGGAAGCA
GCAAATGGAAAGAAATCCCACTGGGCAGAGCTTGAAATAAGTGGAAAAGTAAGAAGCTTAAGCGCATCT
TTGTGGTCACTAACTCACCTGACAGCTTTGCATTTGAGTGACAATTCCCTGTCCCGAATTCCTTCAGAC
ATTGCCAAGCTTCACAATCTGGTGTATTTGGACCTGTCATCTAATAAAATTCGTAGCTTACCCGCAGAA
CTCGGAAACATGGTATCACTCAGGGAGCTCCATTTAAATAACAACCTGTTACGAGTTCTACCTTTTGAG
CTGGGAAAACCTGTTTCAGTTGCAGACTTTAGGCCTGAAAGGAAATCCCCTTACCCAGGATATATTGAAC
CTTTATCAGGAACCAGATGGAACAAGACGGCTGCTGAACATTTTGCTTGATAATTTGTCAGTACTGCA
AAAAGAATTACAACAGAACAACCACCTCCAAGGCTTGGATTATGTTACAAGAACCAGATAGGACAAGG
CCAACCTGCCTGTTTTCTGTCATGTGCTATAATGTTCTTTGTGATAAATATGCGACCCGGCAGTTATAC
GGCTACTGTCCATCATGGGCGCTAAACTGGGACTACAGGAAAAAGGCCATTATTCAGAAATCTTGAGC
TGCAATGCTGATATCGTAAGTCTTCAGGAGGTTGAAACGGAACAGTATTACAGTTTTTTTTCTGGTAGAG
CTGAAAGAACGTGGCTATAATGGATTCTTCAGTCTAAGTCTAGAGCTAGGACAATGTCAGAAACAAGAA
AGGAAACATGTTGATGGCTGTGCAATATCTTCAAGACAGAAAAATTTACTTTGGTTTCAGAAACACACT
GTTGAATTTAATCAGCTAGCCATGGCAAATTCAGGGGTCTGAAGCTATGCTGAACAGAGTCATGACA
AAAGATAACATTGGGGTTCAGTACTGCTAGAACTTCGGAAGGAATCGATTGAAATGCCGTCCGGAAG
CCACATCTTGAACAGAAAAACAATTTCTTGTGGCTAACGCCACATGCATTGGGACCCTGAATAC
TCTGATGTGAAGTTGGTACAAACTATGATGTTCCCTCTCAGAAGTGAAGAACATTATTGATAAAGCCTCT
CGCAACCTCAAATCCAGTGTGTTGGGAGAATTTGGAATTTCCACTTGTGTTATGTGCAGATCTTAAT
TCTTTGCCAGACTCTGGTGTGTAAGTATTTGAGCACAGGTGGAGTAGAAACAAATCACAAAGACTTT
AAGGAGTTGAGGTATAATGAAAGTCTCACAAACTTCAGCTGTCATGGGAAGAATGGAACCAATGGA
AGGATCACTCATGGTTTCAAGTTACAGAGTGCCTATGAGAGTGGCCTGATGCCTTACACGAATTACACA
TTTGATTTCAAGGTATAATAGACTACATTTTCTATTCTAAACCTCAGCTGAACACCTTAGGCATCCTG
GGCCCTCTGGACCACCACTGGCTGGTTGAGAATAACATCAGTGGCTGCCCGCACCCCTCATCCCCTCT
GACCCTTCTCACTTTTTGCACAACCTGGAGCTTACTGCCTTCTGCCCAAGTCAACGGCATCCAC
CTTCTGGCAGGAGGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_001303241

Insert Size: 1674 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001303241.1](#)

RefSeq Size: 6347 bp

RefSeq ORF: 1674 bp

Locus ID: 57472

Cytogenetics: 5q35.3

Protein Families: Transcription Factors

Protein Pathways: RNA degradation

MW: 63.3 kDa

Gene Summary:

This gene encodes the catalytic component of the CCR4-NOT core transcriptional regulation complex. The encoded protein has a 3'-5' RNase activity and prefers polyadenylated substrates. The CCR4-NOT complex plays a role in many cellular processes, including miRNA-mediated repression, mRNA degradation, and transcriptional regulation. [provided by RefSeq, Dec 2014]