

## Product datasheet for **SC336797**

### FAM222B (NM\_001288632) Human Untagged Clone

#### Product data:


**Product Type:** Expression Plasmids  
**Product Name:** FAM222B (NM\_001288632) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** FAM222B  
**Synonyms:** C17orf63  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC336797 representing NM\_001288632.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCTAGCCTGTCTACCAGGGCCAGGTGACCTGTCTTTTCAGCTTCTTTCTCACACGCAGATGAACACT
GGACTTCAGAAATGGGACACTACACAGAAAATGAGAACTGCTCACTATCTACCCAGCCGAATTGGAT
GCGTATGCTAAGAAGGTCGAAACAACCCACTGACTATAAAAATCTTCCCAACAGTGTGAAGGTTCCC
CAGCGGAAACACGTTTCGTCTGACTGTGAACGGCCTCGACACATCAGCCCAGCGCTACAGCCCTACCCG
ACACAGGCTGCCACCAAGGCAGGCTGCTTGCATTGTCAAAGTGCCAGCCAAAAGCATACTCAAGGAC
TTTGACGGCACCCGAGCCGGTGTCTCCCTGAGGCCATCATGAACCCCCAGTGGCACCCATATGCTACT
GTGGCACCCAGCACTTTAGCCACCCAGGCCAGGCTCTGGCCGCCAGCAGGCCCTGCAGCATGCA
CAGACCCTGGCCCATGCCCTCCCAGACGCTGCAGCACCCCTCAGGGTATCCCGCCACCCAGGCACTG
TCCCACCCTCAGAGCCTCCAGCAGCCTCAGGGCCTGGGCCACCCTCAGCCCATGGCCAAAACCCAGGGC
TTGGTCCACCCTCAGGCCCTGGCTCACCAGGGTCTCCAGCACCCCACAATCCCTTGCTGCATGGAGGC
CGGAAGATGCCAGACTCAGATGCCCCCCGAATGTGACCGTGTCTACCTCAACTATCCCCCTTCAATG
GCGGCCACTCTGCAGCACAGCCAGCCTCCGGACCTGAGTAGCATCGTGACCAGATCAACCAGTTTTGC
CAGACGAGGGCAGGCATCAGCACTACCTCAGTGTGTGAGGGCCAGATCGCAACCCAGCCCCATTAGT
CGCAGCTGCTCATCAATGCAAGCACCCGGGTGTCGACCCACAGGTCACCCACCAATGCCTTCATGT
GTGGTCAATCCCATGGAGCACACCACGGGCCACCGCCGTTGCTGCTGCAGGTCCTGTCAACCTG
CCCACAGGCATCTCTCGGTCCCCTACTGGCTACCTAGCGACCTCAAGCCAGTCACTGGAACAGCAC
CAGCTGGCCACCTACAGCAGATGTGCAGCGAGGCTAGTGGGACGCCAGCCCTGGCCTGACAGGCAAG
CATGCGGCAGGACGCGAGTTGGCAGGGCCTGGCTTTGTGGCAAGGCCCTGCCTACCCGCAAGAACTC
TGCCTGGCGCAGTCTTCCATCTGAAGCCACCCTGAAAAGCCGACCCCATCCCACCAGTCAACGGC
ATGGCAGCCCCATTGGCCTACCCCAATGGTCACTACTTCAACCCCTGTGGAACAACATTCTGCCACT
CCCAATAGCGACAGCTCGGGTCTCAGGACCTTGCATGCCGTTCCACGGTGGGACGCCACAGGTGCA
CCCCTCGACTGTGCGGCAGCTCCCGGGGCCACTACCGAGCAGGGACCGGGGGCGTCCAGTGGCAAGC
CAGAACAGCTTGATGCAAACAGTGGATTACCTAAGTGGGATTTCCAACAGGCCTGCTTCCGAGAACAG
AGCCTGGCCATGCTGAGCAAGGCCACCGAGCCCTGGCAACCGAGCCCCGATCCCACAGAGAGTCGA
AGTCTTCATATTCAGCACCCAGGTATAGTAG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001288632

**Insert Size:** 1689 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001288632.1](#)

**RefSeq Size:** 4261 bp

**RefSeq ORF:** 1689 bp

**Locus ID:** 55731

**UniProt ID:** [Q8WU58](#)

**Cytogenetics:** 17q11.2

**MW:** 59.7 kDa