

Product datasheet for **SC336799**

FAM222B (NM_001288634) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: FAM222B (NM_001288634) Human Untagged Clone
Tag: Tag Free
Symbol: FAM222B
Synonyms: C17orf63
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC336799 representing NM_001288634.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCTAGCCTGTCTACCAGGGCCAGGTGACCTGTCTTTTCAGCTTCTTTCTCACACGCAGATGAACACT
GGACTTCAGAAATGGGACACTACACAGAAAATGAGAACTGCTCACTATCTACCCAGCCGAATTGGAT
GCGTATGCTAAGAAGGTCGAAACAACCCACTGACTATAAAAATCTTCCCAACAGTGTGAAGGTTCCC
CAGCGGAAACACGTTTCGTCTGACTGTGAACGGCCTCGACACATCAGCCCAGCGCTACAGCCCCTACCCG
ACACAGGCTGCCACCAAGCAGGCTGCTTGCATTGTCAAAGTGCCAGCCAAAAGCATACTCAAGGAC
TTTGACGGCACCCGAGCCGGTGTCTCCCTGAGGCCATCATGAACCCCCAGTGGCACCCATATGCTACT
GTGGCACCCAGCACTTTAGCCACCCAGGCCAGGCTCTGGCCGCCAGCAGGCCCTGCAGCATGCA
CAGACCCTGGCCCATGCCCTCCCAGACGCTGCAGCACCCCTCAGGGTATCCCGCCACCCAGGCACTG
TCCCACCCTCAGAGCCTCCAGCAGCCTCAGGGCCTGGGCCACCCTCAGCCATGGCCAAAACCCAGGGC
TTGGTCCACCCTCAGGCCCTGGCTCACCAGGGTCTCCAGCACCCCAACAATCCCTTGCTGCATGGAGGC
CGGAAGATGCCAGACTCAGATGCCCCCCGAATGTGACCGTGTCTACCTCAACTATCCCCTTTCAATG
GCGGCCACTCTGCAGCACAGCCAGCCTCCGGACCTGAGTAGCATCGTGACCAGATCAACCAGTTTTGC
CAGACGAGGGCAGGCATCAGCACTACCTCAGTGTGTGAGGGCCAGATCGCAACCCAGCCCCATTAGT
CGCAGCTGCTCATCAATGCAAGCACCCGGGTGTCGACCCACAGGTCACCCACCAATGCCTTCATGT
GTGGTCAATCCCATGGAGCACACCACGGGCCACCGCCGTTGCTGCTGCAGGTCTGTCAACCTG
CCCACAGGCATCTCTCGGTCCCCTACTGGCTACCTAGCGACCTCAAGCCAGTCACTGGAACAGCAC
CAGCTGGCCACCTACAGCAGATGTGCAGCGAGGCTAGTGGGACGCCAGCCCCTGGCCTGACAGGCAAG
CATGCGGCAGGACGCGAGTTGGCAGGGCCTGGCTTTGTGGCAAGGCCCTGCCTACCCGCAAGAACTC
TGCCTGGCGCAGTCTTCCATCTGAAGCCACCCTGAAAAGCCGACCCCATCCCACCAGTCAACGGC
ATGGCAGCCCCATTGGCCTACCCCAATGGTCACTACTTCAACCCTGTGGAACAACATTCTGCCACT
CCCAATAGCGACAGCTCGGGTCTCAGGACCTTGGCATGCCGTTCCACGGTGGGACGCCACAGGTGCA
CCCCTCGACTGTGCGGCAGCTCCCGGGGCCACTACCGAGCAGGGACCGGGGGCGTCCAGTGGCAAGC
CAGAACAGCTTGATGCAACAGTGGATTACCTAAGTGGGATTTCCAACAGGCCTGCTTCCGAGAACAG
AGCCTGGCCATGCTGAGCAAGGCCACCGAGCCCTGGCAACCGAGCCCCGATCCCACAGAGAGTCCA
AGTCTTCATATTCAGCACCCAGGTATAGTAG
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001288634

Insert Size: 1689 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288634.1](#)

RefSeq Size: 4248 bp

RefSeq ORF: 1689 bp

Locus ID: 55731

UniProt ID: [Q8WU58](#)

Cytogenetics: 17q11.2

MW: 59.7 kDa