

Product datasheet for **SC336804**

CRNKL1 (NM_001278626) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CRNKL1 (NM_001278626) Human Untagged Clone
Tag:	Tag Free
Symbol:	CRNKL1
Synonyms:	CLF; Clf1; CRN; HCRN; MSTP021; SYF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336804 representing NM_001278626.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAAATGAAGAATCGCCAAGTCAACCATGCTCGAAATATCTGGGACCGGGCCATAACAACGCTGCCT
CGAGTTAATCAGTCTGGTACAAGTACACGTACATGGAGGAAATGTTGGGAAACGTTGCCGGTGCCCGG
CAGGTGTTTGAGCGCTGGATGGAGTGGCAGCCTGAGGAGCAAGCCTGGCACTCCTACATCAACTTTGAG
CTGAGATACAAAGAGGTGGATCGGGCCCGCACCATTTATGAGCGATTTGCTCGTGACCCCTGATGTT
AAGAAGTGGATCAAGTATGCCCGCTTTGAAGAAAAACATGCTTATTTGCCATGCACGAAAGTGTAT
GAGAGAGCTGTGGAATCTTTGGAGATGAACATATGGATGAGCACCTTTATGTTGCCTTTGCCAAGTTT
GAAGAAAAACAGAAAGAGTTTGAAGGGTACGAGTGATTTACAAGTATGCCCTGGACAGAATTTCAAAA
CAAGATGCCCAAGAAGCTTTAAAAATTATACCATCTTTGAGAAGAAGTTTGGTGTATAGGCGGGGTATT
GAAGATATCATTGTGAGCAAACGGAGATCCAGTACGAAGAAGAAGTGAAGGCGAATCCACACAATTAT
GATGCATGGTTTGATTACTTGCCTTGGTGAAGAGTACGCAGAAGCTGAAGCCGTGAGAGAAGTCTAT
GAAAGGGCCATTGCCAATGTCCACCCATTGAGGAGAAGAGGCACTGGAAGCGCTACATTTATCTTTGG
ATCAACTATGCACTCTATGAAGAATTGGAGGCAAAGGATCCTGAGAGGACAAGACAGGTGTATCAAGCC
TCTTTGGAATAATTCTCACAAAAAGTTACATTTGCCAAAATGTGGTACTGTATGCACAGTTTGAA
ATACGACAGAAGAATCTGTCATTAGCCAGAAGAGCATTGGGAAGTCCATAGGCAAAATGTCCAAAGAAC
AAATTATTTAAAGTTTACATAGAATTGGAGCTACAGCTTCGAGAATTTGACAGATGCCGGAAGCTTTAT
GAAAAGTTCTGGAATTTGGACCTGAAAATTGTACCTCATGGATTAATTCGCTGAATTAGAGACAATC
CTTGGTGTATTGACAGAGCACGGGCAATCTATGAATTAGCCATCAGTCAGCCACGTTTAGACATGCCA
GAGGTGCTTTGGAATCATATATTGATTTTGAATTGAGCAGGAAGAAACAGAAAGAACCGAAACCTT
TACCGGGCGGTTGCTTCAACGGACGAGCATGTCAAGGTATGGATCAGCTTTGCTCAGTTTGAGTTGTCT
TCAGGAAAAGAAGGAAGTTTGAATAATGCAGACAAATTTATGAAGAAGCTAACAAAACCATGCGAAAC
TGTGAAGAAAAGGAGAGAGACTTATGCTGCTGGAATCTTGGCGAAGTTTGAAGAAGAATTTGGAACA
GCTTCAGATAAGGAGAGAGTAGACAAACTCATGCCAGAGAAAGTCAAGAAGAGAAGAAAGGTCCAGACT
GATGATGGGTCTGATGCAGGCTGGGAAGAATACTTTGATTACATCTTTCCAGAAGATGCTGCCAACCA
CCTAACCTCAAACCTCTGGCCATGGCCAAACTGTGGAAGAAACAGCAGCAGGAAAAGGAGGATGCTGAG
CACCATCCAGATGAGGACGTCGATGAGAGTGAATCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001278626

Insert Size: 1695 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278626.1](#)

RefSeq Size: 4190 bp

RefSeq ORF: 1695 bp

Locus ID: 51340

UniProt ID: [Q9BZJ0](#)

Cytogenetics: 20p11.23

Protein Pathways: Spliceosome

MW: 68.6 kDa

Gene Summary:

The crooked neck (crn) gene of *Drosophila* is essential for embryogenesis and is thought to be involved in cell cycle progression and pre-mRNA splicing. A protein encoded by this human locus has been found to localize to pre-mRNA splicing complexes in the nucleus and is necessary for pre-mRNA splicing. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jul 2013]

Transcript Variant: This variant (3) uses two alternate splice sites in the 5' UTR, compared to variant 1. These differences result in predicted translation initiation at a downstream in-frame start site through ribosomal re-initiation and/or leaky scanning. The encoded isoform (c) is shorter than isoform a. Both variants 3 and 4 encode the same isoform (c).