

## Product datasheet for **SC336806**

### **FAM222B (NM\_001288631) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	FAM222B (NM_001288631) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM222B
Synonyms:	C17orf63
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336806 representing NM\_001288631.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGTCAAGCTTCCCCCTTGCCTCTGCTAAAAAAACCTCAGAGTTGGTACTAACTTGAACAAAGG
GACTCTTCTACTCATCAAGGGGACACTACACAGAAAATGAGAAGTCTCACTATCCTACCCAGCCGAA
TTGGATGCGTATGCTAAGAAGGTCGAAACAACCCACTGACTATAAAAAATCTTCCCAACAGTGTGAAG
GTTCCCAAGCGGAAACACGTTCTGCTGACTGTGAACGGCCTCGACACATCAGCCAGCGCTACAGCCCC
TACCCGACACAGGCTGCCACCAAGGCAGGCCTGCTTGCCATTGTCAAAGTGCCAGCCAAAAGCATACTC
AAGGACTTTGACGGCACCCGAGCCCGTTGCTCCCTGAGGCCATCATGAACCCCCAGTGGCACCCTAT
GCTACTGTGGCACCCAGCACTTTAGCCACCCAGGCCAGGCTCTGGCCCGCCAGCAGGCCCTGCAG
CATGCACAGACCTGGCCATGCCCTCCCGACGCTGCAGCACCTCAGGGTATCCCGCCACCCAG
GCACTGTCCACCCTCAGAGCCTCCAGCAGCCTCAGGGCCTGGGCCACCCTCAGCCCATGGCCAAAAC
CAGGGCTTGGTCCACCCTCAGGCCTGGCTCACCAGGGTCTCCAGCACCCCAACAATCCCTTGCTGCAT
GGAGGCCGGAAGATGCCAGACTCAGATGCCCCCCGAATGTGACCGTGTCTACCTCAACTATCCCCCTT
TCAATGGCGGCCACTCTGCAGCACAGCCAGCCTCCGGACCTGAGTAGCATCGTGACCCAGATCAACCAG
TTTTGCCAGACGAGGGCAGGCATCAGCACTACCTCAGTGTGTGAGGGCCAGATCGCAACCCAGCCCC
ATTAGTCGAGTCTGCTCATCAATGCAAGCACCCGGGTGTGACCCACAGCGTCCCAACCAATGCCT
TCATGTGTGGTCAATCCCATGGAGCACACCCACGGGCCACCGCCGCTTGCCTGCTGCAGGTCTGTG
AACCTGCCACAGGCATCTCTCGGTCCCCTGGCTACCCTAGCGACCTCAAGCCAGTCACTGGAAC
CAGCACCAGTGGCCACCTACAGCAGATGTGCAGCGAGGCTAGTGGGACGCCAGCCCTGGCCTGACA
GGCAAGCATGCGCAGGACGCGAGTTGGCAGGGCCTGGCTTTGTGGCAAGGCCCTGCTACCCGCGAG
GAACTCTGCCTGGCGCAGTCTTCCATCTGAAGCCACCCCTGGAAAAGCCGACCCATCCCAACAGTC
AACGGCATGGCAGCCCCATTGGCCTACCCCAATGGTCACTACTTCCAACCCCTGTGGAACAACATTCTG
CCCCTCCCAATAGCGACAGCTCGGGTCTCAGGACCTTGCCATGCCGTTCCACGGTGGCAGCCACA
GGTGCACCCCTCGACTGTGCGGAGCTCCCGGGGCCACTACCGAGCAGGGACCGGGGGCGGTCCAGTG
GCAAGCCAGAACAGCTTGATGCAACAGTGGATTACCTAAGTGGGATTTCACAGGCCTGCTCCGA
GAACAGAGCCTGGCCATGCTGAGCAAGGCCACCGAGCCCTGGCAACCGAGCCCCGATCCACAGAG
AGTCGAAGTCTTCATATTCAGCACCCAGGGTATAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001288631

**Insert Size:** 1695 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001288631.1](#)

**RefSeq Size:** 4350 bp

**RefSeq ORF:** 1695 bp

**Locus ID:** 55731

**UniProt ID:** [Q8WU58](#)

**Cytogenetics:** 17q11.2

**MW:** 59.9 kDa