

Product datasheet for **SC336808**

ZNF286 (ZNF286A) (NM_001288642) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF286 (ZNF286A) (NM_001288642) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF286A
Synonyms:	ZNF286
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336808 representing NM_001288642.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGAGAGCCCAGGAGGGAAGGCAGGAAGTGGGACCAGAGGCCGGGGGGCTTCTAACTCGTCTGGC
CAGGGAGATCTGAATTGGGGTGAAGAGCAGAATCTCCAGAACAAGGAGGAGGTGGTGATCATGGAGACA
GATTTGGCTGAAATGCCTGAGAAAAGGAGCTCTGTCTTCCAGGATTCTCCCCATTTCCAAGAGAAGAGC
ACAGAAGAGGGAGAAGTGGCTGTCTGCGCCTCACGGCCAGATCCAGGAAAACAGTGACATTCAAGGAT
GTGGCCATGGACTTTACACCAGAGGAGTGGGGGAAGCTGGATCCTGCACAAAGGGATGTGATGCTGGAG
AACTATAGGAACCTAGTCTCACTTTGGCTTCCAGTTTCAAACCTGAGAGCTACAACCTGGAGAATGGA
AAAGAACCATTGAAGCTTGAGAGAAAAGCCCCAAAAGCAGCTATTCAGACATGGAGACTAGACCACAG
AGCAAGGATTCAACTTCAGTCAAGATTTTCAAAGCAGAATCATGCAAAGTTGCAATAATAGACAGA
CTGACACGGAATAGTGTCTATGACTCTAACTTGAAGCAGCCCTTGAATGTGAAAATTGGTTAGAGAAT
CAGCAAGGAAATCAGGAGAGACATTTGAGAGAAATGTTCACTCACATGAATTCACCTCTGAGGAAACA
GACCATAAGCATGATGTATACTGGAAAAGCTTCAATCAGAAATCTGTCTTATCACTGAAGACAGAGTT
CCCAAAGGATCTTATGCCTTCCATACACTTGAAAAAGCTTGAAACAAAATCAAACCTAATGAAAAAG
CAGAGAACTTATAAAGAGAAAAACCTCATAAATGTAATGATTGTGGTGAACCTTTCACCTACCATTCA
GTGCTTATTCGACACCAGAGAGTCCATACTGGAGAGAAACCCTATACCTGCAATGAATGTGGGAAATCT
TTAGCCACAGAGCTAATTTAACTAAACACCAGAGAACTCATACTAGAATTCTTTTGAATGCAGTGAA
TGCAAGAAAACCTTCACAGAAAGCTCATCCCTTGCAACACATCAGAGAATTCACGTTGGAGAGAGACCT
TATGAATGCAATGAATGTGGGAAAGTTTTAATCGAAGTACACATCTTGTGCAGCATCAGTTGATTCAT
ACTGGAGTGAAGCCTTATGAATGAATGAATGTGATAAAGCTTTTATTCATTCATCAGCACTCATTAAA
CATCAAAGAACTCATACTGGAGAGAAACCTTATAAATGTCAAGAATGTGGGAAAGCCTTTAGCCATTGC
TCATCCCTAACTAAGCATCAGAGAGTTCATACTGGAGAAAAGCCATATGAATGCAGTGAATGTGGAAAA
ACTTTTAGTCAGAGCACACATCTTGTTCAACATCAGAGAATTCACACTGGAGAGAAAACCTATGAGTGT
AATGAATGTGGGAAAACCTTCAGCCGAGCTCCAATTTTGCTAAACATCAAAGAATTCATATTGGAAAG
AAACCGTACAAATGTAGCGAGTGTGGAAAAGCCTTCATTCATTCATCAGCTCTCATTCAACATCAGAGA
ACTCATAACCGGAGAGAAAACCTTAGATGTAATGAGTGTGGGAAAAGCTTAAAGTGCAGTTCATCTCTC
ATCAGACATCAAAGAGTTCACACTGAAGAGCAACCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001288642

Insert Size: 1695 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288642.1](#)

RefSeq Size: 5399 bp

RefSeq ORF: 1695 bp

Locus ID: 57335

UniProt ID: [Q9HBT8](#)

Cytogenetics: 17p12

Protein Families: Transcription Factors

MW: 65.1 kDa

Gene Summary:

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) encodes the longest isoform (1).