

Product datasheet for **SC336812**

ZNF26 (NM_001256279) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF26 (NM_001256279) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF26
Synonyms:	HEL-179; KOX20
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336812 representing NM_001256279.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCACCAGTTTCCGGACAGCTTCGTGCTGGGTAGTGCTGCTGAGGAGGAGCCTGCTAATGAGCTCA
AGTGTCAAGTTCGGGAATCTTCTGGGTGTTCTTGGGCCAGGGAAACAGTGCCCTGCCTGAGGGATTA
TTGTCATTCAAGGATATATCTATGGAGTTCACCTGGGATGAATGGCAGCTACTGGATTCTACACAGAAG
TACCTGTACAGAGATGTGATATTGGAACCTATCATAACCTGATATCAGTGGGGTATCATGGTACCAAG
CCTGACTTAATCTTCAAGTTGGAACAAGGAGAAGATCCATGGATAATAAATGCCAAATTTCCAGGCAG
AGCTGTCCAGATGGCTGGGAAGATGGTACCAGAACAATCAAGATGAGCTTGAGAGTATTGAAAGAAGC
TATGCTTGTAGTGTGGGAAGACTTAATCTGAGCAAAACCCATGATTCTTCAAGACAGAGACTCTAT
AACACACGTGGAAAAAGTTTGACACAAACTCAGCTCCAAGCAGAAGTTATTTAAGAAAGAATCCTGAT
AAGTTTCATGGTTATGAAGAACCATAATTTCTTAAGCATCAAAGAGCTCATAGCATAGAAAAAACTGT
GTGTGTAGTGAATGTGGGAAAGCTTTTCTGTTGAAGTCACAGCTCATTGTACATCTCAGAATTCATACA
GGAGAGAGACCTTATGAATGCAGTAAATGTGAAAGAGCCTTCAGTGCCAAGTCAAACCTTAATGCTCAT
CAGAGAGTTCATACAGGAGAAAAACCTACTCATGTAGTGAGTGCAGAGAAGGTCTTCTTTTCAGGTCA
CAGCTCATTGTCCATCAGGAAATTCACACAGGAGGAAACCTATGGCTGCAGTGAATGTGGGAAAGCC
TACAGTTGGAAATCACAGCTTCTTTACACCAGAGAAGTCACACAGGAGTGAAACCGTATGAATGCAGC
GAATGTGGGAAAGCCTTTAGTTTGAAGTCTCCATTCGTTGTACACCAGAGAAGTCATACAGGAGTGAAA
CCCCATAAATGCAGTGAATGTGGGAAAGCCTTTAGGAGTAAGTCCTATCTCCTTGTTACATCCGAATG
CATACAGGAGAAAAACCTATCAATGCAGTGATTGTGGGAAAGCCTTCAATATGAAGACACAACCTATT
GTACATCAGGGAGTTCACACAGGAAATAATCCTTATCAATGCGGTGAATGTGGGAAAGCCTTTGGTAGG
AAGGAACAGCTCACTGCACATCTGAGAGCTCATGCAGGAGAGAAGCCCTATGGATGCAGTGAATGTGGG
AAGGCTTTTCAGCAGCAAGTCATACCTTGTATACATAGGAGAACACACACCGGAGAGAGACCTATGAA
TGAGTTTGTGTGAGAGAGCCTTTTGTGGAATACACAGCTGATTATACATCAGAGAAGTCACTCAACT
GAGAAGCCCTATGAATGCAATGAATGTGAAAAAGCCTACCCTAGGAAGGCATCACTTCAGATACACCAG
AAAACCTATTCCGGGAGAGAAACCTTTAAATGCAGTGAATGTGGGAAAGCCTTCACTCAGAAGTCATCT
CTCAGTGAACATCAGAGAGTTCACACCGGAGAGAAACCATGGAATGCTCTGAATGTGGGAAATCCTTC
TGTTGGAATTCAGGGCTTCGTATACATCGGAAGACTCATAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001256279

Insert Size: 1701 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256279.1</u>
RefSeq Size:	3297 bp
RefSeq ORF:	1701 bp
Locus ID:	7574
Cytogenetics:	12q24.33
Protein Families:	Transcription Factors
MW:	64.8 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.