

Product datasheet for **SC336813**

HSZFP36 (ZNF823) (NM_017507) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HSZFP36 (ZNF823) (NM_017507) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF823
Synonyms:	HSZFP36
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336813 representing NM_017507.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAATGGGAGGACCAGAACATTGGAGATCAGTGCCAAAATGCCAAGAGAAATCTAAGAAGTCATACA
TGTGAAATTAAGATGACAGTCAATGTGGAGAACTTTTGGCCAGATTCCAGATAGTATTGTGAACAAG
AACACTCCTCGAGTAAATCCATGTGACAGTGGTGTGAGTGTGGAGAAGTCGTCTTGGGTCATTTCGTCTCTT
AATTGCAACATCAGAGTTGACACTGGACACAAATCATGTGAGCATCAGGAATATGGAGAGAAGCCATAT
ACACATAAAACAACGTGGGAAAGCCATCAGTCATCAGCACTCCTTCCAGACACATGAAAGGCCCCACC
GGAAAGAAACCCTTCGATTGTAAAGAATGTGCAAAAACCTTTAGTTCTCTTGGAAACCTCCGAAGACAC
ATGGCGGCACACCATGGAGATGGACCTTATAAATGTAAGTTGTGTGGGAAAGCCTTTGTTTGGCCAGT
TTATTTCAATTTGCACGAAAGAACACACTGGAGAGAAACCGTATGAATGTAAGCAGTGTCTAAAGCC
TTTCCTTTTACAGTTCCTATCTAAGACATGAAAGAATCCACACGGGAGAGAAAGCGTATGAATGTAAG
CAGTGTCCAAAGCCTTTCCTGATTACAGTACCTATCTAAGACATGAGAGAACTCACACCGGAGAGAAA
CCCTATAAATGTACACAATGTGGGAAAGCCTTCAGCTGTACTATTACACTCGACTACATGAAAGGACT
CACACGGGAGAACAACCTATGCATGTAAGCAATGTGGGAAAACGTTTTATCATCACACAAGCTTTCGA
AGACACATGATAAGGCACACTGGAGACGGACCACATAAATGTAAGATATGTGGGAAAGGCTTTGATTGT
CCTAGTTCAGTTCGAAATCATGAACTACTCACACTGGAGAGAAACCTATGAATGTAAGCAGTGTGGG
AAAGTGTATCTCATAGCTCGAGCTTTCGAAAGTACATGATAACACACACAGGAGATGGACCCAGAAA
TGCAAGATATGTGGGAAAGCCTTTGGTTGTCCAGTTTATTTCAAAGACATGAAAGGACTCACACTGGA
GAGAAACCCTATCAATGTAACAATGTGGTAAAGCCTTCAGTCTTGCCGTTCCCTTCGAAGACATGAA
GCAACTCACACTGGAGTAAACCTATAAATGTCAGTGTGGGAAAGCCTTTAGTGATCTCTCTTCCCTT
CAAAATCATGAGACAACCTCACACTGGAGAGAAAGCCATATGAGTGTAAAGGAATGTGGGAAAGCATTAGT
TGTTTCAAATACCTTTCTCAACATAAAAGGACCCACACAGTAGAAAAACCTTATGAGTGTAAAAATGT
AGAAAAGCCTTCAGTCATTTTCAGTAACTTAAAAGTCCATGAAAGGATTCCTCTGGAGAGAAGCCATAT
GAATGTAAGGAATGTGGAAAAGCATTCTTGGCTCACTTGCCTTCTACGACATGAAAGAATTCACACT
GGAGAGAAACCCTATGAATGTCTACAATGTGGTAAAGCCTTCAGTGTCCCGTTTCCCTTCGAGGACAT
GAAAAAATCACACTGGAGAGAAGCTGTATGAATGTAAGGAATGTGGGAAAGCATTGAGTCTCTCCGT
TCCTTGCATAGACATAAAAGGACTCACTGGAAAGATACTCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_017507

Insert Size: 1701 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017507.1](#)

RefSeq Size: 2351 bp

RefSeq ORF: 1701 bp

Locus ID: 55552

UniProt ID: [P16415](#)

Cytogenetics: 19p13.2

MW: 65.2 kDa

Gene Summary:

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream start codon, compared to variant 1. It encodes isoform 2, which has a shorter N-terminus, compared to isoform 1.