

Product datasheet for **SC336824**

MPP2 (NM_001278376) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MPP2 (NM_001278376) Human Untagged Clone
Tag:	Tag Free
Symbol:	MPP2
Synonyms:	DLG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336824 representing NM_001278376.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGGCAGCCAGGCAGCGGGTCTCCTTGGAGGGTATATCCCTGGAGTCTTCTGAGGAAGCGGAG
CTCCAGAGGGAAGCCATGCAGCAAGTCTGGACAACCTGGGATCCCTCCCAGTGCCACGGGGCTGCA
GAGCTGGACCTGATCTTCTTCGAGGCATTATGGAAAGTCCCATAGTAAGATCCCTGGCCAAGGCCAT
GAGAGGCTGGAGGAGACGAAGCTGGAGGCCGTGAGAGACAACAACCTGGAGCTGGTGCAGGAGATCTG
CGGGACCTGGCGCAGCTGGCTGAGCAGAGCAGCACAGCCGCCGAGCTGGCCACATCCTCCAGGAGCCC
CACTTCCAGTCCCTCCTGGAGACGCAGACTCTGTGGCCTCAAAGACCTATGAGACACCACCCCCAGC
CCTGGCCTGGACCCTACATTAGCAACCAGCCTGTACCTCCCGATGCTGTGCGCATGGTGGGCATCCGC
AAGACAGCCGGAGAACATCTGGGTGTACGTTCCGCGTGGAGGGCGGCGAGCTGGTATCGCGCGCATT
CTGCATGGGGCATGGTGGCTCAACAAGGCCTGCTGCATGTGGGTGACATCATCAAGGAGGTGAACGGG
CAGCCAGTGGGCAGTGACCCCCGCGCACTGCAGGAGCTCCTGCGCAATGCCAGTGGCAGTGTATCCTC
AAGATCCTGCCAGCTACCAGGAGCCCCATCTGCCCGCCAGGATTTTGTGAAATGTCACTTTGACTAT
GACCCGGCCCGAGACAGCCTCATCCCCTGCAAGGAAGCAGGCCTGCGCTTCAACGCCGGGACTTGCTC
CAGATCGTAAACCAGGATGATGCCAACTGGTGGCAGGCATGCCATGTCAAGGGGGCAGTGTGGGCTC
ATTCAGCCAGCAGCTGCTGGAGGAGAAGCGGAAAGCATTGTCAAGAGGGACCTGGAGCTGACACCAAAC
TCAGGGACCCTATGCGGCAGCCTTTCAGGAAAGAAAAAGAAGCGAATGATGTATTTGACCACCAAGAA
GCAGAGTTTGACCGTATGAGCTGCTCATTTATGAGGAGGTGGCCCGCATGCCCCGTTCCGCGGAAA
ACCCTGGTACTGATTGGGGCTCAGGGCGTGGGACGGCGCAGCCTGAAGAACAAGCTCATCATGTGGGAT
CCAGATCGCTATGGCACCACGGTGCCCTACACCTCCCGCGGCCGAAAGACTCAGAGCGGGAAGGTGAG
GGTTACAGCTTTGTGCCGTGGGGAGATGGAGGCTGACGTCCGTGCTGGGCGCTACCTGGAGCATGGC
GAATACGAGGGCAACCTGTATGGCACACGTATTGACTCCATCCGGGGCGTGGTCCGCTGCTGGGAAGGTG
TGCGTGTGGATGTCAACCCCCAGGCGGTGAAGGTGCTACGAACGGCCGAGTTTGTCCCTTACGTGGTG
TTCATCGAGGCCCCAGACTTCGAGACCCTGCGGGCCATGAACAGGGCTGCGCTGGAGAGTGGAAATATCC
ACCAAGCAGCTCACGGAGGCGGACCTGAGACGGACAGTGGAGGAGAGCAGCCGCATCCAGCGGGGTAC
GGGCACTACTTTGACCTGCTGCTGGTCAATAGCAACCTGGAGAGGACCTCCGCGAGCTCCAGACAGCC
ATGGAGAAGCTACGGACAGAGCCCCAGTGGGTGCCTGTCAGCTGGGTGTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001278376

Insert Size: 1710 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278376.2](#)

RefSeq Size: 4275 bp

RefSeq ORF: 1710 bp

Locus ID: 4355

UniProt ID: [Q14168](#)

Cytogenetics: 17q21.31

Protein Families: Druggable Genome

MW: 63.3 kDa

Gene Summary:

Palmitoylated membrane protein 2 is a member of a family of membrane-associated proteins termed MAGUKs (membrane-associated guanylate kinase homologs). MAGUKs interact with the cytoskeleton and regulate cell proliferation, signaling pathways, and intracellular junctions. Palmitoylated membrane protein 2 contains a conserved sequence, called the SH3 (src homology 3) motif, found in several other proteins that associate with the cytoskeleton and are suspected to play important roles in signal transduction. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (8) differs in the 5' UTR and represents use of an alternate promoter, compared to variant 1. It encodes isoform 8 which has a distinct N-terminus and is shorter, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.